



Yes, you can:

Tel: 0800 043 0881
Email: salesUK@dhcare.com
Online Spares available at: www.dhcare.co.uk

Power Retail Prescription Form

September 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 NB Modulite

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 20°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- Standard 80° hangers with fixed footplate
- LiNX DLX-REM110 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Stability Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 50 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,995



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code **Description**

BASE MODEL*

ARM0000	☐	TDX SP2 Narrow Base	£6,995
ARZ1300	☐	Special TDX SP2 Narrow Base (Extended Lead Time). Please provide details in comments section.	POA

USER WEIGHT*

ARM1920	☐	User weight up to 125kg	☐
ARM1921	☐	User weight 126kg to 136kg	☐
ARM1922	☐	User weight 137kg to 150kg (NC 1092, 2052, 2053, 2055)	☐

SEATING TYPE*				
ARM2100		Modulite Seat Unit	(NC 2400, 2401, 2402, 2403, 4520, 4521, 2410, 2411, 2412, 2328, 2329, 2330, 2331, 4510, 4530, 4540, 4550)	◆
ARM2103		Modulite Support+	(NC 2132, 2140, 2142, 2143, 2144, 2147, 2148, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2163, 2170, 2171, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2179, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2155, 2195, 2196, 2197, 2199, 2184, 2186, 2188, 2185, 2187, 2189, 1350, 1427, 1432, 1495, 1496, 2324, 2325, 2327)	○
SEAT FRAME*				
ARM2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 550 mm. Additional 50 mm width by arm rest adjustment.		◆
SEAT WIDTH*				
ARZ1001		Seat width 330 mm. Seat width achieved on a 380 mm seat frame (2121) with inrigged armrests (2183)	(NC 2181, 2193, 2192, 2194, 2148)	○
ARM2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.	(NC 2143, 2701, 2703, 2705, 2148)	○
ARM2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.	(NC 2703, 2705, 2148)	○
ARM2123		480 mm Seat Plate	(NC 2702, 2704, 2706, 2148)	○
ARM2124		530 mm Seat Plate	(NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706, 2148)	○
ARM2125		550 mm Seat Plate	(NC 2132, 2134, 2142, 2143, 2147, 2144, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2170, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2152, 2155, 2183, 2197, 2184, 2186, 2188, 2185, 2187, 2189, 1496, 0556, 0589, 2701, 2702, 2703, 2704, 2705, 1404, 1405)	○
BODY FIT KIT				
UNZ2106		Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)	(NC 2124, 2192, 2194)	○
SEAT DEPTH RAIL LENGTH*				
ARM2129		Standard Depth Package (410mm - 510mm. Simple stepless adjustment from the side)	(NC 2138)	○
ARM2130		Extended Depth Package (510mm - 560mm. Simple stepless adjustment from the side)	(NC 2132, 2134)	○
SEAT DEPTH*				
ARM2132		410 mm	(NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706)	○
ARM2134		460 mm		○
ARM2136		510 mm		○
ARM2138		560 mm	(NC 2143)	○
SEAT OPTION				
ARZSWO		"H" Frame Interface Kit		£495
SEAT CUSHION*				
ARM2142		Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£120
ARM2143		Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm, up to 125kg)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£169
ARM2147		Matrx Libra (90mm - standard up to 136kg or HD up to 227kg). Extended lead time.		£410
ARM2144		Matrx Vi Cushion (Mesh fabric)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	£400
ARM2148		Standard Cushion (90mm thick) - For 550mm Seat Plate Only	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○
ARM2149		Support+ Cushion		£270
ARM2140		Without seat cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○
1st POSITIONING ACCESORIES FOR MATRX LIBRA				
ARM1980		Tapered Wedges (Pair)		£30
ARM1981		Abductor		£23
ARM1982		Lateral Wedge		£20
ARM1983		Anterior Seat Wedges		£27

2nd POSITIONING ACCESORIES FOR MATRX LIBRA			
ARM1980		Tapered Wedges (Pair)	£30
ARM1981		Abductor	£23
ARM1982		Lateral Wedge	£20
ARM1983		Anterior Seat Wedges	£27
POSITIONING ACCESSORIES FOR MATRX VI			
ARM1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	£54
BACKREST SUPPORT*			
ARM2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	(NC 2172, 2173, 2052, 2053) ○
ARM2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	◆
PUSH HANDLES FOR BACKPOSTS*			
ARM2107		Backposts WITH push handles	◆
ARM2108		Backposts WITHOUT push handles. Push handles will be cut down on the vertical section.	(NC 1332, 2219, 2220, 2227, 2228, 2245) ○
BACKREST UPHOLSTERY / CUSHION*			
ARM2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 1445, 1496, 2155, 2052, 2053) ◆
ARM2179		Flex3 Tension Adjustable (540 mm or 480 mm)	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2128) £155
ARM2172		Flex3 adjustable Backrest Plates (without Backrest Cushion)	(NC 2163, 2128, 2151, 2152, 2155) £250
ARM2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion	(NC 2163, 2128) £400
ARM2400		Support+ Standard Backrest (Spacer fabric)	(NC 2152, 2155) £330
ARM2401		Support+ Active Strapped Backrest (Spacer fabric)	(NC 2152, 2155) £400
ARM2402		Support+ Passive Strapped Backrest (Spacer fabric)	(NC 2152, 2155) £400
ARM2403		Support+ Extra Deep Backrest (Spacer fabric)	(NC 2152, 2155) £450
ARM2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £353
ARM2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £353
ARM2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £494
ARM2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £496
ARM2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £546
ARM2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445) £475
ARM2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445) £575
ARM2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 1496, 2151, 2152, 2155) ○
SUPPORT+ BACKREST HEIGHT (Only applicable for Support+ Back options)			
ARM4520		Tall - 570 mm (Measured from Seat Plate)	○
ARM4521		Short - 500 mm (Measured from Seat Plate)	○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)			
ARM2380		405 mm (16 inch) MAC BH	(NC 1445) ○
ARM2381		455 mm (18 inch) MAC BH	(NC 1445) ○
ARM2382		505 mm (20 inch) MAC BH	(NC 1445) ○
ARM2383		560 mm (22 inch) MAC BH	(NC 1445) ○

MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)				
ARM2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445)	£100
ARM2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445)	£100
ARM2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445)	£100
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)				
ARM2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445)	£100
ARM2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445)	£100
ARM2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445)	£100
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)				
ARM2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445)	○
ARM2394		Backrest matches seat width	(NC 1445)	○
ARM2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445)	○
UPHOLSTERY COLOUR*				
ARM2151		Black Nylon Fabric TR18		◆
ARM2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26	(NC 2171)	○
ARM2155		Black Spacer Fabric TR36	(NC 2171, 2179)	○
ARX1002		No Upholstery Colour required - No backrest upholstery required	(NC 2179, 2173, 2174, 2175, 2176, 2177, 2178)	○
ARMRESTS*				
ARM2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm		◆
ARM2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm		£40
ARM2183		Inrigged Armrests. Reduces SW by 25 mm each side	(NC 2185, 2186, 2187, 2188, 2189)	£30
ARM2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).	(NC 1445, 2179)	£65
ARM2192		Following Armrest (height adjustable from 230mm - 300mm).	(NC 1445, 2179)	£185
ARM2194		Following Armrest (height adjustable from 300mm - 360mm).	(NC 1445, 2179)	£185
ARMPADS*				
ARM2195		Short Armpad (325mm)		◆
ARM2196		Long Armpad (375mm)		○
ARM2197		Long Ergonomic Armpad (397mm)		○
ARM2410		Support+ Standard Armpad Short (80 mm x 300 mm)		£100
ARM2411		Support+ Standard Armpad Long (80 mm x 400 mm)		£110
ARM2412		Support+ P-Shaped Armpad (120 mm x 400 mm)		£130
ARM2199		Soft Armpad (REA Azalea Wide Pad)	(NC 2179)	£95
HIP SUPPORTS LH				
ARM2184		Hip Support LH, pad size 140 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
ARM2186		Hip Support LH, pad size 230 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
ARM2188		Hip Support LH, pad size 330 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
HIP SUPPORTS RH				
ARM2185		Hip Support RH, pad size 140 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
ARM2187		Hip Support RH, pad size 230 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
ARM2189		Hip Support RH, pad size 330 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
MOUNTING HARDWARE FOR HIP SUPPORTS				
ARM2164		Mounting hardware for hip supports (6.1" / 155mm high)		○
ARM2165		Mounting hardware for hip supports (8.1" / 206mm high)		○

HEADREST / NECKREST (Adequate headrest option must be selected)*				
ARM1350		Neckrest depth adjustable		£80
ARM1432		Headrest sideways adjustable		£225
ARM1427		Headrest with cheek support		£225
ARM2328		Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).		£250
ARM2329		Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem (_2332).		£270
ARM2330		Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).		£230
ARM2331		Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem (_2332).		£250
ARM2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware		£325
ARM2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware		£325
ARM2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware		£440
ARM1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	(NC 2171, 2174, 2175, 2176, 2177, 2178)	○
ARM1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts		○
SUPPORT+ ACCESSORIES				
ARM4510		Support+ Wipedown Cushion Cover (PU Coated Type)		£30
ARM4530		Support+ Wipedown Backrest Cover (PU Coated Type)		£30
ARM4540		Support+ Wipedown Armpad Covers (PU Coated Type)		£30
ARM4550		Support+ Wipedown Headrest Cover (PU Coated Type)		£30
LATERAL SUPPORT				
ARM1445		Lateral Support with extended arm, multi-adjustable pad	(NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£165
MANUAL / POWERED SEATING OPTIONS*				
ARM1091		Manually settable Seat Tilt 0° - 20°		◆
ARM1092		Powered "fixed pivot" Tilt 20° (136kg max user weight)	(NC 2201)	£400
ARM2055		Powered "fixed pivot" Tilt 20° & Sit Smarter (136kg max user weight)	(NC 2201, 2054)	£540
ARM1095		Powered "fixed pivot" Tilt 20°, heavy duty (180kg max user weight)	(NC 2201)	£400
ARM2056		Powered "fixed pivot" Tilt 20° & Sit Smarter, heavy duty (180kg max user weight)	(NC 2201, 2054)	£540
ARM1093		Powered Tilt (0° - 30°) with CoG shift (160kg max user weight)	(NC 0354, 2201)	£1,100
ARM2057		Powered Tilt (0° - 30°) with CoG shift & Sit Smarter (160kg max user weight)	(NC 2201)	£1,490
ARM2052		Powered Tilt (0°-50°) with CoG shift (136kg max user weight)	(NC 2201)	£1,450
ARM2059		Powered Tilt (0° - 50°) with CoG shift & Sit Smarter (136kg max user weight)	(NC 2201, 2171)	£2,000
ARM1094		Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (160kg max user weight)	(NC 0354, 2201)	£1,900
ARM2058		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift & Sit Smarter (160kg max user weight)	(NC 2201, 2220)	£2,570
ARM2053		Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (136kg max user weight)	(NC 2201)	£2,500
ARM2060		Seat Lifter with Powered Tilt (300mm, 45°) with CoG & Sit Smarter (136kg max user weight)	(NC 2201, 2220, 2171)	£3,370

SEAT TO GROUND HEIGHT*				
ARM0354		450mm (17.75") (from floor to front of seat plate at 0° tilt)	(NC 1094, 2052, 2053)	○
ARM0364		475mm (18.75") (from floor to front of seat plate at 0° tilt)		○
SEAT TILT OPTION				
ARM2054		5° Fixed Anterior Tilt, Assist mount. Will reduce max tilt angle by 5°	(NC 1092, 1095)	£320
BACKREST RECLINER*				
ARM0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable		◆
ARM0977		Powered Recliner (0° - 30°) Foldable	(NC 2201)	£500
ARM0979		Powered Recliner 30° & Sit Smarter		£670
LEG REST TYPE*				
ARM0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0517, 0512, 0504, 0513, 0501)	£370
ARM0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 2201, 0504, 0512, 0504, 0513, 0501)	£1,620
ARM0600		Centre Mounted: LNX Powered Centre Mounted Legrest & Sit Smarter, length compensating. Length 340 - 410 mm	(NC 2201, 0512, 0504, 0513, 0501)	£2,190
UNZ0592		Centre Mounted: LNX Powered Centre Mounted Legrest with Telescoping Foot Platform (lowers to the ground). Length: 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting.	(NC 2201, 0504, 0512, 0504, 0513, 0501)	POA
ARM0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm	(NC 2365, 0559, 0516, 0517, 0501)	◆
ARM0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	(NC 0559, 0516, 0517, 0554, 0555, 0501)	£120
ARM0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£235
ARM0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£660
ARM0530		Side Mounted: ADE pair & Sit Smarter. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£900
ARM0501		Without Leg Rests	(NC 0504, 0513, 2365, 0559, 0516, 0517, 0559, 0516, 0517)	○
FOOTPLATE TYPE*				
ARM2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	(NC 0503, 0533, 0531)	£105
ARM0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	(NC 0503, 0533, 0531)	○
ARM0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	(NC 0503, 0533, 0531)	£135
ARM0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	(NC 0526, 0503, 0533, 0531)	£270
ARM0512		Fixed Footplate pair	(NC 0533, 0531)	◆
ARM0504		Angle Adjustable Footplate pair	(NC 0526)	£75
ARM0513		Angle and depth adjustable Footplate pair		£120
ARM0501		Without footplates for option 0501 without Leg Rests	(NC 0503, 0533, 0531)	○
AMPUTEE SUPPORT				
ARM0554		Amputee Support RH	(NC 0556, 0589, 0592)	£195
ARM0555		Amputee Support LH	(NC 0556, 0589, 0592)	£195

LiNX CONTROLS*				
ARM2201		DLX-REM110 (Drive Only)	(NC 1437, 1092, 1095, 1093, 1094, 0977, 0589, 0531, 2231, 2232, 2219, 2245, 2270, 1908, 2209, 2210)	◆
ARM2203		DLX-REM211 (Drive & Actuators)	(NC 1437, 2231, 2232, 2270, 1908, 2209, 2210)	○
ARM2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2231, 2232, 2270, 2209, 2210)	£200
ARM2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2209, 2210)	£425
ARM2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1208, 1707, 1708, 1709, 1710, 2213, 2214, 2227, 2228, 2220, 2270)	£765
REM400 TOGGLE OPTION				
ARM2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206)	£315
CONTROL BRACKETRY*				
ARM1207		Standard Bracket RH		○
ARM1208		Standard Bracket LH		○
ARM1707		Swing Away RH Bracket		£135
ARM1708		Swing Away LH Bracket		£135
ARM1709		Matrx Resolve Swing in/out Bracket RH	(NC 2213, 2213, 1404)	£245
ARM1710		Matrx Resolve Swing in/out Bracket LH	(NC 2214, 2214, 1405)	£245
ARM2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 1405, 1415)	£495
ARM2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 1405, 1339)	£495
ARM2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 1415)	○
ARM2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 1339)	○
ARM2227		Bracket RH mounted for Attendant	(NC 2220, 2108, 2128)	£65
ARM2228		Bracket LH mounted for Attendant	(NC 2220, 2108, 2128)	£65
ARZ2239		Individual mounting option per request (Special). Requires additional information/specification.		POA
DUAL CONTROLS				
ARM2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 2108, 2205, 2206, 1221, 1222, 9260)	£395
ARM2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1222, 2108, 2201, 9260)	£395
ARM2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2201, 9260)	£745
ARM2420		mo-vis Scoot Control w/ dedicated bracket included.	(NC 2108, 1221, 1222)	£3,420
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
ARM1207		Dual Control Standard Bracket RH		○
ARM1208		Dual Control Standard Bracket LH		○
ARM1221		IDC Dual Control Detachable Bracket RH	(NC 2220, 2227)	○
ARM1222		IDC Dual Control Detachable Bracket LH	(NC 2220, 2228)	○
ARM9260		Bracket for mo-vis Scoot	(NC 2219, 2245)	○
COMPACT REMOTES				
ARM2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2220)	£825
ARM2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2220)	£1,075
ARM3000		mo-vis Micro LiNX	(NC 2201, 2203, 2204)	£2,960
ARM3001		mo-vis Multi LiNX	(NC 4401, 4402, 2201, 2203, 2204)	£2,640
ARM3110		mo-vis All Round Low Force Omni	(NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
ARM3111		mo-vis All Round Omni	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
ARM3112		mo-vis Heavy Duty Omni	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4300, 4301, 4302, 4400, 4401, 4402, 2201, 2203, 2204)	£3,000

COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
ARM2207		RH Mount	(NC 2201, 2203, 2204)	○
ARM2208		LH Mount	(NC 2201, 2203, 2204)	○
ARM2213		Centre Mount Swing Away to RH (Nucleus)	(NC 1405, 1415, 2201, 2203, 2204)	£495
ARM2214		Centre Mount Swing Away to LH (Nucleus)	(NC 1405, 1339, 2210, 2203, 2204)	£495
ARZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
ARM3300		Maxx Resolve RH	(NC 4301, 4302)	£270
ARM3301		Maxx Resolve LH	(NC 4301, 4302)	£270
ARM3310		BodyPoint Trilock Midline RH	(NC 4301, 4302)	£480
ARM3311		BodyPoint Trilock Midline LH	(NC 4301, 4302)	£480
ARM3320		Tray Mounted	(NC 4301, 4302)	○
ARZ3330		Individual mounting option per request (Special). Requires additional information/specification.		POA
CHIN CONTROLS				
ARM3000		mo-vis Micro LiNX	(NC 2201, 2203, 2204)	£2,960
ARM3001		mo-vis Multi LiNX	(NC 4401, 4402, 2201, 2203, 2204)	£2,640
ARM3110		mo-vis All Round Low Force Omni (consider option 3810 for switch access)	(NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
ARM3111		mo-vis All Round Omni (consider option 3810 for switch access)	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
ARM3400		mo-vis Chin Control Harness (One Size)	(NC 4301, 4302)	£1,000
ARM3410		mo-vis Manual Lateral Swing Away RH Mounted	(NC 4301, 4302)	£1,520
ARM3411		mo-vis Manual Lateral Swing Away LH Mounted	(NC 4301, 4302)	£1,520
ARM3412		mo-vis Manual Vertical Swing Away RH Mounted	(NC 4301, 4302)	£1,520
ARM3413		mo-vis Manual Vertical Swing Away LH Mounted	(NC 4301, 4302)	£1,520
ARM3420		mo-vis Powered Lateral Swing Away RH Mounted	(NC 4301, 4302)	£2,400
ARM3421		mo-vis Powered Lateral Swing Away LH Mounted	(NC 4301, 4302)	£2,400
ARM3422		mo-vis Powered Vertical Swing Away RH Mounted	(NC 4301, 4302)	£2,400
ARM3423		mo-vis Powered Vertical Swing Away LH Mounted	(NC 4301, 4302)	£2,400
ARZ3330		Individual mounting option per request (Special). Requires additional information/specification.		POA
HEAD CONTROLS				
ARM3200		mo-vis Proportional Head Control Omni	(NC 2201, 2203, 2204)	£4,500
ARM3221		mo-vis Switched Head Control Omni	(NC 2201, 2203, 2204)	£4,500
ARM3222		mo-vis Switched Head Control w/ Sip & Puff Omni	(NC 2201, 2203, 2204)	£4,750
FOOT CONTROL				
ARM3010		mo-vis Foot Control LiNX RH Side	(NC 2203, 2204, 2220)	£3,300
ARM3011		mo-vis Foot Control LiNX LH Side	(NC 2203, 2204, 2220)	£3,300
LiNX 10-WAY SWITCH				
ARM2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		£500
LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)				
ARM2217		Bracket RH, mounted for user	(NC 2213, 1415)	○
ARM2218		Bracket LH, mounted for user	(NC 2214, 1339)	○
ARM2227		Bracket RH, mounted for attendant	(NC 2128, 1221)	£65
ARM2228		Bracket LH, mounted for attendant	(NC 2128, 1222)	£65
LiNX ACCESSORY CABLE				
ARM2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)		£100
LiNX USB CHARGER PORT				
ARM2253		USB Charger Port DLX-USB02-A	(NC 2701, 2702, 2703, 2704, 2705, 2706)	£160

LINX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
ARM2207		RH Mount	○
ARM2208		LH Mount	○
FURTHER LINX / MO-VIS COMPONENTS / ACCESSORIES			
ARM2302		Not Necessary (NC 2254, 1711, 1712, 1713, 1714, 1715, 1716)	◆
ARM2300		To be enclosed only. Standard 300 mm Cable Length for mo-vis Switches	○
ARM2305		To be enclosed only. Extended 1500 mm Cable Length for mo-vis Switches	○
ARZ2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
ARM1685		Buddy Button. Red / Yellow / Green / Blue / Black provided. Interchangeable colours.	£90
EGG SWITCH			
ARM2254		Egg Switch, black	£90
PIKO BUTTONS			
ARM1711		Red Piko Button	£160
ARM1712		Yellow Piko Button	£160
ARM1713		Green Piko Button	£160
ARM1714		Orange Piko Button	£160
ARM1715		Blue Piko Button	£160
ARM1716		Black Piko Button	£160
TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle - direct control for Tilt only	£270
ARZ2262		Stereo Button - direct control for Tilt only	£255
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	£330
ARZ2264		Four-Way Button - direct control for Tilt & Recline	£355
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
LINX INPUT MODULE			
UNZ2273		DLX-IN200 Module (Jack port for external power button) (NC 2220)	£270
UNZ9249		DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	£785
MO-VIS TWISTER BASIC			
ARM3500		Black	£140
ARM3501		Red	£140
ARM3502		Yellow	£140
ARM3503		Blue	£140
ARM3504		Green	£140
MO-VIS TWISTER BASIC MOUNTING (Only applicable for mo-vis Twister Basic)			
ARM3520		Mounting Set	£50
ARM3525		Table Mounting Set	£50

MO-VIS TWISTER			
ARM3600		Black	£160
ARM3601		Red	£160
ARM3602		Yellow	£160
ARM3603		Blue	£160
ARM3604		Green	£160
MO-VIS TWISTER ON BENDED TUBE INC. D6 CLAMP			
ARM3700		Black	£220
ARM3701		Red	£220
ARM3702		Yellow	£220
ARM3703		Blue	£220
ARM3704		Green	£220
MO-VIS SATELLITE TWISTER (Only applicable for mo-vis All Round Remotes)			
ARM3800		Black	£200
ARM3801		Red	£200
ARM3802		Yellow	£200
ARM3803		Blue	£200
ARM3804		Green	£200
ARM3810		All Round Chin Control Set (includes All Round Joystick Ball and 2x Satellite Switches)	£420
MO-VIS TWISTER GOOSENECK			
ARM3900		Black	£260
ARM3901		Red	£260
ARM3902		Green	£260
MO-VIS Q2M TWIST & SWING			
ARM4000		Black	£500
ARM4001		Red	£500
ARM4002		Yellow	£500
ARM4003		Blue	£500
ARM4004		Green	£500
MO-VIS TWISTER PRO			
ARM4100		Black	£160
ARM4101		Red	£160
ARM4102		Green	£160
MO-VIS MULTI SWITCH			
ARM4200		Multi Switch	£620
ARM4201		Multi Switch Proximity Sensors	£860
ARM4202		Multi Switch Proximity Set	£840
MO-VIS HANDWARMER			
ARM4300		Hand Warmer	£800
ARM4301		Hand Warmer 2 Dome Hand Hood	£400
ARM4302		Hand Warmer 1 Dome Hand Hood	£300

MO-VIS JOYSTICK KNOBS			
ARM4400		Micro Joystick Ball	£60
ARM4401		Micro Joystick Cup 10.1mm	£60
ARM4402		Micro Joystick Cup 15.5mm	£60
ARM4403		Multi Joystick Ball 4mm	£60
ARM4404		Multi Joystick Cup 24mm	£60
ARM4405		All Round Joystick Ball 6.35mm	£60
ARM4406		All Round Joystick Cup 37mm	£60
BODYPOINT JOYSTICK KNOBS			
ARM2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£75
ARM2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£75
ARM2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
ARM2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
SHROUD COLOUR*			
ARM2640		Obsidian Black Frame Shrouds	○
ARM2641		Ruby Red Frame Shrouds	○
ARM2642		Purple Passion Frame Shrouds	○
ARM2643		Peacock Blue Frame Shrouds	○
ARM2644		Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARM2650		Obsidian Black Rim Inserts	○
ARM2651		Ruby Red Rim Inserts	○
ARM2652		Purple Passion Rim Inserts	○
ARM2653		Peacock Blue Rim Inserts	○
ARM2654		Arctic Blue Rim Inserts	○
FORKS*			
ARM2631		Single sided fork	£300
ARM2632		Double sided fork	◆
CENTRE WHEELS (DRIVE WHEELS)*			
ARM0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
ARM0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	○
ARM0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected, gel-sealed	○
MOTOR TYPE*			
ARM1294		3 km/h (2mph)	○
ARM1282		6 km/h (4mph)	◆
ARM1908		10 km/h (6mph) (includes light system and indicators (NC 2203, 1437, 1254)	£515
BATTERY*			
ARM1257		50 Ah Gel	◆
G TRAC			
ARM1215		G-Trac (G-trac ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£400
CHARGER*			
ARM1343		OptiCharge Charger 8A (200 - 250V)	◆
ARM1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£155

LIGHTS*				
ARM1437		LED Light System with Indicators (30 LUX)	(NC 2203, 1908)	£450
ARX1006		Lights Included with Motor Package 10 & 12kmh	(NC 1282)	○
ARM1254		Without Light system		◆
TRANSPORTATION KIT*				
ARM1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.		£580
ARM1493		4-Point tie-down (For wheelchair transportation without user)		○
ARM1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.		◆
FIXATION BELT*				
ARM1329		Lap Belt with Auto Lock	(NC 2702, 2703, 2704, 2705, 2706)	◆
ARM2113		Lap belt with auto lock (flexible mount in side rails)	(NC 2701, 2702, 2703, 2704, 2705, 2706)	○
ARM1332		H-Belt	(NC 2108)	£75
ARM2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£95
ARM2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£95
ARM2703		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£90
ARM2704		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£90
ARM2705		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£90
ARM2706		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£90
REAR VIEW MIRROR				
ARM1415		Rear view Mirror RH	(NC 2213, 2217)	£45
ARM1339		Rear view Mirror LH	(NC 2214, 2218)	£45
BAG HOOKS				
ARM2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	(NC 2173, 2227, 2228)	£80
CANE HOLDER				
ARM1312		Cane holder RH		£45
ARM1363		Cane holder LH		£45
TRAY				
ARM1404		Tray, Swing Away RH	(NC 1707, 2209, 2213, 2214, 2207)	£200
ARM1405		Tray, Swing Away LH	(NC 1708, 2210, 2213, 2214, 2208)	£200
TOTAL PRICE				

Non Conforming Options:

ARM2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
ARM2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
ARM2123	480 mm Seat Plate
ARM2124	530 mm Seat Plate
ARM2132	410 mm
ARM2134	460 mm
ARM2138	560 mm
ARM2142	Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg)
ARM2143	Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm, up to 125kg)
ARM2147	Matrx Libra (90mm - standard up to 136kg or HD up to 227kg). Extended lead time.
ARM2144	Matrx Vi Cushion (Mesh fabric)
ARM2148	Standard Cushion (90mm thick) - For 550mm Seat Plate Only
ARM2140	Without seat cushion
ARM1980	Tapered Wedges (Pair)
ARM1981	Abductor
ARM1982	Lateral Wedge
ARM1983	Anterior Seat Wedges
ARM1980	Tapered Wedges (Pair)
ARM1981	Abductor
ARM1982	Lateral Wedge
ARM1983	Anterior Seat Wedges
ARM1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
ARM2163	Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.
ARM2108	Backposts WITHOUT push handles. Push handles will be cut down on the vertical section.
ARM2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
ARM2179	Flex3 Tension Adjustable (540 mm or 480 mm)
ARM2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ARM2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ARM2400	Support+ Standard Backrest (Spacer fabric)
ARM2401	Support+ Active Strapped Backrest (Spacer fabric)
ARM2402	Support+ Passive Strapped Backrest (Spacer fabric)
ARM2403	Support+ Extra Deep Backrest (Spacer fabric)
ARM2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARM2168	Matrx MAC Back - Foam Cushion [Mesh fabric]
ARM2169	Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]

ARM2170	Without Backrest Upholstery (Back Posts/Back Frame only)
ARM4520	Tall - 570 mm (Measured from Seat Plate)
ARM4521	Short - 500 mm (Measured from Seat Plate)
ARM2380	405 mm (16 inch) MAC BH
ARM2381	455 mm (18 inch) MAC BH
ARM2382	505 mm (20 inch) MAC BH
ARM2383	560 mm (22 inch) MAC BH
ARM2385	90 mm (3.5 inch) MAC Contour Depth LH
ARM2386	150 mm (6 inch) MAC Contour Depth LH
ARM2387	230 mm (9 inch) MAC Contour Depth LH
ARM2389	90 mm (3.5 inch) MAC Contour Depth RH
ARM2390	150 mm (6 inch) MAC Contour Depth RH
ARM2391	230 mm (9 inch) MAC Contour Depth RH
ARM2393	Backrest -25 mm (-1 inch) to seat width
ARM2394	Backrest matches seat width
ARM2395	Backrest +25 mm (+1 inch) to seat width
ARM2151	Black Nylon Fabric TR18
ARM2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
ARM2155	Black Spacer Fabric TR36
ARM2181	T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm
ARM2183	Inrigged Armrests. Reduces SW by 25 mm each side
ARM2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).
ARM2192	Following Armrest (height adjustable from 230mm - 300mm).
ARM2194	Following Armrest (height adjustable from 300mm - 360mm).
ARM2195	Short Armpad (325mm)
ARM2196	Long Armpad (375mm)
ARM2197	Long Ergonomic Armpad (397mm)
ARM2410	Support+ Standard Armpad Short (80 mm x 300 mm)
ARM2411	Support+ Standard Armpad Long (80 mm x 400 mm)
ARM2412	Support+ P-Shaped Armpad (120 mm x 400 mm)
ARM2199	Soft Armpad (REA Azalea Wide Pad)
ARM2184	Hip Support LH, pad size 140 x 100 mm
ARM2186	Hip Support LH, pad size 230 x 100 mm
ARM2188	Hip Support LH, pad size 330 x 100 mm
ARM2185	Hip Support RH, pad size 140 x 100 mm
ARM2187	Hip Support RH, pad size 230 x 100 mm
ARM2189	Hip Support RH, pad size 330 x 100 mm
ARM1350	Neckrest depth adjustable
ARM1432	Headrest sideways adjustable
ARM1427	Headrest with cheek support
ARM2328	Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).
ARM2329	Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem (_2332).
ARM2330	Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).

ARM2331	Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem (2332).
ARM2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
ARM2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
ARM2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
ARM1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
ARM1495	Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts
ARM4510	Support+ Wipedown Cushion Cover (PU Coated Type)
ARM4530	Support+ Wipedown Backrest Cover (PU Coated Type)
ARM4540	Support+ Wipedown Armpad Covers (PU Coated Type)
ARM4550	Support+ Wipedown Headrest Cover (PU Coated Type)
ARM1445	Lateral Support with extended arm, multi-adjustable pad
ARM1092	Powered "fixed pivot" Tilt 20° (136kg max user weight)
ARM2055	Powered "fixed pivot" Tilt 20° & Sit Smarter (136kg max user weight)
ARM1095	Powered "fixed pivot" Tilt 20°, heavy duty (180kg max user weight)
ARM1093	Powered Tilt (0° - 30°) with CoG shift (160kg max user weight)
ARM2052	Powered Tilt (0°-50°) with CoG shift (136kg max user weight)
ARM1094	Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (160kg max user weight)
ARM2053	Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (136kg max user weight)
ARM0354	450mm (17.75") (from floor to front of seat plate at 0° tilt)
ARM2054	5° Fixed Anterior Tilt, Assist mount. Will reduce max tilt angle by 5°
ARM0977	Powered Recliner (0° - 30°) Foldable
ARM0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
ARM0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
UNZ0592	Centre Mounted: LNX Powered Centre Mounted Legrest with Telescoping Foot Platform (lowers to the ground). Length: 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting.
ARM0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm
ARM0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
ARM0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
ARM0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
ARM0501	Without Leg Rests
ARM2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
ARM0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
ARM0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)
ARM0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
ARM0512	Fixed Footplate pair
ARM0504	Angle Adjustable Footplate pair
ARM0513	Angle and depth adjustable Footplate pair
ARM0501	Without footplates for option 0501 without Leg Rests
ARM0554	Amputee Support RH
ARM0555	Amputee Support LH
ARM2201	DLX-REM110 (Drive Only)

ARM2203	DLX-REM211 (Drive & Actuators)
ARM2204	DLX-REM216 (Drive, Actuators & Lights)
ARM2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
ARM2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
ARM2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
ARM1208	Standard Bracket LH
ARM1707	Swing Away RH Bracket
ARM1708	Swing Away LH Bracket
ARM1709	Matrx Resolve Swing in/out Bracket RH
ARM1710	Matrx Resolve Swing in/out Bracket LH
ARM2213	Centre Mount Swing Away to RH Nucleus Bracket
ARM2214	Centre Mount Swing Away to LH Nucleus Bracket
ARM2209	REM500 Swing Away RH Multi-Adjustable Bracket
ARM2210	REM500 Swing Away LH Multi-Adjustable Bracket
ARM2227	Bracket RH mounted for Attendant
ARM2228	Bracket LH mounted for Attendant
ARM2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
ARM2219	LiNX DLX-ACU200 (Compact Dual Control)
ARM2245	LiNX DLX-IDC (Intuitive Dual Control)
ARM1208	Dual Control Standard Bracket LH
ARM1221	IDC Dual Control Detachable Bracket RH
ARM1222	IDC Dual Control Detachable Bracket LH
ARM9260	Bracket for mo-vis Scoot
ARM2231	LiNX DLX-CR400 (Compact Remote)
ARM2232	LiNX DLX-CR400LF (Compact Remote Low Force)
ARM2207	RH Mount
ARM2208	LH Mount
ARM2213	Centre Mount Swing Away to RH (Nucleus)
ARM2214	Centre Mount Swing Away to LH (Nucleus)
ARM3410	mo-vis Manual Lateral Swing Away RH Mounted
ARM3411	mo-vis Manual Lateral Swing Away LH Mounted
ARM3412	mo-vis Manual Vertical Swing Away RH Mounted
ARM3413	mo-vis Manual Vertical Swing Away LH Mounted
ARM3420	mo-vis Powered Lateral Swing Away RH Mounted
ARM3421	mo-vis Powered Lateral Swing Away LH Mounted
ARM3422	mo-vis Powered Vertical Swing Away RH Mounted
ARM3423	mo-vis Powered Vertical Swing Away LH Mounted
ARM3200	mo-vis Proportional Head Control Omni
ARM3221	mo-vis Switched Head Control Omni
ARM3222	mo-vis Switched Head Control w/ Sip & Puff Omni
ARM2217	Bracket RH, mounted for user
ARM2218	Bracket LH, mounted for user
ARM2227	Bracket RH, mounted for attendant
ARM2228	Bracket LH, mounted for attendant

ARM2253	USB Charger Port DLX-USB02-A
ARM2207	RH Mount
ARM2208	LH Mount
ARM2254	Egg Switch, black
ARM1711	Red Piko Button
ARM1712	Yellow Piko Button
ARM1713	Green Piko Button
ARM1714	Orange Piko Button
ARM1715	Blue Piko Button
ARM1716	Black Piko Button
ARZ2207	Toggle / Button RH Mount
ARZ2208	Toggle / Button LH Mount
ARM4300	Hand Warmer
ARM4301	Hand Warmer 2 Dome Hand Hood
ARM4302	Hand Warmer 1 Dome Hand Hood
ARM4400	Micro Joystick Ball
ARM4401	Micro Joystick Cup 10.1mm
ARM4402	Micro Joystick Cup 15.5mm
ARM1282	6 km/h (4mph)
ARM1908	10 km/h (6mph) (includes light system and indicators)
ARM1437	LED Light System with Indicators (30 LUX)
ARM1254	Without Light system
ARM1329	Lap Belt with Auto Lock
ARM2113	Lap belt with auto lock (flexible mount in side rails)
ARM1332	H-Belt
ARM2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARM2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARM2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ARM2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ARM2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ARM2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ARM1415	Rear view Mirror RH
ARM1339	Rear view Mirror LH
ARM2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
ARM1404	Tray, Swing Away RH
ARM1405	Tray, Swing Away LH



For a more complex customer requirement,
please email: salesUK@dhcare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):