



**Yes, you can:**

# Powerchairs Ireland Prescription Form

Please refer to [www.invacare.ie](http://www.invacare.ie) for Local Distributor

NOVEMBER 2023  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## TDX SP2 NB

### Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- Compact, manoeuvrable, CWD TDX SP2 Narrow Base
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable backrest Recline 0° - 30°, in 3.75° steps
- Standard 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 355mm x 76mm (14" x 3") Black tyres with Gel inserts (Solid)
- 6 km/h (4mph) 4 Pole motors and 50Ah batteries
- 8A (200 - 250V) charger
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: €7,624**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
ARM2100	<input type="checkbox"/> TDX SP2 NARROW BASE	€7,624
<b>USER WEIGHT*</b>		
ARM1920	<input type="checkbox"/> User weight up to 125kg	○
ARM1921	<input type="checkbox"/> User weight 126kg to 136kg	○
ARM1922	<input type="checkbox"/> User weight 137kg to 150kg (NC 1092)	○
<b>BASE MODEL*</b>		
ARM0000	<input type="checkbox"/> TDX SP2 Narrow Base	◆
ARZ0000	<input type="checkbox"/> Special TDX SP2 Narrow Base	POA

SEATING TYPE*			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	◆
SEAT FRAME*			
ARM2111	<input type="checkbox"/>	Telescopic Seat Frame	◆
SEAT WIDTH*			
ARZ1001	<input type="checkbox"/>	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	(NC 2181, 2193, 2192, 2194) ○
ARM2121	<input type="checkbox"/>	380mm (15")	(NC 2143, 2701, 2703, 2705) ○
ARM2122	<input type="checkbox"/>	430mm (17")	(NC 2703, 2705) ○
ARM2123	<input type="checkbox"/>	480mm (19")	(NC 2702, 2704, 2706) ○
ARM2124	<input type="checkbox"/>	530mm (21")	(NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706) ○
SEAT DEPTH RAIL LENGTH*			
ARM2129	<input type="checkbox"/>	Standard 410mm - 510mm	(NC 2138) ○
ARM2190	<input type="checkbox"/>	Extended 510mm - 560mm	(NC 2132, 2134) ○
SEAT DEPTH*			
ARM2132	<input type="checkbox"/>	410mm (16")	(NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706) ○
ARM2134	<input type="checkbox"/>	460mm (18")	○
ARM2136	<input type="checkbox"/>	510mm (20")	○
ARM2138	<input type="checkbox"/>	560mm (22")	(NC 2143) ○
SEAT OPTION			
ARZSWO	<input type="checkbox"/>	H Frame Interface Kit	€507
SEAT CUSHION*			
ARM2142	<input type="checkbox"/>	Matrx Flo-Tech Lite	€126
ARM2143	<input type="checkbox"/>	Matrx Flo-tech Contour Visco (430mm x 430mm only - up to 125kg)	(NC 2136, 2121, 2124, 2138) €126
ARM2147	<input type="checkbox"/>	Matrx Libra	€469
ARM2140	<input type="checkbox"/>	Without seat cushion	○
BACKREST SUPPORT*			
ARM2163	<input type="checkbox"/>	Short 480mm	(NC 2172, 2173) ○
ARM2162	<input type="checkbox"/>	Tall 540mm	◆
PUSH HANDLES FOR BACKPOSTS*			
ARM2107	<input type="checkbox"/>	WITH push handles	◆
ARM2108	<input type="checkbox"/>	WITHOUT push handles	(NC 1332, 2219, 2220, 2227, 2228, 2245) €48
BACKREST UPHOLSTERY / CUSHION*			
ARM2171	<input type="checkbox"/>	Tension Adjustable Backrest	(NC 1445, 1496, 2155, 2230, 2247, 1288, 1289) ◆
ARM2179	<input type="checkbox"/>	Flex3 Tension Adjustable	(NC 2151, 2192, 2193, 2194, 2199, 2229, 2155, 2230, 2247, 2128) €182
ARM2172	<input type="checkbox"/>	Flex3 Backrest plates without cushion	(NC 2163, 2128) €252
ARM2173	<input type="checkbox"/>	Flex3 Backrest plates with Laguna cushion	(NC 2163, 2128) €411
ARM2174	<input type="checkbox"/>	Matrx PB 400mm	(NC 1496, 2124) €572
ARM2175	<input type="checkbox"/>	Matrx PB 500mm	(NC 1496, 2124) €572
ARM2176	<input type="checkbox"/>	Matrx Elite	(NC 1496, 2124) €572
ARM2177	<input type="checkbox"/>	Matrx Elite TR	(NC 1496, 2124) €572
ARM2178	<input type="checkbox"/>	Matrx Elite Deep	(NC 1496, 2124) €741
ARM2170	<input type="checkbox"/>	Without Backrest	(NC 1496, 2230, 2247) ○
UPHOLSTERY TYPE*			
ARM2151	<input type="checkbox"/>	Black fabric	◆
ARM2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	(NC 2171) ○
ARM2155	<input type="checkbox"/>	Black Spacer fabric	(NC 2171, 2179) ○
ARZ1002	<input type="checkbox"/>	Without Backrest Upholstery	(NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178) ○

ARMREST*			
ARM2191		T-Armrest standard	◆
ARM2181		T-Armrest extended height	€48
ARM2183		In-rigged armrests (NC 2185, 2186, 2187, 2188, 2189)	€28
ARM2193		Flip-up armrests (NC 1445, 2179)	€41
ARM2192		Following 230mm - 300mm (NC 1445, 2179)	€200
ARM2194		Following 300mm - 360mm (NC 1445, 2179)	€200
ARMPAD*			
ARM2195		Short 325mm Armpad	◆
ARM2196		Long 375mm Armpad	○
ARM2197		Long ergonomic 397mm Armpad (NC 2230, 2247)	○
ARM2199		Soft Armpad (NC 2179)	€103
HIP SUPPORTS LH			
ARM2184		Small 140mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARM2186		Medium 230mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARM2188		Large 330mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	€272
HIP SUPPORTS RH			
ARM2185		Small 140mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARM2187		Medium 230mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	€272
ARM2189		Large 330mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	€272
MOUNTING HARDWARE FOR HIP SUPPORTS			
ARM2165		Extended height (+2")	○
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARM1350		Neckrest height and depth adjustable (NC 2247)	€38
ARM1432		Headrest, height, depth and sideways adjustable (NC 2247)	€228
ARM1427		Headrest with cheek support, height and depth adjustable (NC 2247)	€239
ARM2324		Elan standard pad 255mm (10") (NC 2247)	€406
ARM2325		Elan occipital pad 230mm x 180mm (9" X 7") (NC 2247)	€406
ARM2327		Elan 4 point headrest pad 280mm x 255mm (11" X 10") (NC 2247)	€548
ARM1496		All bracketry as per option 1432 (NC 2171, 2174, 2175, 2176, 2177, 2178)	○
ARM1495		Push bar 22mm	○
LATERAL SUPPORT			
ARM1445		Multi-adjustable pad (NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	€170
SEAT LIFTER / SEAT TILT*			
ARM1091		Manually settable Seat Tilt (0° - 20°) Seat Height 450mm	○
ARM1092		Electric Tilt (0° - 20°) (NC 2201)	€372
ARM1095		Electric Tilt (0° - 20°) Max User weight 150kg (NC 2201)	€372
ARM1093		Electric Tilt (0° - 30°) with CoG shift Seat Height 460mm (NC 0354, 0364, 2201)	€374
ARM1094		Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 460mm - 490mm (NC 0354, 0364, 2201)	€1,419
SEAT HEIGHT SETTING*			
ARM0354		Seat height 450mm (17.75")	○
ARZ1003		Seat height set at 460mm (18") (from floor to front of seat plate at 0° tilt) For 30° tilt and Lifter (NC 1092, 1095)	○
ARM0364		Seat height 475mm (18.75")	○
SEAT TILT OPTION			
ARM2054		5° fixed anterior Tilt assist mount (NC 1092, 1095)	€476

<b>BACKREST RECLINER*</b>			
ARM0975		Manual settable Recliner	◆
ARM0977		Powered Recliner (NC 2201)	€338
<b>FOOTREST HANGER TYPE*</b>			
ARM0556		Centre mounted : Manual (NC 0517, 0512, 0504, 0513)	€321
ARM0589		Centre mounted : LNX Powered (NC 2201, 0504, 0512, 0504, 0513)	€1,914
ARM0526		Side mounted : Standard 80° hanger	◆
ARM0503		Side mounted : VARI F	○
ARM0533		Side mounted : VARI A	€133
ARM0531		Side mounted : ADE Powered (NC 2201)	€686
ARM0501		Without Hanger (NC 0504, 0513)	○
<b>AMPUTEE SUPPORT</b>			
ARM0555		LH Amputee Support	€202
ARM0554		RH Amputee Support	€202
<b>FOOTPLATE TYPE*</b>			
ARM2365		Angle Adjustable Footplates Small 127mm x 190mm (5" x 7.5") for Centre Mounted Hanger only (NC 0503, 0533, 0531, 0501)	€129
ARM0559		Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger only (NC 0503, 0533, 0531, 0501)	○
ARM0516		Angle Adjustable Footplates Large 150mm x 290mm (6" x 11.5") for Centre Mounted Hanger only (NC 0503, 0533, 0531, 0501)	€167
ARM0517		Angle Adjustable foot platform Medium 305mm x 275mm (12" x 10.75") for Centre Mounted Hanger only (NC 0526, 0503, 0533, 0531, 0501)	€376
ARM0512		Fixed Footplate (NC 0533, 0531)	◆
ARM0504		Angle adjustable Footplate (NC 0526)	€75
ARM0513		Angle and Depth adjustable	€122
ARZ1004		Without Footplates (NC 0589, 0526, 0503, 0533, 0531)	○
<b>CONTROLS*</b>			
ARM2201		REM110 (NC 1437, 1092, 1095, 1093, 1094, 0977, 0589, 0531, 2229, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2219, 2245, 2270, 1908, 2209, 2210)	◆
ARM2203		REM211 (NC 1437, 2229, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2270, 1908, 2209, 2210)	○
ARM2204		REM216 (NC 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2270, 2209, 2210)	€159
ARM2205		REM400 (NC 2209, 2210)	€473
ARM2206		REM500 (NC 1708, 2213, 2214, 2227, 2228, 2220, 2270)	€909
<b>CONTROL BRACKETRY*</b>			
ARM1207		Standard Bracket RH - Not Compatible with REM500	○
ARM1208		Standard Bracket LH - Not Compatible with REM500	○
ARM1707		Swing away RH Bracket	€151
ARM1708		Swing away LH Bracket	€151
ARM2213		Centre mount, swing away to RH Nucleus (NC 1405, 1415)	€739
ARM2214		Centre mount, swing away to LH Nucleus (NC 1405, 1339)	€739
ARM2227		RH Bracket mounted for Attendant (NC 2220, 2108, 1286, 2128)	€71
ARM2228		LH Bracket mounted for Attendant (NC 2220, 2108, 1286, 2128)	€71
ARM2209		REM500 Swing away RH multi-adjustable Bracket (NC 1415)	○
ARM2210		REM500 Swing away LH multi-adjustable Bracket (NC 1339)	○
ARZ2239		Individual mounting option per request (Special)	POA
<b>REM400 TOGGLE OPTION</b>			
ARM2270		Assignable toggle switch module (NC 2203, 2204, 2206)	€329

DUAL CONTROL				
ARM2220		Dual control REM050	(NC 2108, 2205, 2206, 1221, 1222, 2230, 2247, 2257, 2233, 2240)	€489
ARM2219		Compact (ACU200)	(NC 1222, 2108, 2201)	€489
ARM2245		IDC Dual control	(NC 2201)	€748
DUAL CONTROL BRACKETRY				
ARM1207		Standard Bracket RH		○
ARM1208		Standard Bracket LH		○
ARM1221		IDC Dual control detachable Bracket RH	(NC 2220, 2227)	○
ARM1222		IDC Dual control detachable Bracket LH	(NC 2220, 2228)	○
LINX COMPACT REMOTE & BRACKETS (Requires REM400 or REM500)				
ARM2231		Proportional Compact Remote	(NC 2203, 2204)	€845
ARM2232		Proportional Compact Remote low force	(NC 2203, 2204)	€1,199
ARM2207		RH mount Bracket	(NC 2201, 2203, 2204)	○
ARM2208		LH mount Bracket	(NC 2201, 2203, 2204)	○
ARM2213		Centre mount, swing away to RH Nucleus Bracket	(NC 1405, 1415, 2201, 2203, 2204)	€739
ARM2214		Centre mount, swing away to LH Nucleus Bracket	(NC 1405, 1339, 2210, 2203, 2204)	€739
ARZ2215		Individual mounting option per request (Special)		POA
CHIN CONTROL & BRACKETS (Requires REM400 or REM500)				
ARM2229		Proportional Chin joystick (ASL138)	(NC 2203, 2204, 2128, 2179)	€3,106
ARM1289		Swing away RH Manual Bracket	(NC 1432, 1427, 1496, 1495, 2171, 2201, 2203, 2204)	○
ARM1288		Swing away LH Manual Bracket	(NC 1432, 1427, 1496, 1495, 2171, 2201, 2203, 2204)	○
ARM1287		Swing away RH Powered Bracket	(NC 1404, 1415, 1312, 2201, 2203, 2204)	€2,500
ARM1286		Swing away LH Powered Bracket	(NC 1339, 1363, 2201, 2203, 2204)	€2,500
LINX HEAD CONTROL ATOM TECHNOLOGY (Requires REM400 or REM500)				
ARM2230		ATOM with ultrasonic proximity sensor	(NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2229, 2233, 2324, 2325, 2327, 1350, 1432, 1427, 2257, 2240)	€4,497
ARM2247		PROTON with ultrasonic proximity sensor	(NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2229, 2233, 2324, 2325, 2327, 1350, 1432, 1427, 2257, 2240)	€4,497
FOUR-SWITCH PROXIMITY ARRAY (with ATOM TECHNOLOGY)				
ARM2257		Four-switch Array (ASL106)	(NC 2203, 2204, 2220, 2229, 2230, 2247, 2233, 2240, 2128)	€8,205
ARM2275		Centre mount, swing-away to RH nucleus	(NC 2203, 2204, 1404, 1405, 1415)	€586
ARM2276		Centre mount, swing-away to LH nucleus	(NC 2203, 2204, 1404, 1405, 1339)	€586
ARZ2215		Individual mounting option per request (Special)		POA
WIRELESS ACCESSORIES FOR ATOM				
ARM2249		Wireless triple switch receiver (ASL557-3)		€422
ARM2250		Wireless mouse emulator (ASL558)		€1,013
LINX PROPORTIONAL REMOTE / INPUTS & BRACKETS (Requires REM400 or REM500)				
ARM2233		MEC proportional joystick	(NC 2203, 2204, 2220, 2229, 2230, 2247, 2257, 2240, 2128)	€4,230
ARM2240		Proportional heavy duty with mode switch (ASL133)	(NC 2203, 2204, 2220, 2229, 2230, 2247, 2257, 2233, 2128)	€2,683
ARM2235		MEC tray mount RH (ASL603)	(NC 1415, 2201, 2203, 2204)	○
ARM2236		MEC tray mount LH (ASL603)	(NC 1339, 2201, 2203, 2204)	○
ARM2213		Centre mount, swing away to RH nucleus	(NC 1405, 1415, 2201, 2203, 2204)	€739
ARM2214		Centre mount, swing away to LH nucleus	(NC 1405, 1339, 2201, 2203, 2204)	€739
ARM2241		Centre mount, swing away to RH nucleus	(NC 1405, 1415, 2201, 2203, 2204)	€739
ARM2242		Centre mount, swing away to LH nucleus	(NC 1405, 1339, 2201, 2203, 2204)	€739
ARZ2215		Individual mounting option per request (Special)		POA

<b>LINX 10-WAY SWITCH</b>			
ARM2251		DLX-SW10-A	€513
ARM2217		RH Bracket	(NC 2213, 2275, 2235, 2241, 1415) ○
ARM2218		LH Bracket	(NC 2214, 2276, 2236, 2242, 1339) ○
ARM2227		RH Bracket	(NC 2128, 1287, 1221) ○
ARM2228		LH Bracket	(NC 2128, 1286, 1222) ○
<b>LINX ACCESSORY</b>			
ARM2252		Loom providing jack sockets for drive and mode	€105
<b>LINX USB CHARGER PORT</b>			
ARM2253		USB charger port DLX-USB02-A	(NC 2701, 2702, 2703, 2704, 2705, 2706) €162
ARM2207		RH mount	○
ARM2208		LH mount	○
<b>REMOTE STOP SWITCH</b>			
ARM2258		Remote Stop Switch (ASL504)	(NC 2203, 2204) €1,019
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
ARM2300		To be enclosed only	○
ARZ2301		To be connected and programmed (special)	POA
<b>BUDDY BUTTONS</b>			
ARM1685		Red	€90
ARM1686		Yellow	€90
ARM1687		Green	€90
ARM1690		Blue	€90
<b>EGG SWITCH</b>			
ARM2254		Egg Switch, Black	€90
<b>PIKO BUTTONS</b>			
ARM1711		Red	€161
ARM1712		Yellow	€161
ARM1713		Green	€161
ARM1714		Orange	€161
ARM1715		Blue	€161
ARM1716		Black	€161
<b>TOGGLES / BUTTONS</b>			
ARZ2261		Stereo Toggle	€273
ARZ2262		Stereo Button	€263
ARZ2263		Four-way Toggle	€342
ARZ2264		Four-way Button	€365
ARZ2207		RH mount	○
ARZ2208		LH mount	○
ARZ2225		Without mounting	○
<b>BODYPOINT JOYSTICK KNOBS</b>			
ARM2721		3" U-shaped, joystick handle	€75
ARM2722		4" U-shaped, joystick handle	€75
ARM2723		Domed rubber, joystick handle	€75
ARM2724		Mushroom grooved, joystick handle	€75

<b>SHROUD COLOUR*</b>			
ARM0228	<input type="checkbox"/>	Wet Black	○
ARM2601	<input type="checkbox"/>	SpaceStation Silver	○
ARM2602	<input type="checkbox"/>	Island Blue	○
ARM2603	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARM2607	<input type="checkbox"/>	Rockstar Red	○
ARM2608	<input type="checkbox"/>	Lights Out Black - Matt	○
ARM2609	<input type="checkbox"/>	Invacare Blue	○
<b>RIM COLOUR*</b>			
ARM2611	<input type="checkbox"/>	SpaceStation Silver	○
ARM2612	<input type="checkbox"/>	Island Blue	○
ARM2613	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARM2617	<input type="checkbox"/>	Rockstar Red	○
ARM2618	<input type="checkbox"/>	Lights Out Black - Matt	○
ARM2619	<input type="checkbox"/>	Invacare Blue	○
ARM2620	<input type="checkbox"/>	Wet Black	○
ARM2621	<input type="checkbox"/>	Without rim insert	○
<b>FORKS*</b>			
ARM2631	<input type="checkbox"/>	Single sided	€358
ARM2632	<input type="checkbox"/>	Double sided	◆
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>			
ARM0704	<input type="checkbox"/>	Black Tyres, Pneumatic	◆
ARM0733	<input type="checkbox"/>	Black tyres with gel foam inserts (solid)	◆
<b>MOTOR TYPE*</b>			
ARM1294	<input type="checkbox"/>	3 km/h (2mph)	○
ARM1282	<input type="checkbox"/>	6 km/h (4mph)	◆
ARM1908	<input type="checkbox"/>	10 km/h (6mph)	(NC 2203, 1437, 1254) €527
<b>BATTERY*</b>			
ARM1257	<input type="checkbox"/>	50 Ah Gel	◆
<b>PERFORMANCE ACCESSORIES</b>			
ARM1215	<input type="checkbox"/>	G-Trac	€865
<b>CHARGER*</b>			
ARM1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V)	◆
ARM1313	<input type="checkbox"/>	OptiCharge 10A (200 - 250V)	€174
<b>LIGHTS*</b>			
ARM1437	<input type="checkbox"/>	LED headlights & indicators	(NC 2203, 1908) €394
ARZ1006	<input type="checkbox"/>	Lights Included with Motor Package 10 & 12kmh	(NC 1282) ○
ARM1254	<input type="checkbox"/>	Without light system	◆
<b>TRANSPORTATION KIT*</b>			
ARM1450	<input type="checkbox"/>	Preparation for DAHL docking station	€593
ARM1494	<input type="checkbox"/>	4-Point tie-down	◆

FIXATION BELT*			
ARM1329		Lap Belt fitted to base of backrest	(NC 2702, 2703, 2704, 2705, 2706) ◆
ARM2113		Lap Belt adjustable mounting on seat rail extrusion	(NC 2701, 2702, 2703, 2704, 2705, 2706) ○
ARM1332		H-Belt	(NC 2108) €75
ARM2701		Evoflex pelvic stabilizer, size L	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €117
ARM2702		Evoflex pelvic stabilizer, size M	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €117
ARM2703		2-Point hip belt centre pull, size L	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €108
ARM2704		2-Point hip belt centre pull, size M	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €108
ARM2705		2-Point hip belt dual pull, size L	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €108
ARM2706		2-Point Hip Belt Dual pull, size M	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €108
REAR VIEW MIRROR			
ARM1415		RH Mirror	(NC 2213, 1287, 2217, 2275, 2235, 2241) ○
ARM1339		LH Mirror	(NC 2214, 1286, 2218, 2276, 2236, 2242) ○
OTHER OPTIONS			
ARM2128		Bag hooks	(NC 2173, 2227, 2228, 2229, 2257, 2233, 2240) €95
ARM1312		RH Cane holder	€25
ARM1363		LH Cane holder	€25
ARM1404		Swing away RH Tray	(NC 1707, 2209, 2213, 2214, 2235, 2236, 2241, 2242, 1287, 2207, 2275, 2276) €373
ARM1405		Swing away LH Tray	(NC 1708, 2210, 2213, 2214, 2235, 2236, 2241, 2242, 1286, 2208, 2275, 2276) €373

<b>TOTAL PRICE</b>	
--------------------	--



**Non Conforming Options:**

ARM2121	380mm (15")
ARM2122	430mm (17")
ARM2123	480mm (19")
ARM2124	530mm (21")
ARM2132	410mm (16")
ARM2134	460mm (18")
ARM2136	510mm (20")
ARM2138	560mm (22")
ARM2143	Matrx Flo-tech Contour Visco (430mm x 430mm only - up to 125kg)
ARM2163	Short 480mm
ARM2108	WITHOUT push handles
ARM2171	Tension Adjustable Backrest
ARM2179	Flex3 Tension Adjustable
ARM2172	Flex3 Backrest plates without cushion
ARM2173	Flex3 Backrest plates with Laguna cushion
ARM2174	Matrx PB 400mm
ARM2175	Matrx PB 500mm
ARM2176	Matrx Elite
ARM2177	Matrx Elite TR
ARM2178	Matrx Elite Deep
ARM2151	Black fabric
ARM2155	Black Spacer fabric
ARM2181	T-Armrest extended height
ARM2183	In-rigged armrests
ARM2193	Flip-up armrests
ARM2192	Following 230mm - 300mm
ARM2194	Following 300mm - 360mm
ARM2196	Long 375mm Armpad
ARM2197	Long ergonomic 397mm Armpad
ARM2199	Soft Armpad
ARM2184	Small 140mm x 100mm
ARM2186	Medium 230mm x 100mm
ARM2188	Large 330mm x 100mm
ARM2185	Small 140mm x 100mm
ARM2187	Medium 230mm x 100mm
ARM2189	Large 330mm x 100mm
ARM1350	Neckrest height and depth adjustable
ARM1432	Headrest, height, depth and sideways adjustable
ARM1427	Headrest with cheek support, height and depth adjustable
ARM2324	Elan standard pad 255mm (10")
ARM2325	Elan occipital pad 230mm x 180mm (9" X 7")
ARM2327	Elan 4 point headrest pad 280mm x 255mm (11" X 10")
ARM1496	All bracketry as per option 1432

ARM1495	Push bar 22mm
ARM1445	Multi-adjustable pad
ARM1092	Electric Tilt (0° - 20°)
ARM1095	Electric Tilt (0° - 20°) Max User weight 150kg
ARM1093	Electric Tilt (0° - 30°) with CoG shift Seat Height 460mm
ARM1094	Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 460mm - 490mm
ARM0354	Seat height 450mm (17.75")
ARM0364	Seat height 475mm (18.75")
ARM0977	Powered Recliner
ARM0589	Centre mounted : LNX Powered
ARM0526	Side mounted : Standard 80° hanger
ARM0503	Side mounted : VARI F
ARM0533	Side mounted : VARI A
ARM0531	Side mounted : ADE Powered
ARM0501	Without Hanger
ARM0517	Angle Adjustable foot platform Medium 305mm x 275mm (12" x 10.75") for Centre Mounted Hanger only
ARM0512	Fixed Footplate
ARM0504	Angle adjustable Footplate
ARM0513	Angle and Depth adjustable
ARM2201	REM110
ARM2203	REM211
ARM2204	REM216
ARM2205	REM400
ARM2206	REM500
ARM1707	Swing away RH Bracket
ARM1708	Swing away LH Bracket
ARM2213	Centre mount, swing away to RH Nucleus
ARM2214	Centre mount, swing away to LH Nucleus
ARM2227	RH Bracket mounted for Attendant
ARM2228	LH Bracket mounted for Attendant
ARM2209	REM500 Swing away RH multi-adjustable Bracket
ARM2210	REM500 Swing away LH multi-adjustable Bracket
ARM2270	Assignable toggle switch module
ARM2220	Dual control REM050
ARM2219	Compact (ACU200)
ARM2245	IDC Dual control
ARM1221	IDC Dual control detachable Bracket RH
ARM1222	IDC Dual control detachable Bracket LH
ARM2231	Proportional Compact Remote
ARM2232	Proportional Compact Remote low force
ARM2207	RH mount Bracket
ARM2208	LH mount Bracket
ARM2213	Centre mount, swing away to RH Nucleus Bracket
ARM2214	Centre mount, swing away to LH Nucleus Bracket

ARM2229	Proportional Chin joystick (ASL138)
ARM1289	Swing away RH Manual Bracket
ARM1288	Swing away LH Manual Bracket
ARM1287	Swing away RH Powered Bracket
ARM1286	Swing away LH Powered Bracket
ARM2230	ATOM with ultrasonic proximity sensor
ARM2247	PROTON with ultrasonic proximity sensor
ARM2257	Four-switch Array (ASL106)
ARM2275	Centre mount, swing-away to RH nucleus
ARM2276	Centre mount, swing-away to LH nucleus
ARM2233	MEC proportional joystick
ARM2240	Proportional heavy duty with mode switch (ASL133)
ARM2235	MEC tray mount RH (ASL603)
ARM2236	MEC tray mount LH (ASL603)
ARM2213	Centre mount, swing away to RH nucleus
ARM2214	Centre mount, swing away to LH nucleus
ARM2241	Centre mount, swing away to RH nucleus
ARM2242	Centre mount, swing away to LH nucleus
ARM2217	RH Bracket
ARM2218	LH Bracket
ARM2227	RH Bracket
ARM2228	LH Bracket
ARM2253	USB charger port DLX-USB02-A
ARM2207	RH mount
ARM2208	LH mount
ARZ2207	RH mount
ARZ2208	LH mount
ARM1282	6 km/h (4mph)
ARM1908	10 km/h (6mph)
ARM1437	LED headlights & indicators
ARM1254	Without light system
ARM1329	Lap Belt fitted to base of backrest
ARM2113	Lap Belt adjustable mounting on seat rail extrusion
ARM1332	H-Belt
ARM2701	Evoflex pelvic stabilizer, size L
ARM2702	Evoflex pelvic stabilizer, size M
ARM2703	2-Point hip belt centre pull, size L
ARM2704	2-Point hip belt centre pull, size M
ARM2705	2-Point hip belt dual pull, size L
ARM2706	2-Point Hip Belt Dual pull, size M
ARM1415	RH Mirror
ARM1339	LH Mirror
ARM2128	Bag hooks
ARM1312	RH Cane holder

---

ARM1363	LH Cane holder
ARM1404	Swing away RH Tray
ARM1405	Swing away LH Tray

---



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):