



Yes, you can:

Tel: 0800 043 0881
Email: salesUK@dhcare.com
Online Spares available at: www.dhcare.co.uk

Power Retail Prescription Form

September 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 NB Teen

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Highly adjustable Modulite Seating System
- Seat widths from 305mm - 380mm (12" - 15")
- Seat depths 305mm - 560mm (12" - 22"), stepless adjustment (12" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 20°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- Standard 80° hangers with fixed footplate
- LiNX DLX-REM110 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Stability Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 50 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,995



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARZ0000	TDX SP2 NB Modulite Teen	£6,995
ARZ1300	Special TDX SP2 NB Modulite Teen (Extended Lead Time). Please provide details in comments section.	◆
SEATING TYPE*		
ARM2100	Modulite Seat Unit	◆
SEAT FRAME*		
ARM2111	Telescopic Seat Frame + Plate with Side Rails.	◆

SEAT WIDTH*				
ARZ0001		305mm		○
ARZ0002		330mm		○
ARM2121		380mm		○
SEAT DEPTH RAIL LENGTH*				
ARM2129		Standard Depth Package (410mm - 510mm. Simple stepless adjustment from the side)		○
ARM2190		Extended Depth Package (510mm - 560mm. Simple stepless adjustment from the side)	(NC 2134)	○
SEAT DEPTH*				
ARZ0003		305mm		○
ARZ0004		330mm		○
ARZ0005		355mm		○
ARZ0006		380mm		○
ARM2132		410mm		○
ARM2134		460mm		○
ARM2136		510mm		○
ARM2138		560mm		○
SEAT OPTION				
ARZSWO		"H" Frame Interface Kit	(NC 0002)	£495
SEAT CUSHION*				
ARM2142		Flo-tech Cushion (Lite - 75mm up to 127kg)	(NC 0001, 0002)	£120
ARM2147		Matrx Libra (90mm - standard). Extened lead time.	(NC 0001, 0002)	£410
1626906		Matrx Mini series - PS cushion with Kidabra print 310mm x 310mm, 12" x 12" (max user weight 68kg)		£232
1626907		Matrx Mini series - PS cushion with Kidabra print 310mm x 360mm, 12" x 14" (max user weight 68kg)		£232
1626908		Matrx Mini series - PS cushion with Kidabra print 310mm x 410mm, 12" x 16" (max user weight 68kg)		£232
1627456		Matrx Mini series - PS cushion with Kidabra print 330mm x 310mm, 13" x 12" (max user weight 68kg)		£232
1627457		Matrx Mini series - PS cushion with Kidabra print 330mm x 360mm, 13" x 14" (max user weight 68kg)		£232
1627458		Matrx Mini series - PS cushion with Kidabra print 330mm x 410mm, 13" x 16" (max user weight 68kg)		£232
1626909		Matrx Mini series - PS cushion with Kidabra print 360mm x 310mm, 14" x 12" (max user weight 68kg)		£232
1626910		Matrx Mini series - PS cushion with Kidabra print 360mm x 360mm, 14" x 14" (max user weight 68kg)		£232
1626911		Matrx Mini series - PS cushion with Kidabra print 360mm x 410mm, 14" x 16" (max user weight 68kg)		£232
ARM2140		Without seat cushion		◆
1st POSITIONING ACCESORIES FOR MATRX LIBRA				
ARM1980		Tapered Wedges (Pair)		£30
ARM1981		Abductor		£23
ARM1982		Lateral Wedge		£20
ARM1983		Anterior Seat Wedges		£27
2nd POSITIONING ACCESORIES FOR MATRX LIBRA				
ARM1980		Tapered Wedges (Pair)		£30
ARM1981		Abductor		£23
ARM1982		Lateral Wedge		£20
ARM1983		Anterior Seat Wedges		£27

BACKREST SUPPORT*			
ARM2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrux Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	◆
ARM2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrux Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. (NC 2173)	○
PUSH HANDLES FOR BACKPOSTS*			
ARM2107		Backposts WITH push handles	◆
ARM2108		Backposts WITHOUT push handles. Push handles will be cut down on the vertical section. (NC 2219, 2227, 2228, 2245, 2171, 1495)	○
BACKREST UPHOLSTERY / CUSHION*			
ARM2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery. (NC 0002, 2155, 1496, 1445)	◆
ARM2172		Flex3 adjustable Backrest Plates (without Backrest Cushion)	£250
ARM2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion	£400
SP1562143		Mistral Flex3 Backrest shell - shorter & narrower than the Laguna Flex3 (NC 0002)	£159
SP1431668-TR26		Mistral Backrest, multi stretch PU coated fabric Black Cushion (NC 0002)	£145
SP1431668-TR35		Mistral Backrest, Black Plush Cushion (NC 0002)	£163
EMP		Extended Mounting Pins	£519
MINI1207		Matrx Mini series - Mini Back with Kidabra print 310mm x 180mm, 12" x 7" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1210		Matrx Mini series - Mini Back with Kidabra print 310mm x 260mm, 12" x 10" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1212		Matrx Mini series - Mini Back with Kidabra print 310mm x 310mm, 12" x 12" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1214		Matrx Mini series - Mini Back with Kidabra print 310mm x 360mm, 12" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1307		Matrx Mini series - Mini Back with Kidabra print 330mm x 180mm, 13" x 7" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1310		Matrx Mini series - Mini Back with Kidabra print 330mm x 260mm, 13" x 10" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1312		Matrx Mini series - Mini Back with Kidabra print 330mm x 310mm, 13" x 12" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1314		Matrx Mini series - Mini Back with Kidabra print 330mm x 360mm, 13" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1407		Matrx Mini series - Mini Back with Kidabra print 360mm x 180mm, 13" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1410		Matrx Mini series - Mini Back with Kidabra print 360mm x 260mm, 13" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1412		Matrx Mini series - Mini Back with Kidabra print 360mm x 310mm, 13" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1414		Matrx Mini series - Mini Back with Kidabra print 360mm x 360mm, 14" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£44
ARM2170		Without Backrest Upholstery (Back Posts/Back Frame only) (NC 2151, 2152, 2155)	○

BACKREST OPTIONS			
SP1571542		Angle adjustable backrest support brackets RH	£44
SP1571553		Angle adjustable backrest support brackets LH	£36
UPHOLSTERY COLOUR*			
ARM2151		Black Nylon Fabric TR18	◆
ARM2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26 (NC 2170)	○
ARM2155		Black Spacer Fabric TR36 (NC 2147, 2171)	○
ARZ1004		Matrx Mini In Kidabra fabric (NC 2172, 2173)	○
ARZ1003		Matrx Mini In black fabric (NC 2172, 2173)	○
ARX1002		No Upholstery Colour required - No backrest upholstery required (NC 2172, 2173)	○
ARMRESTS*			
ARM2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm	◆
ARM2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm	£40
ARM2183		Inrigged Armrests. Reduces SW by 25 mm each side	£30
SP1582159		Inrigged Armrests for 12" seat width, Standard inrigged arm option MUST BE SELECTED with additional brackets (must be selected with 12" seat width)	£45
ARMPADS*			
ARM2195		Short Armpad (325mm)	◆
ARM2196		Long Armpad (375mm)	○
HIP SUPPORTS LH			
ARM2184		Hip Support LH, pad size 140 x 100 mm (NC 0556, 0589)	£250
ARM2186		Hip Support LH, pad size 230 x 100 mm (NC 0556, 0589)	£250
HIP SUPPORTS RH			
ARM2185		Hip Support RH, pad size 140 x 100 mm (NC 0556, 0589)	£250
ARM2187		Hip Support RH, pad size 230 x 100 mm (NC 0556, 0589)	£250
MOUNTING HARDWARE FOR HIP SUPPORTS			
ARM2164		Mounting hardware for hip supports (6.1" / 155 mm high)	○
ARM2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
HEADREST / NECKREST - Adequate headrest option must be selected*			
1603255		Elan Standard Pad - Small (160mm x 90mm) Standard cover (please select mounting Hardware type for Elan Headrest)	£124
1603254		Elan Standard Pad - Medium (260mm x 130mm) Standard cover (please select mounting Hardware type for Elan Headrest)	£124
1603291		Elan Occipital Pad (230mm x 180mm) Standard cover (please select mounting Hardware type for Elan Headrest)	£232
1629822		Elan 4 Point Headrest Pad (Small) Standard cover (please select mounting Hardware type for Elan Headrest)	£184
ARM1432		Headrest sideways adjustable	£225
ARM1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	○
ARM1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○

ELAN MOUNTING HARDWARE				
1603250		Matrx Elan 6-point Hardware Infinite adjustments		£184
1603251		Matrx Elan 6-point Mini (paediatric 100mm shorter post)		£166
1603249		Matrx Standard multi-axis hardware Limited adjustability		£166
1603252		Elan adapter for Flex3 back		£140
LATERAL SUPPORT				
ARM1445		Lateral Support with extended arm, multi-adjustable pad	(NC 2171)	£150
MATRX MINI LATERAL SUPPORTS				
1627466		Swing-away lateral Hardware (each)		£86
1544292		Extra small lateral Pad (each)		£59
1627465		Chest harness interface Kit		£70
MANUAL / POWERED SEATING OPTIONS*				
ARM1091		Manually settable Seat Tilt 0° - 20°		◆
ARM1092		Powered "fixed pivot" Tilt 20° (136kg max user weight)		£400
ARM1093		Powered Tilt (0° - 30°) with CoG shift (160kg max user weight)	(NC 2201)	£1,100
ARM1094		Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (160kg max user weight)	(NC 2201)	£1,900
SEAT TO GROUND HEIGHT*				
ARM0354		450mm (17.75") (from floor to front of seat plate at 0° tilt)	(NC 1094)	○
ARM0364		475mm (18.75") (from floor to front of seat plate at 0° tilt)		○
SEAT TILT OPTION				
ARM2054		5° Fixed Anterior Tilt, Assist mount. Will reduce max tilt angle by 5°		£320
BACKREST RECLINER*				
ARM0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable		◆
ARM0977		Powered Recliner (0° - 30°) Foldable		£500
LEG REST TYPE*				
ARM0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0504, 0513, 0526, 0555, 0554)	£370
ARM0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 0512, 0504, 0513, 0526, 0555, 0554)	£1,620
ARM0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm	(NC 0559, 0504)	◆
1606173-0601		Side Mounted: Black short 90° Legrest Hanger x 2 - Action style (please select legrest support carriers below)	(NC 0512, 0504, 0513)	£136
1606172-0601		Side Mounted: Black short 80° Legrest Hanger x 2 - Action style (please select legrest support carriers below)	(NC 0512, 0504, 0513)	£136
ARM0501		Without Leg Rests	(NC 0559, 0504, 0513)	○
LEGREST SUPPORT (Only applicable for Action Style Hangers)				
SP1531195		Legrest support carrier LH (to be selected for the Action style hangers, 1606173-0601 and 1606172-0601)		£39
SP1531196		Legrest support carrier RH (to be selected for the Action style hangers, 1606173-0601 and 1606172-0601)		£39

FOOTPLATE TYPE*				
ARM0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	(NC 0501)	○
ARM0512		Fixed Footplate pair	(NC 0556, 0589, 0501)	◆
ARM0504		Angle Adjustable Footplate pair	(NC 0589, 0501, 0526)	£75
ARM0513		Angle and depth adjustable Footplate pair	(NC 0589, 0501)	○
ARM0501		Without footplates for option 0501 without Leg Rests	(NC 0589, 0526)	○
FOOTPLATE TYPE (Only applicable for Action Style Hangers)				
1428660		Angle and depth adjustable Footplate LH, including stem (for use with short 80° & 90° Action style Legrest)		£37
1428654		Angle and depth adjustable Footplate RH, including stem (for use with short 80° & 90° Action style Legrest)		£37
1530620		Fixed Footplate LH, including stem (for use with short 80° & 90° Action style Legrest)		£19
1530614		Fixed Footplate RH, including stem (for use with short 80° & 90° Action style Legrest)		£19
1492150		Upmounted Footplate LH (for use with short 80° & 90° Action style Legrest)		£74
1492149		Upmounted Footplate RH(for use with short 80° & 90° Action style Legrest)		£63
AMPUTEE SUPPORT				
ARM0554		Amputee Support RH	(NC 0556, 0589)	£195
ARM0555		Amputee Support LH	(NC 0556, 0589)	£195
LiNX CONTROLS*				
ARM2201		DLX-REM110 (Drive Only)	(NC 1093, 1094, 0977, 0589, 1437, 2209, 2210, 2219, 2245, 2270)	◆
ARM2203		DLX-REM211 (Drive & Actuators)	(NC 2209, 2210, 2270)	○
ARM2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2210, 2270)	£200
ARM2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2210)	£425
ARM2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1208, 1707, 1708, 1709, 1710, 2213, 2214, 2227, 2228, 2220, 2270)	£765
REM400 TOGGLE OPTION				
ARM2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204)	£315
CONTROL BRACKETRY*				
ARM1207		Standard Bracket RH	(NC 2206)	○
ARM1208		Standard Bracket LH	(NC 2206)	○
ARM1707		Swing Away RH Bracket	(NC 2206)	£135
ARM1708		Swing Away LH Bracket	(NC 2206)	£135
ARM1709		Matrx Resolve Swing in/out Bracket RH	(NC 2213, 2213, 1404)	£245
ARM1710		Matrx Resolve Swing in/out Bracket LH	(NC 2214, 2214, 1405)	£245
ARM2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 1405)	£495
ARM2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 1405)	£495
ARM2209		REM500 Swing Away RH Multi-Adjustable Bracket		○
ARM2210		REM500 Swing Away LH Multi-Adjustable Bracket		○
ARM2227		Bracket RH mounted for Attendant		£65
ARM2228		Bracket LH mounted for Attendant		£65
ARZ2239		Individual mounting option per request (Special). Requires additional information/specification.		POA
DUAL CONTROLS				
ARM2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 1222, 2227, 2228, 9260)	£395
ARM2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1222, 2227, 2228, 2201, 9260)	£395
ARM2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2220, 2227, 2228, 9260)	£745
ARM2420		mo-vis Scoot Control w/ dedicated bracket included.	(NC 2108, 1221, 1222)	£3,420

DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
ARM1207		Dual Control Standard Bracket RH		○
ARM1208		Dual Control Standard Bracket LH		○
ARM1221		IDC Dual Control Detachable Bracket RH	(NC 2220, 2227)	○
ARM1222		IDC Dual Control Detachable Bracket LH	(NC 2220, 2228)	○
ARM9260		Bracket for mo-vis Scoot	(NC 2219, 2245)	○
COMPACT REMOTES				
ARM2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2220)	£825
ARM2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2220)	£1,075
ARM3000		mo-vis Micro LiNX	(NC 2201, 2203, 2204)	£2,960
ARM3001		mo-vis Multi LiNX	(NC 4401, 4402, 2201, 2203, 2204)	£2,640
ARM3110		mo-vis All Round Low Force Omni	(NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
ARM3111		mo-vis All Round Omni	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
ARM3112		mo-vis Heavy Duty Omni	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4300, 4301, 4302, 4400, 4401, 4402, 2201, 2203, 2204)	£3,000
COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
ARM2207		RH Mount	(NC 2201, 2203, 2204)	○
ARM2208		LH Mount	(NC 2201, 2203, 2204)	○
ARM2213		Centre Mount Swing Away to RH (Nucleus)	(NC 1405, 2201, 2203, 2204)	£495
ARM2214		Centre Mount Swing Away to LH (Nucleus)	(NC 1405, 2210, 2203, 2204)	£495
ARZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
ARM3300		Maxx Resolve RH	(NC 4301, 4302)	£270
ARM3301		Maxx Resolve LH	(NC 4301, 4302)	£270
ARM3310		BodyPoint Trilock Midline RH	(NC 4301, 4302)	£480
ARM3311		BodyPoint Trilock Midline LH	(NC 4301, 4302)	£480
ARM3320		Tray Mounted	(NC 4301, 4302)	○
ARZ3330		Individual mounting option per request (Special). Requires additional information/specification.		POA
CHIN CONTROLS				
ARM3000		mo-vis Micro LiNX	(NC 2201, 2203, 2204)	£2,960
ARM3001		mo-vis Multi LiNX	(NC 4401, 4402, 2201, 2203, 2204)	£2,640
ARM3110		mo-vis All Round Low Force Omni (consider option 3810 for switch access)	(NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
ARM3111		mo-vis All Round Omni (consider option 3810 for switch access)	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
ARM3400		mo-vis Chin Control Harness (One Size)	(NC 4301, 4302)	£1,000
ARM3410		mo-vis Manual Lateral Swing Away RH Mounted	(NC 4301, 4302)	£1,520
ARM3411		mo-vis Manual Lateral Swing Away LH Mounted	(NC 4301, 4302)	£1,520
ARM3412		mo-vis Manual Vertical Swing Away RH Mounted	(NC 4301, 4302)	£1,520
ARM3413		mo-vis Manual Vertical Swing Away LH Mounted	(NC 4301, 4302)	£1,520
ARM3420		mo-vis Powered Lateral Swing Away RH Mounted	(NC 4301, 4302)	£2,400
ARM3421		mo-vis Powered Lateral Swing Away LH Mounted	(NC 4301, 4302)	£2,400
ARM3422		mo-vis Powered Vertical Swing Away RH Mounted	(NC 4301, 4302)	£2,400
ARM3423		mo-vis Powered Vertical Swing Away LH Mounted	(NC 4301, 4302)	£2,400
ARZ3330		Individual mounting option per request (Special). Requires additional information/specification.		POA
HEAD CONTROLS				
ARM3200		mo-vis Proportional Head Control Omni	(NC 2201, 2203, 2204)	£4,500
ARM3221		mo-vis Switched Head Control Omni	(NC 2201, 2203, 2204)	£4,500
ARM3222		mo-vis Switched Head Control w/ Sip & Puff Omni	(NC 2201, 2203, 2204)	£4,750

FOOT CONTROL			
ARM3010	☐	mo-vis Foot Control LiNX RH Side (NC 2203, 2204, 2220)	£3,300
ARM3011	☐	mo-vis Foot Control LiNX LH Side (NC 2203, 2204, 2220)	£3,300
LiNX 10-WAY SWITCH			
ARM2251	☐	10-Way Switch DLX-SW10-A for direct Powered Seating Access	£500
LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
ARM2217	☐	Bracket RH, mounted for user (NC 2213)	☐
ARM2218	☐	Bracket LH, mounted for user (NC 2214)	☐
ARM2227	☐	Bracket RH, mounted for attendant (NC 1221)	£65
ARM2228	☐	Bracket LH, mounted for attendant (NC 1222)	£65
LiNX ACCESSORY CABLE			
ARM2252	☐	Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£100
LiNX USB CHARGER PORT			
ARM2253	☐	USB Charger Port DLX-USB02-A	£160
LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
ARM2207	☐	RH Mount	☐
ARM2208	☐	LH Mount	☐
FURTHER LiNX / MO-VIS COMPONENTS / ACCESSORIES			
ARM2302	☐	Not Necessary (NC 2254, 1711, 1712, 1713, 1714, 1715, 1716)	☒
ARM2300	☐	To be enclosed only. Standard 300 mm Cable Length for mo-vis Switches	☐
ARM2305	☐	To be enclosed only. Extended 1500 mm Cable Length for mo-vis Switches	☐
ARZ2301	☐	To be connected and programmed (Special)	POA
BUDDY BUTTONS			
ARM1685	☐	Buddy Button. Red / Yellow / Green / Blue / Black provided. Interchangeable colours.	£90
EGG SWITCH			
ARM2254	☐	Egg Switch, black	£90
PIKO BUTTONS			
ARM1711	☐	Red Piko Button	£160
ARM1712	☐	Yellow Piko Button	£160
ARM1713	☐	Green Piko Button	£160
ARM1714	☐	Orange Piko Button	£160
ARM1715	☐	Blue Piko Button	£160
ARM1716	☐	Black Piko Button	£160

TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle - direct control for Tilt only	£270
ARZ2262		Stereo Button - direct control for Tilt only	£255
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	£330
ARZ2264		Four-Way Button - direct control for Tilt & Recline	£355
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
ARZ2261		Stereo Toggle - direct control for Tilt only	£270
ARZ2262		Stereo Button - direct control for Tilt only	£255
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	£330
ARZ2264		Four-Way Button - direct control for Tilt & Recline	£355
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
LINX INPUT MODULE			
UNZ2273		DLX-IN200 Module (Jack port for external power button) (NC 2220)	£270
UNZ9249		DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	£785
MO-VIS TWISTER BASIC			
ARM3500		Black	£140
ARM3501		Red	£140
ARM3502		Yellow	£140
ARM3503		Blue	£140
ARM3504		Green	£140
MO-VIS TWISTER BASIC MOUNTING (Only applicable for mo-vis Twister Basic)			
ARM3520		Mounting Set	£50
ARM3525		Table Mounting Set	£50
MO-VIS TWISTER			
ARM3600		Black	£160
ARM3601		Red	£160
ARM3602		Yellow	£160
ARM3603		Blue	£160
ARM3604		Green	£160
MO-VIS TWISTER ON BENDED TUBE INC. D6 CLAMP			
ARM3700		Black	£220
ARM3701		Red	£220
ARM3702		Yellow	£220
ARM3703		Blue	£220
ARM3704		Green	£220

MO-VIS SATELLITE TWISTER (Only applicable for mo-vis All Round Remotes)			
ARM3800		Black	£200
ARM3801		Red	£200
ARM3802		Yellow	£200
ARM3803		Blue	£200
ARM3804		Green	£200
ARM3810		All Round Chin Control Set (includes All Round Joystick Ball and 2x Satellite Switches)	£420
MO-VIS TWISTER GOOSENECK			
ARM3900		Black	£260
ARM3901		Red	£260
ARM3902		Green	£260
MO-VIS Q2M TWIST & SWING			
ARM4000		Black	£500
ARM4001		Red	£500
ARM4002		Yellow	£500
ARM4003		Blue	£500
ARM4004		Green	£500
MO-VIS TWISTER PRO			
ARM4100		Black	£160
ARM4101		Red	£160
ARM4102		Green	£160
MO-VIS MULTI SWITCH			
ARM4200		Multi Switch	£620
ARM4201		Multi Switch Proximity Sensors	£860
ARM4202		Multi Switch Proximity Set	£840
MO-VIS HANDWARMER			
ARM4300		Hand Warmer	£800
ARM4301		Hand Warmer 2 Dome Hand Hood	£400
ARM4302		Hand Warmer 1 Dome Hand Hood	£300
MO-VIS JOYSTICK KNOBS			
ARM4400		Micro Joystick Ball	£60
ARM4401		Micro Joystick Cup 10.1mm	£60
ARM4402		Micro Joystick Cup 15.5mm	£60
ARM4403		Multi Joystick Ball 4mm	£60
ARM4404		Multi Joystick Cup 24mm	£60
ARM4405		All Round Joystick Ball 6.35mm	£60
ARM4406		All Round Joystick Cup 37mm	£60
BODYPOINT JOYSTICK KNOBS			
ARM2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£75
ARM2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£75
ARM2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
ARM2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75

SHROUD COLOUR*			
ARM2640		Obsidian Black Frame Shrouds	○
ARM2641		Ruby Red Frame Shrouds	○
ARM2642		Purple Passion Frame Shrouds	○
ARM2643		Peacock Blue Frame Shrouds	○
ARM2644		Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARM2650		Obsidian Black Rim Inserts	○
ARM2651		Ruby Red Rim Inserts	○
ARM2652		Purple Passion Rim Inserts	○
ARM2653		Peacock Blue Rim Inserts	○
ARM2654		Arctic Blue Rim Inserts	○
FORKS*			
ARM2632		Double sided fork	◆
CENTRE WHEELS (DRIVING WHEELS)*			
ARM0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	○
MOTOR TYPE*			
ARM1294		3 km/h (2mph)	○
ARM1282		6 km/h (4mph)	◆
ARM1908		10 km/h (6mph) (includes LED light system and indicators)	£515
BATTERY*			
ARM1257		50 Ah Gel	◆
G TRAC			
ARM1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£400
CHARGER*			
ARM1343		OptiCharge Charger 8A (200 - 250V)	◆
ARM1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£155
LIGHTS*			
ARM1437		LED Light System with Indicators (30 LUX)	£450
ARX1006		Lights Included with Motor Package 10 & 12kmh (NC 1282)	○
ARM1254		Without Light system	◆
TRANSPORTATION KIT*			
ARM1450		Preparation for DAHL Docking Station: Docking station for car implementation, please order docking station separately. (Adequate headrest is compulsory.	£580
ARM1493		4-Point tie-down (For wheelchair transportation without user)	○
ARM1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
FIXATION BELT*			
ARM1329		Lap Belt with Auto Lock	◆
ARM2113		Lap belt with auto lock (flexible mount in side rails) (not available in combination with 410mm seat depth (2132) and T-armrest (2191)	○

BODYPOINT BELTS		
HB216-S38-B2		HB216 2-point belt Small x 380mm Pads, cinch mount £63
HB216-M46-B2		HB216 2-point belt Medium x 460mm Pads, cinch mount £75
HB219-S38-B2		HB219 2-point belt Small x 380mm Pads, cinch mount £75
HB219-M46-B2		HB219 2-point belt Medium x 460mm Pads, cinch mount £75
HB235-S32-B2		HB235 2-point belt Extra small x, 320mm Pads, cinch mount £75
HB235-S38-B2		HB235 2-point belt Small x 380mm Pads, cinch mount £81
EB235S		Evoflex Small pelvic stabilizer, push button £81
EB225M		Evoflex Medium pelvic stabilizer, Latch £81
EB225L		Evoflex Large pelvic stabilizer, Latch £110
BODYPOINT PIVOTFIT		
SH290XXS-B6		X-Extra small standard PivotFit with multi-directional design and cinch mount £110
SH290XS-B6		Extra small standard PivotFit with multi-directional design and cinch mount £110
SH290S-B6		Small standard PivotFit with multi-directional design and cinch mount £110
SH290M-B6		Medium standard PivotFit with multi-directional design and cinch mount £114
SH295XXS-B6		X-Extra small dynamic PivotFit with multi-directional design and cinch mount £114
SH295XS-B6		Extra small dynamic PivotFit with multi-directional design and cinch mount £114
SH295S-B6		Small dynamic PivotFit with multi-directional design and cinch mount £114
SH295M-B6		Medium dynamic PivotFit with multi-directional design and cinch mount £92
BODYPOINT TRIMLINE		
SH235S-B6		Small Trimeline shoulder harness with multi-directional straps and cinch mount £92
SH235M-B6		Medium Trimeline shoulder harness with multi-directional straps and cinch mount £133
BODYPOINT STAYFLEX		
SH351XS-B6		Extra small standard Stayflex without zip, with multidirectional straps and plastic cinch mounts £133
SH351S-B6		Small standard Stayflex without zip, with multidirectional straps and plastic cinch mounts £133
SH360M-B6		Medium narrow Stayflex with zip, with multidirectional straps and cinch mounts £133
SH361M-B6		Medium narrow Stayflex without zip, with multidirectional straps and cinch mount £133
SH350M-B6		Medium standard Stayflex with zip, with multifunctional straps and plastic cinch mounts £76

BODYPOINT MONOFLEX			
SH120XS-B3		Extra Small, front pull, grommet strap	£76
SH120S-B3		Small, front pull, grommet Strap	£76
SH120M-B3		Medium, front pull, grommet strap	£76
SH120XS-B1		Extra Small, Monoflex underarm release, flat mount	£76
SH120S-B1		Small, Monoflex underarm release, flat mount	£76
SH120M-B1		Medium, Monoflex underarm release, flat mount	£76
SH122XS-B1		Extra-small Monflex centre release, flat mount	£76
SH122S-B1		Small Monflex centre release, flat mount	£76
SH122M-B1		Medium Monflex centre release, flat mount	£59
BODYPOINT HARNESSES			
1627465		Chest harness interface kit	£70
CANE HOLDER			
ARM1312		Cane holder RH	£45
ARM1363		Cane holder LH	£45
TRAY			
ARM1404		Tray, Swing Away RH (NC 1707, 2213, 2214, 2209)	£200
ARM1405		Tray, Swing Away LH (NC 1708, 2213, 2214, 2210)	£200
TOTAL PRICE			

Non Conforming Options:	
ARZ0001	305mm
ARZ0002	330mm
ARM2134	460mm
ARM2147	Matrx Libra (90mm - standard). Extened lead time.
ARM2108	Backposts WITHOUT push handles. Push handles will be cut down on the vertical section.
ARM2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
ARM2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ARM2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ARM2170	Without Backrest Upholstery (Back Posts/Back Frame only)
ARM2151	Black Nylon Fabric TR18
ARM2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
ARM2155	Black Spacer Fabric TR36
ARM1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
ARM1495	Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts
ARM1445	Lateral Support with extended arm, multi-adjustable pad
ARM1093	Powered Tilt (0° - 30°) with CoG shift (160kg max user weight)
ARM1094	Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (160kg max user weight)
ARM0977	Powered Recliner (0° - 30°) Foldable
ARM0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
ARM0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
ARM0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm
ARM0501	Without Leg Rests
ARM0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
ARM0512	Fixed Footplate pair
ARM0504	Angle Adjustable Footplate pair
ARM0513	Angle and depth adjustable Footplate pair
ARM0501	Without footplates for option 0501 without Leg Rests
ARM0554	Amputee Support RH
ARM0555	Amputee Support LH
ARM2201	DLX-REM110 (Drive Only)
ARM2203	DLX-REM211 (Drive & Actuators)
ARM2204	DLX-REM216 (Drive, Actuators & Lights)
ARM2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
ARM2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
ARM1208	Standard Bracket LH
ARM1707	Swing Away RH Bracket
ARM1708	Swing Away LH Bracket
ARM1709	Matrx Resolve Swing in/out Bracket RH
ARM1710	Matrx Resolve Swing in/out Bracket LH
ARM2213	Centre Mount Swing Away to RH Nucleus Bracket
ARM2214	Centre Mount Swing Away to LH Nucleus Bracket

ARM2209	REM500 Swing Away RH Multi-Adjustable Bracket
ARM2210	REM500 Swing Away LH Multi-Adjustable Bracket
ARM2227	Bracket RH mounted for Attendant
ARM2228	Bracket LH mounted for Attendant
ARM2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
ARM2219	LiNX DLX-ACU200 (Compact Dual Control)
ARM2245	LiNX DLX-IDC (Intuitive Dual Control)
ARM1208	Dual Control Standard Bracket LH
ARM1221	IDC Dual Control Detachable Bracket RH
ARM1222	IDC Dual Control Detachable Bracket LH
ARM9260	Bracket for mo-vis Scoot
ARM2213	Centre Mount Swing Away to RH (Nucleus)
ARM2214	Centre Mount Swing Away to LH (Nucleus)
ARM3410	mo-vis Manual Lateral Swing Away RH Mounted
ARM3411	mo-vis Manual Lateral Swing Away LH Mounted
ARM3412	mo-vis Manual Vertical Swing Away RH Mounted
ARM3413	mo-vis Manual Vertical Swing Away LH Mounted
ARM3420	mo-vis Powered Lateral Swing Away RH Mounted
ARM3421	mo-vis Powered Lateral Swing Away LH Mounted
ARM3422	mo-vis Powered Vertical Swing Away RH Mounted
ARM3423	mo-vis Powered Vertical Swing Away LH Mounted
ARM3200	mo-vis Proportional Head Control Omni
ARM3221	mo-vis Switched Head Control Omni
ARM3222	mo-vis Switched Head Control w/ Sip & Puff Omni
ARM2227	Bracket RH, mounted for attendant
ARM2228	Bracket LH, mounted for attendant
ARM2254	Egg Switch, black
ARM1711	Red Piko Button
ARM1712	Yellow Piko Button
ARM1713	Green Piko Button
ARM1714	Orange Piko Button
ARM1715	Blue Piko Button
ARM1716	Black Piko Button
ARM4300	Hand Warmer
ARM4301	Hand Warmer 2 Dome Hand Hood
ARM4302	Hand Warmer 1 Dome Hand Hood
ARM4400	Micro Joystick Ball
ARM4401	Micro Joystick Cup 10.1mm
ARM4402	Micro Joystick Cup 15.5mm
ARM1282	6 km/h (4mph)
ARM1437	LED Light System with Indicators (30 LUX)
SH122XS-B1	Extra-small Monflex centre release, flat mount
SH122S-B1	Small Monflex centre release, flat mount
SH122M-B1	Medium Monflex centre release, flat mount

ARM1404	Tray, Swing Away RH
ARM1405	Tray, Swing Away LH



For a more complex customer requirement,
please email: salesUK@dhcare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):