



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

July 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 NB Ultra

Standard Features

- Max User Weight 136kg (21.4 Stone)
- Crash tested ISO7176- 19
- Ultra Low Maxx Seating System
- Seat widths of 405mm - 610mm (16" - 24"), adjustable between 405mm - 505mm (16" - 20") and 505mm - 610mm (20" - 24")
- Seat depths of 380mm - 580mm (15" - 23") adjustable between 380mm - 480mm (15" - 19") and 480mm - 580mm (19 - 23")
- Backrest heights of 560mm - 660mm in 25mm increments
- Cantilever flip-back arm, width, height and angle adjustable with desk armpad with startex covered visco foam
- Elan multi positioning headrest
- Ultra Low Maxx CG power tilt (0° to 50°) seat heights of 420mm - 480mm
- Powered Recline 90° to 170° with extended shear reduction (ESR)
- Vari F hangers with fixed footplates
- LiNX DLX-REM211 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (driving wheels) 355mm x 76mm (14" x 3") Black tyres
- 10km/h (6mph) 4 pole motors and 50 Ah batteries
- G-Trac
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARP0080	TDX SP2 Ultra Narrow Base	◆
ARU1300	Special TDX SP2 Ultra Narrow Base (Extended Lead Time). Please provide details in comments section.	POA
SEATING TYPE*		
ARP2020	Ultra Low Maxx Powered Positioning System	£13,500

SEAT ADJUSTMENTS*				
ARP2052		Centre of gravity Tilt (0 to 50°)		◆
ARP2053		300mm seat lifter (increased STF height by 4mm) with centre of gravity tilt (0 to 45°)	(NC 0339)	£2,500
SEAT TO GROUND HEIGHT*				
ARP0337		440mm		○
ARP0331		460mm		○
ARP0339		480mm	(NC 2053)	○
ANTERIOR AND POSTERIOR TILT*				
ARP2057		Tilt Standard Mount (No anterior or posterior tilt)		◆
ARP2054		5° Fixed Anterior Assist Mount (will reduce max tilt angle by 5°)		£320
ARP2055		10° Fixed Anterior Assist Mount (will reduce max tilt angle by 10°)	(NC 2035)	£400
ARP2056		5° Pre-Tilt (will increase max tilt angle by 5°)		£320
BACKREST RECLINER*				
ARP2060		Powered Recliner (90° - 170°) with Extended Shear Reduction (ESR), Standard Mount		◆
ARP2063		Powered Recliner 8° Recline mount 82° to 162° with Extended Shear Reduction (ESR)		£35
ARP2066		Powered Recline, 30° Recline mount, 60° to 140° with Extended Shear Reduction (ESR)		£35
SEAT FRAME WIDTH RANGE*				
ARP2021		405mm - 505mm (increments of 25mm) (16" - 20")	(NC 2030, 2031, 2032, 2033)	○
ARZ2003		480mm - 555mm (Increments of 25mm) (19" - 22") - Additional build time of 15+days	(NC 2025, 2026, 2027, 2032, 2033)	○
ARZ2022		505mm - 610mm (increments of 25mm) (20" - 24") - Additional build time of 15+days	(NC 2025, 2026, 2027, 2174, 2175)	○
SEAT FRAME DEPTH RANGE*				
ARP2024		480mm - 580mm (increments of 25mm) (19" - 23")	(NC 2035, 2036, 2037)	○
ARP2008		405mm - 505mm (increments of 25mm) (16" - 20")	(NC 2040, 2041, 2042)	○
SEAT WIDTH*				
ARP2025		405mm (16")	(NC 2003, 2022, 2041, 2042)	○
ARP2026		430mm (17")	(NC 2003, 2022, 2041, 2042)	○
ARP2027		455mm (18")	(NC 2003, 2022)	○
ARP2028		480mm (19")	(NC 2022, 2174, 2175)	○
ARP2029		505mm (20")	(NC 2174, 2175)	○
ARP2030		530mm (21")	(NC 2021, 2035, 2036, 2174, 2175)	○
ARP2031		555mm (22")	(NC 2021, 2035, 2036, 2174, 2175)	○
ARP2032		580mm (23")	(NC 2021, 2035, 2036, 2174, 2175, 2147)	○
ARP2033		640mm (24")	(NC 2021, 2035, 2036, 2174, 2175, 1404, 1405, 2147)	○
SEAT DEPTH*				
ARP2035		405mm (16")	(NC 2024, 2030, 2031, 2032, 2033, 2055)	○
ARP2036		430mm (17")	(NC 2024, 2030, 2031, 2032, 2033)	○
ARP2037		455mm (18")	(NC 2024)	○
ARP2038		480mm (19")		○
ARP2039		505mm (20")		○
ARP2040		530mm (21")		○
ARP2041		555mm (22")	(NC 2025, 2026)	○
ARP2042		580mm (23")	(NC 2025, 2026, 2146)	○

SEAT CUSHION*				
ARP2144		Matrx Vi Cushion (Mesh fabric)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	£400
ARP2146		Matrx Flo-tech Contour Visco	(NC 2042, 2030, 2031, 2033, 1985)	£169
ARP2147		Matrx Libra 90mm thick		£410
ARP2140		Without seat cushion	(NC 2033, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	◆
1st FIRST POSITIONING ACCESORIES FOR MATRX LIBRA AND FLO-TECH LITE				
ARP1980		Tapered Wedges (Pair)		£30
ARP1981		Abductor		£23
ARP1982		Lateral Wedge		£20
ARP1983		Anterior Seat Wedges		£27
2nd FIRST POSITIONING ACCESORIES FOR MATRX LIBRA AND FLO-TECH LITE				
ARP1980		Tapered Wedges (Pair)		£30
ARP1981		Abductor		£23
ARP1982		Lateral Wedge		£20
ARP1983		Anterior Seat Wedges		£27
POSITIONING ACCESSORIES FOR MATRX VI				
ARP1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)		£54
BACKREST SUPPORT*				
ARP2050		High Back, height 410mm (Mesh fabric)	(NC 2073, 2074)	◆
ARP2051		High Back height 460mm (Mesh fabric)	(NC 2070, 2071)	○
ARP2320		High Back, height 505mm / 20" (Mesh fabric)		○
ARP2550		Matrx Elite E2 Back, height 355mm / 14"	(NC 2028, 2029, 2030, 2032, 2033)	£460
ARP2551		Matrx Elite E2 Back, height 405mm / 15"		£460
ARP2552		Matrx Elite E2 Back, height 455mm / 18"		£460
ARP2553		Matrx Elite E2 Back, height 505mm / 20"		£460
ARP2554		Matrx Elite HD Back, height 405mm / 16"		£509
ARP2555		Matrx Elite HD Back, height 505mm / 20"		£509
ARP2556		Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)		£496
ARP2557		Matrx Elite TR Back, height 455mm / 18" (Mesh fabric)		£496
ARP2177		Matrx Elite TR Back height 505mm (Mesh fabric)		£496
ARP2563		Matrx Elite TR HD Back, height 505mm / 20" (Mesh fabric)		£550
ARP2558		Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)		£550
ARP2559		Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)		£550
ARP2178		Matrx Elite Deep Back, height 500mm (Mesh fabric)		£550
ARP2174		Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)		£496
ARP2175		Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)		£496
ARP2560		Matrx PB Deep Back, height 405mm / 16"		£509
ARP2561		Matrx PB Deep Back, height 505mm / 20"		£509
ARP2562		Matrx Elite Deep HD Back, height 505mm / 20"		£613
ARP2168		Matrx MAC Back - Foam Cushion [Mesh fabric]		£475
ARP2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]		£575

BACKREST HEIGHT*				
ARP2506		505mm (20") measured from seat pan to top of back shell	(NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562)	○
ARP2501		535mm (21") measured from seat pan to top of back shell	(NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562)	○
ARP2070		560mm (22") measured from seat pan to top of back shell	(NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563)	○
ARP2071		585mm (23") measured from seat pan to top of back shell	(NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563)	○
ARP2072		610mm (24") measured from seat pan to top of back shell	(NC 2050, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○
ARP2073		635mm (25") measured from seat pan to top of back shell	(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2558, 2559, 2174, 2560)	○
ARP2074		660mm (26") measured from seat pan to top of back shell	(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○
ARP2502		685mm (27") measured from seat pan to top of back shell	(NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)				
ARP2380		405 mm (16 inch) MAC BH		○
ARP2381		455 mm (18 inch) MAC BH		○
ARP2382		505 mm (20 inch) MAC BH		○
ARP2383		560 mm (22 inch) MAC BH		○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)				
ARP2385		90 mm (3.5 inch) MAC Contour Depth LH		£100
ARP2386		150 mm (6 inch) MAC Contour Depth LH		£100
ARP2387		230 mm (9 inch) MAC Contour Depth LH		£100
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)				
ARP2389		90 mm (3.5 inch) MAC Contour Depth RH		£100
ARP2390		150 mm (6 inch) MAC Contour Depth RH		£100
ARP2391		230 mm (9 inch) MAC Contour Depth RH		£100
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)				
ARP2393		Backrest -25 mm (-1 inch) to seat width		○
ARP2394		Backrest matches seat width		○
ARP2395		Backrest +25 mm (+1 inch) to seat width		○
ARMREST*				
ARP2095		Cantilever Flip-Back Arm, width, height and angle adjustable (230mm - 330mm) (9" - 13)		◆
ARP2096		Cantilever Flip-Back Arm, width, height and angle adjustable (320mm - 405mm) (12.5" - 16")		○
ARP2097		Two-post Recline Arm, Flip-Back, adjustable height (240mm - 330mm)		○
ARP0690		Without Armrest		○
ARMPAD LH*				
ARP2012		Desk arm pad with startex cover visco foam (255mm) 40m thick		◆
ARP2013		Full length armpad with startex cover foam (355mm) 40mm thick		○
ARP2014		Ergonomic arm though 30mm thick	(NC 1404, 1405, 2347, 2348, 2349, 2350)	£110
ARP2015		Without LH arm pad		○

ARMPAD RH*			
ARP2016		Desk arm pad with startex cover visco foam (255mm) 40mm thick	◆
ARP2017		Full length armpad with startex cover foam (355mm) 40mm thick	○
ARP2018		Ergonomic arm trough 30mm thick (NC 1404, 1405, 2351, 2352, 2353, 2354)	£110
ARP2019		Without RH arm pad	○
ARMPAD ACCESSORIES & MOUNTING OPTIONS LH			
ARP2080		Flat Hand Pad Left (in combination with arm trough only): to be mounted on the opposite side of remote	£110
ARMPAD ACCESSORIES & MOUNTING OPTIONS RH			
ARP2083		Flat Hand Pad Right (in combination with arm trough only): to be mounted on the opposite side of remote	£110
MULTI-AXIS UPPER EXTREMITY SUPPORTS LH			
ARP2082		Multi Axis Upper Extremity Support (MACES) Left, raises arm pad by 43mm	£270
MULTI-AXIS UPPER EXTREMITY SUPPORTS RH			
ARP2085		Multi Axis Upper Extremity Support (MACES) Right, raises arm pad by 43mm	£270
ARMREST HEIGHT SETTING LH*			
O001		Factory Setting	○
O002		Customised 230mm	○
O003		Customised 240mm	○
O004		Customised 250mm	○
O005		Customised 260mm	○
O006		Customised 270mm	○
O007		Customised 280mm	○
O008		Customised 290mm	○
O009		Customised 300mm	○
10		Customised 310mm	○
11		Customised 320mm	○
12		Customised 330mm	○
13		Customised 340mm	○
14		Customised 350mm	○
15		Customised 360mm	○
16		Customised 370mm	○
17		Customised 380mm	○
18		Customised 390mm	○
19		Customised 405mm	○

ARMREST HEIGHT SETTING RH*				
20		Factory Setting		○
21		Customised 230mm		○
22		Customised 240mm		○
23		Customised 250mm		○
24		Customised 260mm		○
25		Customised 270mm		○
25		Customised 280mm		○
27		Customised 290mm		○
28		Customised 300mm		○
29		Customised 310mm		○
30		Customised 320mm		○
31		Customised 330mm		○
32		Customised 340mm		○
33		Customised 350mm		○
34		Customised 360mm		○
35		Customised 370mm		○
36		Customised 380mm		○
37		Customised 390mm		○
38		Customised 405mm		○
ELBOW BLOCK LH PADS				
ARP2347		Elbow Block XS (110mm wide x 83mm high)	(NC 2014)	£475
ARP2348		Elbow Block S (135mm wide x 90mm high)	(NC 2014)	£475
ARP2349		Elbow Block M (145mm wide x 95mm high)	(NC 2014)	£475
ARP2350		Elbow Block L (165mm wide x 110mm high)	(NC 2014)	£475
ELBOW BLOCK LH MOUNTING HARDWARE				
ARP2361		Vertical 3"/ Horizontal 5"		○
ARP2362		Vertical 3"/ Horizontal 7"		○
ARP2363		Vertical 5"/ Horizontal 5"		○
ARP2364		Vertical 5"/ Horizontal 7"		○
ELBOW BLOCKS RH PADS				
ARP2351		Elbow Block XS (110mm wide x 83mm high)	(NC 2018)	£475
ARP2352		Elbow Block S (135mm wide x 90mm high)	(NC 2018)	£475
ARP2353		Elbow Block M (145mm wide x 95mm high)	(NC 2018)	£475
ARP2354		Elbow Block L (165mm wide x 110mm high)	(NC 2018)	£475
ELBOW BLOCK RH MOUNTING HARDWARE				
ARP2371		Vertical 3"/ Horizontal 5"		○
ARP2372		Vertical 3"/ Horizontal 7"		○
ARP2373		Vertical 5"/ Horizontal 5"		○
ARP2374		Vertical 5"/ Horizontal 7"		○

HIP SUPPORTS LH (Not available with 410mm seat depth and Manual Centre legrests)			
ARP2184		Hip Support LH, Quick Release, pad size 140 x 100 mm	£250
ARP2114		Hip Support LH, Quick Release, pad size 180 x 150 mm	£250
ARP2186		Hip Support LH, Quick Release, pad size 230 x 100 mm	£250
ARP2188		Hip Support LH, Quick Release, pad size 330 x 100 mm	£250
HIP SUPPORTS RH (Not available with 410mm seat depth and Manual Centre legrests)			
ARP2185		Hip Support RH, Quick release, pad size 140 x 100 mm	£250
ARP2115		Hip Support RH, Quick Release, pad size 180 x 150 mm	£250
ARP2187		Hip Support RH, Quick release, pad size 230 x 100 mm	£250
ARP2189		Hip Support RH, Quick Release, pad size 330 x 100 mm	£250
MOUNTING HARDWARE FOR HIP SUPPORTS			
ARP2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
HEADREST*			
ARP2092		Onyx Auto Style headrest (NC 2255, 2230, 2247)	£120
ARP2324		Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	◆
ARP2325		Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	£325
ARP1495		Head Control or equivalent Headrest, provided independently from this order (ACCORDING TO THE REQUIREMENTS OF ISO 7176-19, A HEADREST IS COMPULSORY FOR USING THE WHEELCHAIR AS VEHICLE SEAT)	○
LATERAL SUPPORTS LH			
ARP2343		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£395
ARP2345		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£395
ARP2043		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£395
ARP2045		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£395
LATERAL SUPPORTS RH			
ARP2344		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£395
ARP2346		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£395
ARP2044		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£395
ARP2046		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£395
LATERAL SUPPORT ACCESSORIES			
ARP2049		Additional telescopic link for Lateral support, 25mm, pair (NC 2558, 2559, 2178, 2560, 2561, 2562)	£160
CENTRE MOUNTED LEGREST			
ARP0586		Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 310mm - 410mm	£370
ARP0589		Hangers, centre mounted, length compensating, powered, LNX (0° - 90°): KNEE TO HEL LENGTH 340mm - 410mm (NC 2261, 2262, 2263, 2264)	£1,620
ARP0592		LNX Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380mm - 419mm (Selectable length depends on STF height) Custom-made, requires customised setting (NC 0559, 0516, 0517, 2056, 2263, 2264, 2365)	£2,450
FOOTPLATES			
ARP2365		Angle & width adjustable footplates Small (127 width x 190 depth / 5" x 7.5") (NC 0592)	£105
ARP0559		Angle & width adjustable footplate (240mm depth) (NC 0592)	○
ARP0516		Angle & width adjustable footplates (290mm depth) (NC 0592)	£135

FOOT PLATFORM				
ARP0517		Angle adjustable foot platform (240mm depth, 305mm width)	(NC 0592, 0586)	£270
SIDE MOUNTED HANGERS LH				
ARP0536		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm	(NC 0505)	◆
ARP0537		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	(NC 0505)	£120
ARP0539		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm	(NC 0505, 2261, 2262, 2263, 2264)	£330
ARP0535		Without Hanger and footplates	(NC 0505, 0509, 0511)	○
SIDE MOUNTED HANGERS RH				
ARP0546		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm	(NC 0502)	◆
ARP0547		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	(NC 0502)	£120
ARP0549		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm	(NC 0502, 2261, 2262, 2263, 2264)	£330
ARP0545		Without Hanger and footplates	(NC 0502, 0506, 0508)	○
INDIVIDUAL FOOTPLATE TYPE LH				
ARP0505		Fixed Footplate	(NC 0535, 0537, 0539)	◆
ARP0509		Angle Adjustable Footplate	(NC 0535)	£75
ARP0511		Angle and depth adjustable Footplate	(NC 0535)	○
INDIVIDUAL FOOTPLATE TYPE RH				
ARP0502		Fixed Footplate	(NC 0545, 0547, 0549)	◆
ARP0506		Angle Adjustable Footplate	(NC 0545)	£75
ARP0508		Angle and depth adjustable Footplate	(NC 0545)	○
CONTROLS*				
ARP2203		DLX-REM211 (Drive & Actuators)	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2209, 2210, 2213, 2214, 2075, 2076)	○
ARP2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2209, 2210, 2213, 2214, 2075, 2076)	◆
ARP2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2209, 2210, 2213, 2215, 2075, 2076)	£425
ARP2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1207, 1208, 1707, 1708, 1709, 1710, 2227, 2228, 2075, 2076)	£765
REM400 TOGGLE OPTION				
ARP2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204)	£315

CONTROL BRACKETRY*				
ARP1207		Standard Bracket RH	(NC 1404)	○
ARP1208		Standard Bracket LH	(NC 1405)	○
ARP1707		Swing Away RH Bracket	(NC 1404)	£135
ARP1708		Swing Away LH Bracket	(NC 1405)	£135
ARP1709		Matrx Resolve Swing in/out Bracket RH		£245
ARP1710		Matrx Resolve Swing in/out Bracket LH		£245
ARP2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 1404)	○
ARP2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 1405)	○
ARP2075		Swing Away + height / width adjust. Mount on RH for user, Quad link design		£225
ARP2076		Swing Away + height / width adjust. Mount on LH for user, Quad link design		£225
ARP2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 1404, 1405)	£495
ARP2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 1404, 1405)	£495
ARP2227		Bracket RH mounted for Attendant	(NC 1221, 2206, 2219)	£65
ARP2228		Bracket LH mounted for Attendant	(NC 1222, 2206, 2219)	£65
ARZ2239		Individual mounting option per request (Special). Requires additional information/specification.		POA
DUAL CONTROLS				
ARP2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1221, 1222, 2227, 2228, 2245)	£395
ARP2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2219, 2227, 2228)	£745
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 0080, 1221, 1222)	£2,315
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
ARP1207		Dual Control Standard Bracket RH	(NC 1221, 1222, 2227, 2228, 2245)	○
ARP1208		Dual Control Standard Bracket LH	(NC 1221, 1222, 2227, 2228, 2245)	○
ARP1221		IDC Dual Control Detachable Bracket RH	(NC 2219, 2227, 2228)	○
ARP1222		IDC Dual Control Detachable Bracket LH	(NC 2219, 2227, 2228)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA
COMPACT REMOTES				
ARP2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2235, 2236, 2241, 2242)	£825
ARP2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2235, 2236, 2241, 2242)	£1,075
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	£2,015
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0080, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242)	£3,190
ARP2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2240)	£4,100
ARP2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2233)	£2,285
ARP2243		ASL Paediatric (Proportional Joystick ASL132, incl. DLX-IN500)	(NC 2271, 2203, 2204, 2255, 2256)	£2,365

COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
ARP2207		RH Mount	(NC 2203, 2204)	○
ARP2208		LH Mount	(NC 2203, 2204)	○
ARP2213		Centre Mount Swing Away to RH (Nucleus)	(NC 2203, 2204)	£495
ARP2214		Centre Mount Swing Away to LH (Nucleus)	(NC 2203, 2204)	£495
ARP2235		Tray Mount RH (ASL603, 100mm x 150mm)		○
ARP2236		Tray Mount LH (ASL603, 100mm x 150mm)		○
ARP2213		Centre Mount Swing Away to RH (Nucleus)	(NC 1404, 1405)	£495
ARP2214		Centre Mount Swing Away to LH (Nucleus)	(NC 1404, 1405)	£495
ARP2241		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405)	£495
ARP2242		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405)	£495
ARZ2215		Individual mounting option per request (Special). Requires additional information/specification.	(NC 2203, 2204)	POA
CHIN CONTROLS				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£2,015
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)		£960
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)		£1,380
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)		£2,280
HEAD CONTROLS				
ARP2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2203, 2204)	£4,220
ARP2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2203, 2204)	£4,945
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2092, 2324, 2325, 2257, 2240, 2175, 0080)	£4,775
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2092, 2324, 2325, 2257, 2240, 2175, 0080)	£3,315
FOOT CONTROL				
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2203, 2204)	£2,880
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£2,740
FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)				
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification.		POA
TRAY CONTROL				
ARP2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2092, 2324, 2325, 2257, 2240)	£6,535
TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)				
ARP2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1404, 1405)	£495
ARP2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1404, 1405)	£495

SIP N PUFF HEAD CONTROL			
ARP2255	☐	Sip'n Puff Head Control (ASL109) includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements)	(NC 2271, 2203, 2204, 2092, 2324, 2325, 2243) £6,120
SIP N PUFF CONTROL & INPUT MODULE			
ARP2256	☐	Sip'n Puff Control (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module)	(NC 2271, 2203, 2204, 2243) £855
WIRELESS ACCESSORIES FOR ASL CONTROLS			
ARP2249	☐	Wireless Triple Switch Receiver (ASL557-3)	£340
ARP2250	☐	Wireless Mouse Emulator (ASL558)	£810
LiNX INTERFACE OPTIONS			
ARP2271	☐	Input Module DLX-IN500-A (sub D, Sip'n Puff)	(NC 2243, 2255, 2256, 2203, 2204) £785
ARP2272	☐	LiNX Output Module (OUT500)	£785
LiNX 10-WAY SWITCH			
ARP2251	☐	10-Way Switch DLX-SW10-A for direct Powered Seating Access	£500
LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
ARP2217	☐	Bracket RH, mounted for user	☐
ARP2218	☐	Bracket LH, mounted for user	☐
ARP2227	☐	For attendant, mounted on back shroud	£65
ARP2228	☐	Without mounting	£65
LiNX USB CHARGER PORT			
ARP2253	☐	USB Charger Port DLX-USB02-A	£160
LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
ARP2207	☐	RH Mount	☐
ARP2208	☐	LH Mount	☐
FURTHER LiNX COMPONENTS / ACCESSORIES			
ARP2302	☐	Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
ARP2300	☐	To be enclosed only	☐
BUDDY BUTTONS			
ARP1685	☐	Red Buddy Button	£90
ARP1686	☐	Yellow Buddy Button	£90
ARP1687	☐	Green Buddy Button	£90
ARP1690	☐	Blue Buddy Button	£90
EGG SWITCH			
ARP2254	☐	Egg Switch, Black	£90
PIKO BUTTONS			
ARP1711	☐	Red Piko Button	£160
ARP1712	☐	Yellow Piko Button	£160
ARP1713	☐	Green Piko Button	£160
ARP1714	☐	Orange Piko Button	£160
ARP1715	☐	Blue Piko Button	£160
ARP1716	☐	Black Piko Button	£160

TOGGLES / BUTTONS			
ARP2261		Stereo Toggle - direct control for Tilt only	£270
ARP2262		Stereo Button - direct control for Tilt only	£255
ARP2263		Four-Way Toggle - direct control for Tilt & Recline	£330
ARP2264		Four-Way Button - direct control for Tilt & Recline	£355
ARP2207		Toggle / Button RH Mount	○
ARP2208		Toggle / Button LH Mount	○
ARP2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
ARP2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£75
ARP2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£75
ARP2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
ARP2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
SHROUD COLOUR*			
ARP2640		Obsidian Black Frame Shrouds	○
ARP2641		Ruby Red Frame Shrouds	○
ARP2642		Purple Passion Frame Shrouds	○
ARP2643		Peacock Blue Frame Shrouds	○
ARP2644		Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARP2650		Obsidian Black Rim Inserts	○
ARP2651		Ruby Red Rim Inserts	○
ARP2652		Purple Passion Rim Inserts	○
ARP2653		Peacock Blue Rim Inserts	○
ARP2654		Arctic Blue Rim Inserts	○
FORKS*			
ARP2631		Single sided fork	£300
ARP2632		Double sided fork	○
CENTRE WHEELS (DRIVING WHEELS)*			
ARP0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARP0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	◆
MOTOR TYPE*			
ARP1294		3 km/h (2mph)	○
ARP1282		6 km/h (4mph)	○
ARP1283		10 km/h (6mph) (includes light system and indicators)	◆
BATTERY*			
ARP1257		50 Ah	◆
G TRAC			
ARP1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	◆
CHARGER*			
ARP1343		OptiCharge Charger 8A (200 - 250V)	◆
ARP1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£155

LIGHTS*			
ARP1437		LED Light System with Indicators (30 LUX) (NC 2203)	£450
ARX1006		Lights Included with Motor Package 10kmh (NC 1282)	◆
ARP1254		Without Light system	◆
TRANSPORTATION KIT*			
ARP1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
ARP1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter. (NC 2053)	◆
ARP1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory. (NC 2052)	○
FIXATION BELT*			
ARP2113		Lap belt with auto lock (flexible mount in side rails)	◆
ARP2701		Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	£95
ARP2702		Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	£95
ARP2703		2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	£90
ARP2704		2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	£90
ARP2705		2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	£90
ARP2706		2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	£90
TRAY			
ARP1404		Tray Swing Away RH (NC 1207, 1707, 2033, 2014, 2018, 2209, 2207, 2213, 2214, 2241, 2242, 2275, 2276)	£200
ARP1405		Tray Swing Away LH (NC 1208, 1708, 2033, 2014, 2018, 2210, 2208, 2213, 2214, 2241, 2242, 2275, 2276)	£200

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:	
ARP0080	TDX SP2 Ultra Narrow Base
ARP2052	Centre of gravity Tilt (0 to 50°)
ARP2053	300mm seat lifter (increased STF height by 4mm) with centre of gravity tilt (0 to 45°)
ARP0339	480mm
ARP2055	10° Fixed Anterior Assist Mount (will reduce max tilt angle by 10°)
ARP2056	5° Pre-Tilt (will increase max tilt angle by 5°)
ARP2021	405mm - 505mm (increments of 25mm) (16" - 20")
ARZ2003	480mm - 555mm (Increments of 25mm) (19" - 22") - Additional build time of 15+days
ARZ2022	505mm - 610mm (increments of 25mm) (20" - 24") - Additional build time of 15+days
ARP2024	480mm - 580mm (increments of 25mm) (19" - 23")
ARP2025	405mm (16")
ARP2026	430mm (17")
ARP2027	455mm (18")
ARP2028	480mm (19")
ARP2029	505mm (20")
ARP2030	530mm (21")
ARP2031	555mm (22")
ARP2032	580mm (23")
ARP2033	640mm (24")
ARP2035	405mm (16")
ARP2036	430mm (17")
ARP2037	455mm (18")
ARP2040	530mm (21")
ARP2041	555mm (22")
ARP2042	580mm (23")
ARP2146	Matrx Flo-tech Contour Visco
ARP2147	Matrx Libra 90mm thick
ARP1980	Tapered Wedges (Pair)
ARP1981	Abductor
ARP1982	Lateral Wedge
ARP1983	Anterior Seat Wedges
ARP1980	Tapered Wedges (Pair)
ARP1981	Abductor
ARP1982	Lateral Wedge
ARP1983	Anterior Seat Wedges
ARP1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
ARP2050	High Back, height 410mm (Mesh fabric)
ARP2051	High Back height 460mm (Mesh fabric)
ARP2320	High Back, height 505mm / 20" (Mesh fabric)
ARP2550	Matrx Elite E2 Back, height 355mm / 14"
ARP2551	Matrx Elite E2 Back, height 405mm / 15"
ARP2552	Matrx Elite E2 Back, height 455mm / 18"
ARP2553	Matrx Elite E2 Back, height 505mm / 20"

ARP2554	Matrx Elite HD Back, height 405mm / 16"
ARP2555	Matrx Elite HD Back, height 505mm / 20"
ARP2556	Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)
ARP2557	Matrx Elite TR Back, height 455mm / 18" (Mesh fabric)
ARP2177	Matrx Elite TR Back height 505mm (Mesh fabric)
ARP2563	Matrx Elite TR HD Back, height 505mm / 20" (Mesh fabric)
ARP2558	Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)
ARP2559	Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)
ARP2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)
ARP2174	Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)
ARP2175	Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)
ARP2560	Matrx PB Deep Back, height 405mm / 16"
ARP2561	Matrx PB Deep Back, height 505mm / 20"
ARP2562	Matrx Elite Deep HD Back, height 505mm / 20"
ARP2070	560mm (22") measured from seat pan to top of back shell
ARP2071	585mm (23") measured from seat pan to top of back shell
ARP2073	635mm (25") measured from seat pan to top of back shell
ARP2074	660mm (26") measured from seat pan to top of back shell
ARP2014	Ergonomic arm though 30mm thick
ARP2018	Ergonomic arm though 30mm thick
ARP2347	Elbow Block XS (110mm wide x 83mm high)
ARP2348	Elbow Block S (135mm wide x 90mm high)
ARP2349	Elbow Block M (145mm wide x 95mm high)
ARP2350	Elbow Block L (165mm wide x 110mm high)
ARP2351	Elbow Block XS (110mm wide x 83mm high)
ARP2352	Elbow Block S (135mm wide x 90mm high)
ARP2353	Elbow Block M (145mm wide x 95mm high)
ARP2354	Elbow Block L (165mm wide x 110mm high)
ARP2092	Onyx Auto Style headrest
ARP2324	Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW)
ARP2325	Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW)
ARP0586	Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 310mm - 410mm
ARP0592	LNx Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380mm - 419mm (Selectable length depends on STF height) Custom-made, requires customised setting
ARP2365	Angle & width adjustable footplates Small (127 width x 190 depth / 5" x 7.5")
ARP0559	Angle & width adjustable footplate (240mm depth)
ARP0516	Angle & width adjustable footplates (290mm depth)
ARP0517	Angle adjustable foot platform (240mm depth, 305mm width)
ARP0537	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARP0539	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm
ARP0535	Without Hanger and footplates
ARP0547	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARP0549	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm
ARP0545	Without Hanger and footplates
ARP0505	Fixed Footplate

ARP0509	Angle Adjustable Footplate
ARP0511	Angle and depth adjustable Footplate
ARP0502	Fixed Footplate
ARP0506	Angle Adjustable Footplate
ARP0508	Angle and depth adjustable Footplate
ARP2203	DLX-REM211 (Drive & Actuators)
ARP2204	DLX-REM216 (Drive, Actuators & Lights)
ARP2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
ARP1207	Standard Bracket RH
ARP1208	Standard Bracket LH
ARP1707	Swing Away RH Bracket
ARP1708	Swing Away LH Bracket
ARP1709	Matrx Resolve Swing in/out Bracket RH
ARP1710	Matrx Resolve Swing in/out Bracket LH
ARP2209	REM500 Swing Away RH Multi-Adjustable Bracket
ARP2210	REM500 Swing Away LH Multi-Adjustable Bracket
ARP2075	Swing Away + height / width adjust. Mount on RH for user, Quad link design
ARP2076	Swing Away + height / width adjust. Mount on LH for user, Quad link design
ARP2213	Centre Mount Swing Away to RH Nucleus Bracket
ARP2214	Centre Mount Swing Away to LH Nucleus Bracket
ARP2227	Bracket RH mounted for Attendant
ARP2228	Bracket LH mounted for Attendant
ARP2219	LiNX DLX-ACU200 (Compact Dual Control)
ARP2245	LiNX DLX-IDC (Intuitive Dual Control)
ARP1207	Dual Control Standard Bracket RH
ARP1208	Dual Control Standard Bracket LH
ARP1221	IDC Dual Control Detachable Bracket RH
ARP1222	IDC Dual Control Detachable Bracket LH
ARP2231	LiNX DLX-CR400 (Compact Remote)
ARP2232	LiNX DLX-CR400LF (Compact Remote Low Force)
ARP2233	ASL MEC (Proportional MEC Joystick ASL130)
ARP2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
ARP2243	ASL Paediatric (Proportional Joystick ASL132, incl. DLX-IN500)
ARP2207	RH Mount
ARP2208	LH Mount
ARP2213	Centre Mount Swing Away to RH (Nucleus)
ARP2214	Centre Mount Swing Away to LH (Nucleus)
ARP2235	Tray Mount RH (ASL603, 100mm x 150mm)
ARP2236	Tray Mount LH (ASL603, 100mm x 150mm)
ARP2213	Centre Mount Swing Away to RH (Nucleus)
ARP2214	Centre Mount Swing Away to LH (Nucleus)
ARP2241	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARP2242	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARZ2215	Individual mounting option per request (Special). Requires additional information/specification.

ARP2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
ARP2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
ARP2257	Four-Switch Proximity Array (ASL106)
ARP2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
ARP2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
ARP2255	Sip'n Puff Head Control (ASL109) includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements)
ARP2256	Sip'n Puff Control (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module)
ARP2271	Input Module DLX-IN500-A (sub D, Sip'n Puff)
ARP2227	For attendant, mounted on back shroud
ARP2228	Without mounting
ARP2207	RH Mount
ARP2208	LH Mount
ARP1686	Yellow Buddy Button
ARP1687	Green Buddy Button
ARP1690	Blue Buddy Button
ARP2254	Egg Switch, Black
ARP1711	Red Piko Button
ARP1712	Yellow Piko Button
ARP1713	Green Piko Button
ARP1714	Orange Piko Button
ARP1715	Blue Piko Button
ARP1716	Black Piko Button
ARP2261	Stereo Toggle - direct control for Tilt only
ARP2262	Stereo Button - direct control for Tilt only
ARP2263	Four-Way Toggle - direct control for Tilt & Recline
ARP2264	Four-Way Button - direct control for Tilt & Recline
ARP2207	Toggle / Button RH Mount
ARP2208	Toggle / Button LH Mount
ARP1282	6 km/h (4mph)
ARP1437	LED Light System with Indicators (30 LUX)
ARP1404	Tray Swing Away RH
ARP1405	Tray Swing Away LH



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):