



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2025
Customer Ref:

Customer Ref\*, Account No\*, Contact No\*, Order No\*, Delivery Address\*, Order Qty\*, Order date\*

TDX2 ULTRA NB

Standard Features

- Max user weight 136Kg (21.4 Stone)
Crash test ISO7176-Part 19
Ultra Low Maxx seating by Motion Concepts
Compact, manoeuvrable, CWD TDX SP2 chassis
Seat widths of 405mm - 610mm (16" - 24"), adjustable between 405mm - 505mm (16" - 20") and 505mm - 610mm (20" - 24")
Seat depth of 380mm - 580mm (15" - 23") adjustable between 380mm - 480mm (15" - 19") and 480mm - 580mm (19 - 23")
Backrest heights of 560mm - 660mm in 25mm increments
Cantilever flip-back arm, width, height and angle adjustable with desk armpad with startex covered visco foam
Elan multi positioning headrest
ULTRA Low Maxx CG power tilt (0° to 50°) seat heights of 420mm - 480mm
Powered Recline 90° to 170° with extended shear reduction (ESR)
Vari F hangers with fixed footplates
LiNX control REM211 with standard bracket
Selection of 10 stunning shroud and Rim colours
Centre (driving wheels) 355mm x 76mm (14" x 3") Black tyres
10km/h (6mph) 4 pole motors and 50 Ah batteries
G-Trac
4 Point tie down system
Lapbelt with autolock (adjustable both sides)



Q'STRAIT
QLK-150
Compatible

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code Description

BASE PRODUCT\*

O0001 TDX SP2 ULTRA NARROW BASE ○

CHAIR TYPE\*

ARP0080 TDX SP2 ULTRA Narrow Base ◆

ARX0080 Special TDX SP2 ULTRA Narrow Base (Please provide details in special section) POA

SEATING TYPE*			
ARP2020		Ultra Low Maxx Powered Positioning System	£11,729
SEAT ADJUSTMENTS*			
ARP2052		Centre of gravity Tilt (0 to 50° )	◆
ARP2053		300mm seat lifter (increased STF height by 4mm) with centre of gravity tilt (0 to 45° ) (NC 0339)	£2,138
SEAT HEIGHT SETTING*			
ARP0337		Seat-to -Floor Height 440mm	○
ARP0331		Seat-to Floor Height 460mm	○
ARP0339		Seat-to Floor Height 480mm (NC 2053)	○
ANTERIOR AND POSTERIOR TILT*			
ARP2057		Tilt Standard Mount (No anterior or posterior tilt)	◆
ARP2054		5° Fixed Anterior Assist Mount (will reduce max tilt angle by 5°)	£290
ARP2055		10° Fixed Anterior Assist Mount (will reduce max tilt angle by 10°) (NC 2035)	£361
ARP2056		5° Pre-Tilt (will increase max tilt angle by 5°)	£290
BACKREST RECLINER*			
ARP2060		Powered Recliner (90° - 170°) with Extended Shear Reduction (ESR), Standard Mount	◆
ARP2063		Powered Recliner 8° Recline mount 82° to 162° with Extended Shear Reduction (ESR)	£30
ARP2066		Powered Recline, 30° Recline mount, 60° to 140° with Extended Shear Reduction (ESR)	£30
SEAT FRAME WIDTH RANGE*			
ARP2021		405mm - 505mm (increments of 25mm) (16" - 20") (NC 2030, 2031, 2032, 2033)	○
ARX2003		480mm - 555mm (Increments of 25mm) (19" - 22") - Additional build time of 15+days (NC 2025, 2026, 2027, 2032, 2033)	○
ARX2022		505mm - 610mm (increments of 25mm) (20" - 24") - Additional build time of 15+days (NC 2025, 2026, 2027, 2174, 2175)	○
SEAT FRAME DEPTH RANGE*			
ARP2024		480mm - 580mm (increments of 25mm) (19" - 23") (NC 2035, 2036, 2037)	○
ARP2008		405mm - 505mm (increments of 25mm) (16" - 20") (NC 2040, 2041, 2042)	○
SEAT WIDTH*			
ARP2025		405mm (16") (NC 2003, 2022, 2041, 2042)	○
ARP2026		430mm (17") (NC 2003, 2022, 2041, 2042)	○
ARP2027		455mm (18") (NC 2003, 2022)	○
ARP2028		480mm (19") (NC 2022, 2174, 2175)	○
ARP2029		505mm (20") (NC 2174, 2175)	○
ARP2030		530mm (21") (NC 2021, 2035, 2036, 2174, 2175)	○
ARP2031		555mm (22") (NC 2021, 2035, 2036, 2174, 2175)	○
ARP2032		580mm (23") (NC 2021, 2035, 2036, 2174, 2175, 2147)	○
ARP2033		640mm (24") (NC 2021, 2035, 2036, 2174, 2175, 1404, 1405, 2147)	○
SEAT DEPTH SETTING*			
ARP2035		405mm (16") (NC 2024, 2030, 2031, 2032, 2033, 2055)	○
ARP2036		430mm (17") (NC 2024, 2030, 2031, 2032, 2033)	○
ARP2037		455mm (18") (NC 2024)	○
ARP2038		480mm (19")	○
ARP2039		505mm (20")	○
ARP2040		530mm (21")	○
ARP2041		555mm (22") (NC 2025, 2026)	○
ARP2042		580mm (23") (NC 2025, 2026, 2146)	○

SEAT CUSHION*			
ARP2144		Matrx Vi Cushion (Mesh fabric)	£410
ARP2146		Matrx Flo-tech Contour Visco (NC 2041, 2042, 2030, 2031, 2033)	£112
ARP2147		Matrx Libra 90mm thick (NC 2042)	£391
ARP2140		Solid seat Base, For alternative cushion options (NC 2032, 2033)	◆
BACKREST SUPPORT*			
ARP2050		High Back, height 410mm (Mesh fabric) (NC 2073, 2074)	◆
ARP2051		High Back height 460mm (Mesh fabric) (NC 2070, 2071)	○
ARP2320		High Back, height 505mm / 20" (Mesh fabric)	○
ARP2550		Matrx Elite E2 Back, height 355mm / 14" (NC 2028, 2029, 2030, 2032, 2033)	£460
ARP2551		Matrx Elite E2 Back, height 405mm / 15"	£460
ARP2552		Matrx Elite E2 Back, height 455mm / 18"	£460
ARP2553		Matrx Elite E2 Back, height 505mm / 20"	£460
ARP2554		Matrx Elite HD Back, height 405mm / 16"	£509
ARP2555		Matrx Elite HD Back, height 505mm / 20"	£509
ARP2556		Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)	£496
ARP2557		Matrx Elite TR Back, height 455mm / 18" (Mesh fabric)	£496
ARP2177		Matrx Elite TR Back height 505mm (Mesh fabric)	£496
ARP2563		Matrx Elite TR HD Back, height 505mm / 20" (Mesh fabric)	£550
ARP2558		Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)	£550
ARP2559		Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)	£550
ARP2178		Matrx Elite Deep Back, height 500mm (Mesh fabric)	£550
ARP2174		Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)	£496
ARP2175		Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)	£496
ARP2560		Matrx PB Deep Back, height 405mm / 16"	£509
ARP2561		Matrx PB Deep Back, height 505mm / 20"	£509
ARP2562		Matrx Elite Deep HD Back, height 505mm / 20"	£613
BACKREST HEIGHT*			
ARP2506		505mm (20") measured from seat pan to top of back shell (NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562)	○
ARP2501		535mm (21") measured from seat pan to top of back shell (NC 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562)	○
ARP2070		560mm (22") measured from seat pan to top of back shell (NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563)	○
ARP2071		585mm (23") measured from seat pan to top of back shell (NC 2050, 2551, 2553, 2554, 2555, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563)	○
ARP2072		610mm (24") measured from seat pan to top of back shell (NC 2050, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○
ARP2073		635mm (25") measured from seat pan to top of back shell (NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2558, 2559, 2174, 2560)	○
ARP2074		660mm (26") measured from seat pan to top of back shell (NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○
ARP2502		685mm (27") measured from seat pan to top of back shell (NC 2050, 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560)	○

<b>ARMREST*</b>			
ARP2095		Cantilever Flip-Back Arm, width, height and angle adjustable (230mm - 330mm) (9" - 13)	◆
ARP2096		Cantilever Flip-Back Arm, width, height and angle adjustable (320mm - 405mm) (12.5" - 16")	○
ARP2097		Two-post Recline Arm, Flip-Back, adjustable height (240mm - 330mm)	○
ARP0690		Without Armrest	○
<b>ARMPAD LH*</b>			
ARP2012		Desk arm pad with startex cover visco foam (255mm) 40m thick	◆
ARP2013		Full length armpad with startex cover foam (355mm) 40mm thick	○
ARP2014		Ergonomic arm trough 30mm thick (NC 1404, 1405, 2347, 2348, 2349, 2350)	£98
ARP2015		Without LH arm pad	○
<b>ARMPAD RH*</b>			
ARP2016		Desk arm pad with startex cover visco foam (255mm) 40mm thick	◆
ARP2017		Full length armpad with startex cover foam (355mm) 40mm thick	○
ARP2018		Ergonomic arm trough 30mm thick (NC 1404, 1405, 2351, 2352, 2353, 2354)	£98
ARP2019		Without RH arm pad	○
<b>ARMPAD ACCESSORIES &amp; MOUNTING OPTIONS LH</b>			
ARP2080		Flat Hand Pad Left (in combination with arm trough only): to be mounted on the opposite side of remote	£98
<b>ARMPAD ACCESSORIES &amp; MOUNTING OPTIONS RH</b>			
ARP2083		Flat Hand Pad Right (in combination with arm trough only): to be mounted on the opposite side of remote	£98
<b>MULTI-AXIX UPPER EXTREMITY SUPPORTS LH</b>			
ARP2082		Multi Axis Upper Extremity Support (MACES) Left, raises arm pad by 43mm	£245
<b>MULTI-AXIX UPPER EXTREMITY SUPPORTS RH</b>			
ARP2085		Multi Axis Upper Extremity Support (MACES) Right, raises arm pad by 43mm	£245
<b>ELBOW BLOCK LH PADS</b>			
ARP2347		Elbow Block XS (110mm wide x 83mm high) (NC 2014)	£459
ARP2348		Elbow Block S (135mm wide x 90mm high) (NC 2014)	£459
ARP2349		Elbow Block M (145mm wide x 95mm high) (NC 2014)	£459
ARP2350		Elbow Block L (165mm wide x 110mm high) (NC 2014)	£459
<b>ELBOW BLOCK LH MOUNTING HARDWARE</b>			
ARP2361		Vertical 3"/ Horizontal 5"	○
ARP2362		Vertical 3"/ Horizontal 7"	○
ARP2363		Vertical 5"/ Horizontal 5"	○
ARP2364		Vertical 5"/ Horizontal 7"	○
<b>ELBOW BLOCKS RH PADS</b>			
ARP2351		Elbow Block XS (110mm wide x 83mm high) (NC 2018)	£459
ARP2352		Elbow Block S (135mm wide x 90mm high) (NC 2018)	£459
ARP2353		Elbow Block M (145mm wide x 95mm high) (NC 2018)	£459
ARP2354		Elbow Block L (165mm wide x 110mm high) (NC 2018)	£459
<b>ELBOW BLOCK RH MOUNTING HARDWARE</b>			
ARP2371		Vertical 3"/ Horizontal 5"	○
ARP2372		Vertical 3"/ Horizontal 7"	○
ARP2373		Vertical 5"/ Horizontal 5"	○
ARP2374		Vertical 5"/ Horizontal 7"	○

<b>HIP SUPPORTS LH (Not available with 410mm seat depth and Manual Centre legrests)</b>			
ARP2184		Hip Support LH, pad size 140 x 100 mm	£226
ARP2114		Hip Support LH, pad size 180 x 150 mm	£226
ARP2186		Hip Support LH, pad size 230 x 100 mm	£226
ARP2188		Hip Support LH, pad size 330 x 100 mm	£226
<b>HIP SUPPORTS RH (Not available with 410mm seat depth and Manual Centre legrests)</b>			
ARP2185		Hip Support RH, pad size 140 x 100 mm	£226
ARP2115		Hip Support RH, pad size 180 x 150 mm	£226
ARP2187		Hip Support RH, pad size 230 x 100 mm	£226
ARP2189		Hip Support RH, pad size 330 x 100 mm	£226
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
ARP2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
<b>HEADREST*</b>			
ARP2092		Onyx Auto Style headrest (NC 2255, 2230, 2247)	£109
ARP2324		Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	◆
ARP2325		Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW) (NC 2255, 2230, 2247)	£294
ARP1495		Head Control or equivalent Headrest, provided independently from this order (ACCORDING TO THE REQUIREMENTS OF ISO 7176-19, A HEADREST IS COMPULSORY FOR USING THE WHEELCHAIR AS VEHICLE SEAT)	○
<b>LATERAL SUPPORTS LH</b>			
ARP2343		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£355
ARP2345		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£355
ARP2043		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£355
ARP2045		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£355
<b>LATERAL SUPPORTS RH</b>			
ARP2344		110mm wide x 83mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£355
ARP2346		135mm wide x 90mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£355
ARP2044		145mm wide x 95mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£355
ARP2046		165mm wide x 110mm (NC 2558, 2559, 2178, 2560, 2561, 2562)	£355
<b>LATERAL SUPPORT ACCESSORIES</b>			
ARP2049		Additional telescopic link for Lateral support, 25mm, pair (NC 2558, 2559, 2178, 2560, 2561, 2562)	£141
<b>CENTRE MOUNTED LEGREST</b>			
ARP0586		Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 310mm - 410mm	£286
ARP0589		Hangers, centre mounted, length compensating, powered, LNX (0° - 90°): KNEE TO HEL LENGTH 340mm - 410mm (NC 2261, 2262, 2263, 2264)	£1,588
ARP0592		LNX Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380mm - 419mm (Selectable length depends on STF height) Custom-made, requires customised setting (NC 0559, 0516, 0517, 2056, 2263, 2264, 2365)	£2,417
<b>FOOTPLATES</b>			
ARP2365		Angle & width adjustable footplates Small (127 width x 190 depth / 5" x 7.5") (NC 0592)	£93
ARP0559		Angle & width adjustable footplate (240mm depth) (NC 0592)	○
ARP0516		Angle & width adjustable footplates (290mm depth) (NC 0592)	£121
<b>FOOT PLATFORM</b>			
ARP0517		Angle adjustable foot platform (240mm depth, 305mm width) (NC 0592, 0586)	£243

SIDE MOUNTED HANGERS LH				
ARP0536		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm	(NC 0505)	◆
ARP0537		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	(NC 0505)	£105
ARP0539		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm	(NC 0505, 2261, 2262, 2263, 2264)	£299
ARP0535		Without Hanger and footplates	(NC 0505, 0509, 0511)	○
SIDE MOUNTED HANGERS RH				
ARP0546		Vari F: Angle adjustable with tools 0° - 70°, length 290mm - 460mm	(NC 0502)	◆
ARP0547		Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	(NC 0502)	£105
ARP0549		ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm	(NC 0502, 2261, 2262, 2263, 2264)	£299
ARP0545		Without Hanger and footplates	(NC 0502, 0506, 0508)	○
INDIVIDUAL FOOTPLATE TYPE LH				
ARP0505		Fixed Footplate	(NC 0535, 0537, 0539)	◆
ARP0509		Angle Adjustable Footplate	(NC 0535)	£64
ARP0511		Angle and depth adjustable Footplate	(NC 0535)	○
INDIVIDUAL FOOTPLATE TYPE RH				
ARP0502		Fixed Footplate	(NC 0545, 0547, 0549)	◆
ARP0506		Angle Adjustable Footplate	(NC 0545)	£64
ARP0508		Angle and depth adjustable Footplate	(NC 0545)	○
CONTROLS*				
ARP2203		REM211 Drive & Actuators	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2209, 2210, 2213, 2214, 2075, 2076)	○
ARP2204		REM216 Drive, Actuators & Lights	(NC 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2209, 2210, 2213, 2214, 2075, 2076)	◆
ARP2205		REM400, touch-screen Interface, built in Bluetooth HID	(NC 2209, 2210, 2213, 2215, 2075, 2076)	£305
ARP2206		REM500, without joystick for alternative access, built in Bluetooth HID	(NC 1207, 1208, 1707, 1708, 1709, 1710, 2227, 2228, 2075, 2076)	£697
CONTROL BRACKETRY*				
ARP1207		Standard Bracket RH	(NC 1404)	○
ARP1208		Standard Bracket LH	(NC 1405)	○
ARP1707		Swing Away + height adjustable Bracket RH High end	(NC 1404)	£119
ARP1708		Swing Away + height adjustable Bracket LH High end	(NC 1405)	£119
ARP1709		Matrx Resolve Swing in/out bracket RH		£219
ARP1710		Matrx Resolve Swing in/out bracket LH		£219
ARP2075		Swing Away + height / width adjust mount on RH for user, Quad link design		£202
ARP2076		Swing Away + height / width adjust mount on LH for user, Quad link design		£202
ARP2227		Centre Mount swing-away Bracket to RH (Nucleus)	(NC 2219, 2245)	£57
ARP2228		Centre Mount swing-away Bracket to LH (Nucleus)	(NC 2219, 2245)	£57
ARP2209		REM500 Swing away RH multi-adjustable	(NC 1404)	○
ARP2210		REM500 Swing away LH multi-adjustable	(NC 1405)	○
ARP2213		REM500 Centre mount swing away to RH Nucleus	(NC 1404, 1405)	£446
ARP2214		REM500 Centre mount swing away to LH Nucleus	(NC 1404, 1405)	£446
ARP2075		Swing Away + height / width adjustable Bracket Mount on RH, Quad link design		£202
ARP2076		Swing Away + height / width adjustable Bracket Mount on RH Quad link design		£202
ARX2239		Individual mounting option per (special) - Requires additional information / specification		POA

<b>REM400 TOGGLE OPTION</b>			
ARP2270		Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)	(NC 2203, 2204) £284
<b>CHIN CONTROL</b>			
UNI9250		Mo-vis Micro joystick (Please specify in Special Comments box mounting request if required)	(NC 0080) POA
UNI9251		Mo-vis Multi joystick (Please specify in Special Comments box mounting request if required)	(NC 0080) POA
<b>DUAL CONTROL</b>			
ARP2219		Compact Dual Control (DLX-ACU200) (Please select Dual control bracket below)	(NC 1221, 1222, 2227, 2228, 2245) £355
<b>DUAL CONTROL BRACKETRY</b>			
ARP1207		Dual Control Standard Bracket RH	(NC 1221, 1222, 2227, 2228, 2245) ○
ARP1208		Dual Control Standard Bracket LH	(NC 1221, 1222, 2227, 2228, 2245) ○
<b>INVACARE LINX IDC DUAL CONTROL</b>			
ARP2245		IDC Dual control (Please select AQE1221 OR AQE1222 for bracket option below)	(NC 2219, 2227, 2228) £673
<b>INVACARE LINX IDC DUAL CONTROL BRACKETRY</b>			
ARP1221		IDC Dual control Detachable Bracket RH	(NC 2219, 2227, 2228) ○
ARP1222		IDC Dual control Detachable Bracket LH	(NC 2219, 2227, 2228) ○
<b>LINX COMPACT REMOTE</b>			
ARP2231		Proportional Compact Remote DLX-CR400 (Please select Proportional Bracket)	(NC 2203, 2204) £740
ARP2232		Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)	(NC 2203, 2204) £990
<b>LINX COMPACT REMOTE BRACKET</b>			
ARP2207		Proportional Compact RH Mount Bracket	(NC 2203, 2204) ○
ARP2208		Proportional Compact LH Mount Bracket	(NC 2203, 2204) ○
ARP2213		Proportional Compact Centre Mount swing away to RH (Nucleus)	(NC 2203, 2204) £446
ARP2214		Proportional Compact Centre Mount swing away to LH (Nucleus)	(NC 2203, 2204) £446
ARX2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2203, 2204) POA
<b>LINX HEAD CONTROL</b>			
ARP2230		Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3	(NC 2203, 2204) £3,835
ARP2247		Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings	(NC 2203, 2204) £4,491
<b>LINX MICRO EXTREMITY CONTROL</b>			
ARP2233		MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry)	(NC 2203, 2204, 2230, 2247, 2257, 2240) £3,723
<b>LINX MICRO EXTREMITY CONTROL BRACKET</b>			
ARP2235		MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please also select MEC Bracketry)	○
ARP2236		MEC Tray Mount LH (ASL603, 100mm x 150mm)	○
ARP2213		Centre Mount, swing away to RH, Nucleus	(NC 1404, 1405) £446
ARP2214		Centre Mount, swing away to LH, Nucleus	(NC 1404, 1405) £446
ARP2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405) £446
ARP2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1404, 1405) £446
<b>LINX HEAVY DUTY REMOTE</b>			
ARP2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LiNX Heavy Duty Bracketry)	(NC 2203, 2204, 2230, 2247, 2257, 2233) £2,076

<b>LINX HEAVY DUTY REMOTE BRACKET</b>			
ARP2235		Tray Mount RH, (ASL603, 100mm x 150mm)	○
ARP2236		Tray Mount LH, (ASL603, 100mm x 150mm)	○
ARP2213		Centre Mount, swing away to RH, Nucleus (NC 1404, 1405)	£446
ARP2214		Centre Mount, swing away to LH, Nucleus (NC 1404, 1405)	£446
ARP2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm) (NC 1404, 1405)	£446
ARP2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm) (NC 1404, 1405)	£446
ARX2215		Individual mounting option per request (Special) Requires additional information / specification	POA
<b>LINX 10-WAY SWITCH</b>			
ARP2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch) (Please select 10 way switch bracket below)	£453
<b>LINX 10-WAY SWITCH BRACKET</b>			
ARP2217		10-Way Switch Bracket RH, Mounted for User	○
ARP2218		10-Way Switch Bracket LH, Mounted for User	○
ARP2227		10-Way Switch Bracket RH, Mounted for Attendant	£57
ARP2228		10-Way Switch Bracket LH, Mounted for Attendant	£57
<b>LINX FOUR -SWITCH PROXIMITY ARRAY</b>			
ARP2257		Four-Switch Proximity Array (ASL106) NEW! (NC 2203, 2204, 2230, 2092, 2324, 2325, 2257, 2240)	£5,938
ARP2275		Centre Mount, swing-away to RH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm) (NC 1404, 1405)	£446
ARP2276		Centre Mount, swing-away to LH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm) (NC 1404, 1405)	£446
<b>LINX WIRELESS ACCESSORIES FOR ATOM</b>			
ARP2249		Wireless Triple Switch Receiver (ASL557-3) (Sub D plug) NEW!	£305
ARP2250		Wireless Mouse Emulator (ASL558) NEW!	£733
<b>LINX PEDIATRIC COMPACT</b>			
ARP2243		Proportional Paediatric Compact (ASK132), Includes input Module NEW! (NC 2271, 2203, 2204, 2255, 2256)	£2,146
ARP2235		Tray Mount RH (ASL603, 100mm x 150mm)	○
ARP2236		Tray Mount LH (ASL603, 100mm x 150mm)	○
ARP2213		Centre Mount, swing-away to RH, Nucleus (NC 1404, 1405)	£446
ARP2214		Centre Mount, swing-away to LH, Nucleus (NC 1404, 1405)	£446
ARP2241		Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm) (NC 1404, 1405)	£446
ARP2242		Centre Mount, swing-away to LH, Nucleus (incl Tray (ASL603, 100mm x 150mm) (NC 1404, 1405)	£446
<b>LINX SIP N PUFF HEAD CONTROL</b>			
ARP2255		Sip'n Puff Head Control (ASL109) NEW! includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements) (NC 2271, 2203, 2204, 2092, 2324, 2325, 2243)	£5,562
<b>LINX SIP N PUFF CONTROL &amp; INPUT MODULE</b>			
ARP2256		Sip'n Puff Control NEW! (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module) (NC 2271, 2203, 2204, 2243)	£773
<b>LINX INTERFACE OPTIONS</b>			
ARP2271		Input Module DLX-IN500-A (sub D, Sip'n Puff) NEW! (NC 2243, 2255, 2256, 2203, 2204)	£785
<b>LINX USB CHARGER PORT</b>			
ARP2253		USB Charger Port DLX-USB02-A	£143



<b>LINX USB CHARGER PORT BRACKET</b>			
ARP2207	<input type="checkbox"/>	USB Charger Port RH Mount	<input type="radio"/>
ARP2208	<input type="checkbox"/>	USB Charger Port LH Mount	<input type="radio"/>
<b>OUTPUT MODULE</b>			
ARP2272	<input type="checkbox"/>	LiNX Output Module (OUT500)	£785
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
ARP2300	<input type="checkbox"/>	To be enclosed only	<input type="radio"/>
<b>BUDDY BUTTONS</b>			
ARP1685	<input type="checkbox"/>	Red Buddy Button	£79
ARP1686	<input type="checkbox"/>	Yellow Buddy Button	£79
ARP1687	<input type="checkbox"/>	Green Buddy Button	£79
ARP1690	<input type="checkbox"/>	Blue Buddy Button	£79
<b>EGG SWITCH</b>			
ARP2254	<input type="checkbox"/>	Egg Switch, Black	£79
<b>PIKO BUTTONS</b>			
ARP1711	<input type="checkbox"/>	Red Piko Button	£142
ARP1712	<input type="checkbox"/>	Yellow Piko Button	£142
ARP1713	<input type="checkbox"/>	Green Piko Button	£142
ARP1714	<input type="checkbox"/>	Orange Piko Button	£142
ARP1715	<input type="checkbox"/>	Blue Piko Button	£142
ARP1716	<input type="checkbox"/>	Black Piko Button	£142
<b>TOGGLES / BUTTONS</b>			
ARP2261	<input type="checkbox"/>	Stereo Toggle - direct control for Tilt only	£242
ARP2262	<input type="checkbox"/>	Stereo Button - direct control for Tilt only	£230
ARP2263	<input type="checkbox"/>	Four-Way Toggle - direct control for Tilt & Recline	£300
ARP2264	<input type="checkbox"/>	Four-Way Button - direct control for Tilt & Recline	£321
ARP2207	<input type="checkbox"/>	Toggle / Button RH Mount	<input type="radio"/>
ARP2208	<input type="checkbox"/>	Toggle / Button LH Mount	<input type="radio"/>
ARP2225	<input type="checkbox"/>	Without Mounting	<input type="radio"/>
<b>BODYPOINT JOYSTICK KNOBS</b>			
ARP2721	<input type="checkbox"/>	3" U-Shaped, Joystick Handle, Bodypoint PC101A	£65
ARP2722	<input type="checkbox"/>	4" U-Shaped, Joystick Handle, Bodypoint PC102A	£65
ARP2723	<input type="checkbox"/>	Domed Rubber, Joystick Handle, Bodypoint PC107A	£65
ARP2724	<input type="checkbox"/>	Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£65
<b>REMOTE STOP SWITCH</b>			
ARP2258	<input type="checkbox"/>	Remote Stop Switch (ASL504)	£942
<b>SHROUD COLOUR*</b>			
ARP2640	<input type="checkbox"/>	Obsidian Black Frame Shrouds	<input type="radio"/>
ARP2641	<input type="checkbox"/>	Ruby Red Frame Shrouds	<input type="radio"/>
ARP2642	<input type="checkbox"/>	Purple Passion Frame Shrouds	<input type="radio"/>
ARP2643	<input type="checkbox"/>	Peacock Blue Frame Shrouds	<input type="radio"/>
ARP2644	<input type="checkbox"/>	Arctic Blue Frame Shrouds	<input type="radio"/>

<b>RIM COLOUR*</b>			
ARP2650		Obsidian Black Rim Inserts	○
ARP2651		Ruby Red Rim Inserts	○
ARP2652		Purple Passion Rim Inserts	○
ARP2653		Peacock Blue Rim Inserts	○
ARP2654		Arctic Blue Rim Inserts	○
<b>FORKS*</b>			
ARP2631		Single sided fork (New) includes Black castor	£316
ARP2632		Double sided fork	○
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>			
ARP0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARP0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	◆
<b>MOTOR TYPE*</b>			
ARP1294		3 km/h (2mph)	○
ARP1282		6 km/h (4mph)	○
ARP1285		10 km/h (6mph) (includes light system and indicators)	◆
<b>BATTERY*</b>			
ARP1257		50 Ah	◆
<b>PERFORMANCE ACCESSORIES</b>			
ARP1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	◆
<b>CHARGER*</b>			
ARP1343		OptiCharge Charger 8A (200 - 250V)	◆
ARP1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£137
<b>LIGHTS*</b>			
ARP1437		LED Light System with Indicators (30 LUX) (NC 2203)	£431
ARZ1006		Lights Included with Motor Package 10 & 12kmh (NC 1282)	◆
ARP1254		Without Light system	◆
<b>TRANSPORTATION KIT*</b>			
ARP1450		Preparation for DAHL Docking Station: Docking station for car implementation, please order docking station separately. (Adequate headrest is compulsory).	£524
ARP1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) (NC 2053)	◆
ARP1401		6-Point-tie--down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19-2001) (Adequate headrest is compulsory) (NC 2052)	○
<b>FIXATION BELT*</b>			
ARP2113		Lap belt with auto lock (flexible mount in side rails)	◆

<b>BODYPOINT PELVIC BELTS (includes mounting kit)</b>			
ARP2701		Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	£84
ARP2702		Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	£84
ARP2703		2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	£78
ARP2704		2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	£78
ARP2705		2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	£78
ARP2706		2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	£78
<b>TRAY</b>			
ARP1404		Tray Swing Away RH	(NC 1207, 1707, 2033, 2014, 2018, 2209, 2207, 2213, 2214, 2241, 2242, 2275, 2276) £174
ARP1405		Tray Swing Away LH	(NC 1208, 1708, 2033, 2014, 2018, 2210, 2208, 2213, 2214, 2241, 2242, 2275, 2276) £174

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ARP0080	TDX SP2 ULTRA Narrow Base
ARX0080	Special TDX SP2 ULTRA Narrow Base (Please provide details in special section)
ARP2052	Centre of gravity Tilt (0 to 50° )
ARP2053	300mm seat lifter (increased STF height by 4mm) with centre of gravity tilt (0 to 45° )
ARP0339	Seat-to Floor Height 480mm
ARP2055	10° Fixed Anterior Assist Mount (will reduce max tilt angle by 10°)
ARP2056	5° Pre-Tilt (will increase max tilt angle by 5°)
ARP2021	405mm - 505mm (increments of 25mm) (16" - 20")
ARX2003	480mm - 555mm (Increments of 25mm) (19" - 22") - Additional build time of 15+days
ARX2022	505mm - 610mm (increments of 25mm) (20" - 24") - Additional build time of 15+days
ARP2024	480mm - 580mm (increments of 25mm) (19" - 23")
ARP2025	405mm (16")
ARP2026	430mm (17")
ARP2027	455mm (18")
ARP2028	480mm (19")
ARP2029	505mm (20")
ARP2030	530mm (21")
ARP2031	555mm (22")
ARP2032	580mm (23")
ARP2033	640mm (24")
ARP2035	405mm (16")
ARP2036	430mm (17")
ARP2037	455mm (18")
ARP2040	530mm (21")
ARP2041	555mm (22")
ARP2042	580mm (23")
ARP2146	Matrx Flo-tech Contour Visco
ARP2147	Matrx Libra 90mm thick
ARP2050	High Back, height 410mm (Mesh fabric)
ARP2051	High Back height 460mm (Mesh fabric)
ARP2320	High Back, height 505mm / 20" (Mesh fabric)
ARP2550	Matrx Elite E2 Back, height 355mm / 14"
ARP2551	Matrx Elite E2 Back, height 405mm / 15"
ARP2552	Matrx Elite E2 Back, height 455mm / 18"
ARP2553	Matrx Elite E2 Back, height 505mm / 20"
ARP2554	Matrx Elite HD Back, height 405mm / 16"
ARP2555	Matrx Elite HD Back, height 505mm / 20"
ARP2556	Matrx Elite TR Back, height 405mm / 16" (Mesh fabric)
ARP2557	Matrx Elite TR Back, height 455mm / 18" (Mesh fabric)
ARP2177	Matrx Elite TR Back height 505mm (Mesh fabric)
ARP2563	Matrx Elite TR HD Back, height 505mm / 20" (Mesh fabric)
ARP2558	Matrx Elite Deep Back, height 405mm / 16" (Mesh fabric)
ARP2559	Matrx Elite Deep Back, height 455mm / 18" (Mesh fabric)

ARP2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)
ARP2174	Matrx Posture Back, height 405mm (Mesh fabric max dimension 455mm wide)
ARP2175	Matrx Posture Back height 505mm (Mesh fabric max dimension 455mm wide)
ARP2560	Matrx PB Deep Back, height 405mm / 16"
ARP2561	Matrx PB Deep Back, height 505mm / 20"
ARP2562	Matrx Elite Deep HD Back, height 505mm / 20"
ARP2070	560mm (22") measured from seat pan to top of back shell
ARP2071	585mm (23") measured from seat pan to top of back shell
ARP2073	635mm (25") measured from seat pan to top of back shell
ARP2074	660mm (26") measured from seat pan to top of back shell
ARP2014	Ergonomic arm though 30mm thick
ARP2018	Ergonomic arm though 30mm thick
ARP2347	Elbow Block XS (110mm wide x 83mm high)
ARP2348	Elbow Block S (135mm wide x 90mm high)
ARP2349	Elbow Block M (145mm wide x 95mm high)
ARP2350	Elbow Block L (165mm wide x 110mm high)
ARP2351	Elbow Block XS (110mm wide x 83mm high)
ARP2352	Elbow Block S (135mm wide x 90mm high)
ARP2353	Elbow Block M (145mm wide x 95mm high)
ARP2354	Elbow Block L (165mm wide x 110mm high)
ARP2092	Onyx Auto Style headrest
ARP2324	Elan Standard Pad 255mm / 10" wide w/standard cover (ESP26) complete with Elan Headrest Hardware (MEHW)
ARP2325	Elan Occipital Pad 230mm x 180mm / 9" x 7" w/standard cover (EOP23) complete with Elan Headrest Hardware (MEHW)
ARP0586	Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 310mm - 410mm
ARP0592	LNx Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380mm - 419mm (Selectable length depends on STF height) Custom-made, requires customised setting
ARP2365	Angle & width adjustable footplates Small (127 width x 190 depth / 5" x 7.5")
ARP0559	Angle & width adjustable footplate (240mm depth)
ARP0516	Angle & width adjustable footplates (290mm depth)
ARP0517	Angle adjustable foot platform (240mm depth, 305mm width)
ARP0537	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARP0539	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm
ARP0535	Without Hanger and footplates
ARP0547	Vari Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARP0549	ADE POWERED : Angle adjustable by joystick 0° - 80°, length 290mm - 460mm
ARP0545	Without Hanger and footplates
ARP0505	Fixed Footplate
ARP0509	Angle Adjustable Footplate
ARP0511	Angle and depth adjustable Footplate
ARP0502	Fixed Footplate
ARP0506	Angle Adjustable Footplate
ARP0508	Angle and depth adjustable Footplate
ARP2203	REM211 Drive & Actuators
ARP2204	REM216 Drive, Actuators & Lights
ARP1207	Standard Bracket RH

ARP1208	Standard Bracket LH
ARP1707	Swing Away + height adjustable Bracket RH High end
ARP1708	Swing Away + height adjustable Bracket LH High end
ARP1709	Matrx Resolve Swing in/out bracket RH
ARP1710	Matrx Resolve Swing in/out bracket LH
ARP2075	Swing Away + height / width adjust mount on RH for user, Quad link design
ARP2076	Swing Away + height / width adjust mount on LH for user, Quad link design
ARP2227	Centre Mount swing-away Bracket to RH (Nucleus)
ARP2228	Centre Mount swing-away Bracket to LH (Nucleus)
ARP2209	REM500 Swing away RH multi-adjustable
ARP2210	REM500 Swing away LH multi-adjustable
ARP2213	REM500 Centre mount swing away to RH Nucleus
ARP2214	REM500 Centre mount swing away to LH Nucleus
ARP2075	Swing Away + height / width adjustable Bracket Mount on RH, Quad link design
ARP2076	Swing Away + height / width adjustable Bracket Mount on RH Quad link design
ARP2219	Compact Dual Control (DLX-ACU200) (Please select Dual control bracket below)
ARP1207	Dual Control Standard Bracket RH
ARP1208	Dual Control Standard Bracket LH
ARP2245	IDC Dual control (Please select AQE1221 OR AQE1222 for bracket option below)
ARP1221	IDC Dual control Detachable Bracket RH
ARP1222	IDC Dual control Detachable Bracket LH
ARP2231	Proportional Compact Remote DLX-CR400 (Please select Proportional Bracket)
ARP2232	Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)
ARP2207	Proportional Compact RH Mount Bracket
ARP2208	Proportional Compact LH Mount Bracket
ARP2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ARP2214	Proportional Compact Centre Mount swing away to LH (Nucleus)
ARX2215	Individual mounting option per request (Special) Requires additional information / specification
ARP2230	Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3
ARP2247	Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings
ARP2233	MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry)
ARP2213	Centre Mount, swing away to RH, Nucleus
ARP2214	Centre Mount, swing away to LH, Nucleus
ARP2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARP2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARP2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LiNX Heavy Duty Bracketry)
ARP2213	Centre Mount, swing away to RH, Nucleus
ARP2214	Centre Mount, swing away to LH, Nucleus
ARP2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARP2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARX2215	Individual mounting option per request (Special) Requires additional information / specification
ARP2227	10-Way Switch Bracket RH, Mounted for Attendant
ARP2228	10-Way Switch Bracket LH, Mounted for Attendant
ARP2257	Four-Switch Proximity Array (ASL106) NEW!

ARP2275	Centre Mount, swing-away to RH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm)
ARP2276	Centre Mount, swing-away to LH, Nucleus (Incl Eclipse Tray ASL619, 360mm x 180mm)
ARP2243	Proportional Paediatric Compact (ASK132), Includes input Module NEW!
ARP2213	Centre Mount, swing-away to RH, Nucleus
ARP2214	Centre Mount, swing-away to LH, Nucleus
ARP2241	Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm)
ARP2242	Centre Mount, swing-away to RH, Nucleus (incl Tray (ASL603, 100mm x 150mm)
ARP2255	Sip'n Puff Head Control (ASL109) NEW! includes Input Module and Lip Switch (ASL 301 LIP) (Combines simple sip and puff controls with head movements)
ARP2256	Sip'n Puff Control NEW! (including Sip'n Puff kit, Lip Switch ASL301LIP and Input module)
ARP2271	Input Module DLX-IN500-A (sub D, Sip'n Puff) NEW!
ARP2207	USB Charger Port RH Mount
ARP2208	USB Charger Port LH Mount
ARP2261	Stereo Toggle - direct control for Tilt only
ARP2262	Stereo Button - direct control for Tilt only
ARP2263	Four-Way Toggle - direct control for Tilt & Recline
ARP2264	Four-Way Button - direct control for Tilt & Recline
ARP2207	Toggle / Button RH Mount
ARP2208	Toggle / Button LH Mount
ARP1282	6 km/h (4mph)
ARP1437	LED Light System with Indicators (30 LUX)
ARP1404	Tray Swing Away RH
ARP1405	Tray Swing Away LH



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

**Special Instructions**

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_

Description of special request (please include relevant option code or part number):  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_