





Technical Data

KÜSCHALL CHAMPION

 Frame angle Ti: 75° / 90° Ca: 85°	 Seat width 360 – 460 mm in increments of 20 mm
---	---



 Seat depth Ti: 340 – 460 mm Ca: 400 – 460 mm in increments of 20 mm	 Seat height front: 450 – 540 mm rear: 390 – 490 mm in increments of 10 mm
--	---

 Knee-to-heel length Ti: 220 - 500 mm Ca: 300 - 340 mm 400 - 500 mm in increments of 10 mm	 Backrest height 300 – 465 mm in increments of 15 mm
---	--



 Total weight, unloaded Ti: approx. 10.4 kg Ca: approx. 9.9 kg	 Transport weight (w/out rear wheels) approx. 7.3 kg
---	--



 Max. user weight 120 kg
--



KÜSCHALL K-SERIES

 Frame angle Ti: 75° / 90° Ca: 75°	 Seat width 340 – 480 mm in increments of 20 mm
--	---

 Seat depth 375 – 450 mm in increments of 20 mm	 Seat height front: 480 – 520 mm rear: 380 – 490 mm in increments of 10 mm
---	---

 Knee-to-heel length 290 - 480 mm in increments of 10 mm	 Backrest height 270 – 480 mm in increments of 15 mm
--	--

 Total weight, unloaded Ti: approx. 8.1 kg Ca: approx. 7.9 kg	 Transport weight (w/out rear wheels) Ti: approx. 5.5 kg Ca: approx. 5.3 kg
---	--

 Max. user weight 120 kg	 Crash Test approved Fulfills the requirements of ISO 7176-19
---	---

Küschall AG | Benkenstrasse 260 | CH-4108 Witterswil | Switzerland | Phone +41 (0) 61 487 70 70 | Fax +41 (0) 61 487 70 79 | www.kuschall.com
UNITED KINGDOM: Invacare Ltd. | Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ | Phone +44 (0) 1656 776222 | Fax +44 (0) 1656 776220 | www.invacare.co.uk
IRELAND: Invacare Ltd. | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin | Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085 | www.invacare.ie
Invacare EASTERN EUROPE, MIDDLE EAST, RUSSIA & CIS | Kleiststraße 49 | D-32457 Porta Westfalica | Tel: +49 (0)57 31 754 540 | Fax: +49 (0)57 31 754 541 | edo@invacare.com

Küschall is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication or of parts thereof as well as any reproduction of images must have the written consent of Küschall AG.

2014 Copyright © 08/2014, Küschall AG, Switzerland - All rights reserved



The Prestige Collection

Contemporary designs and precision
in Titanium & Carbon

küschall[®]
designed for life

Exclusive materials for highly active people

The optimal selection for a perfect design. Experience superior driving performance with exceptional strength and aesthetic appeal. The high-end materials of Titanium and Carbon leave nothing else to be desired.

With sleek, contemporary designs combined with durable, lightweight materials, The Prestige Collection is the right choice for active individuals who have a strong affinity for precision and design, and will enhance daily life or any adventurous activity.

KÜSCHALL CHAMPION | KÜSCHALL K-SERIES

Ti
Titanium



Ca
Carbon

Whether your lifestyle is constantly on the go or moving in stylish circles, the **küschall K-Series** and the **küschall Champion** in Titanium or Carbon offer individuals something truly unique as a means of everyday mobility.

These contemporary designs are available in a foldable or a rigid frame. Versatile and highly configurable to suit a wide range of personal requirements, the **küschall K-Series** is designed for individuals who demand effortless mobility combined with efficient driving capabilities. With its ease of transportation and unique scissor folding mechanism, the **küschall Champion** is perfect for life on the go. Its lightweight, minimalist design and sleek aesthetics will give individuals the freedom to participate wherever they may be.

Navigating life with absolute confidence is the hallmark of **küschall wheelchairs - designed for you.**

