

TurnSafe 2

Sit-to-Stand & Transfer Aid




Safe Moving & Handling

Sit-to-stand and transfer support for users and caregivers



Safe, active and comfortable transfers

TurnSafe 2 is an ergonomic sit-to-stand and transfer aid with integrated turning function. Developed in close cooperation with medical professionals and users, it supports both the user's independence and the caregiver's safety.



Activates & Supports

- > **Large foot plate** with anti-slip surface and adjustable height
- > **Adjustable lower leg supports** – sideways and in height
- > **Sturdy, grip-friendly handles** enable the user to actively participate in transfers



Easy to Use & Manoeuvre








- > Relatively **small and convenient** – easy to have close at hand
- > **Simple to rotate** and manoeuvre in most spaces and situations
- > **Foot brake** on the caregiver's side provides stability and security



Protects the Caregiver

- > **Helps reduce the risk of work-related injuries** for carers
- > **Disassembles into three pieces** for easy transport and storage
- > **Suitable** for a wide range of transfer scenarios: bed, wheelchair, toilet

Technical data ▼

						
Total weight ▼	Foot plate diameter ▼	SWL ▼	Height Of Lower Handle User ▼	Height of upper handle, for the caregiver ▼	Height of lower leg support ▼	Lower leg supports adjustable in steps of ▼
15 kg	451 mm	150 kg	800-1146 mm	912-1259 mm	185-531 mm	50 mm

DHCare Ltd.
4D/E, The Ark, 201 Talgarth Road,
Hammersmith, London, W6 8BJ, UK
www.dhcare.com

©2026 DHCare Ltd. All rights reserved. DHCare Ltd reserves the right to modify existing models at any time. Specifications may vary across regions. TurnSafe 2_EU_29052026

