

TurnTable

Safe Moving & Handling



Makes it easier to turn when sitting down

For seated and standing transfers – simple, effective, close at hand



TurnTableSoft

Makes it easier to turn when sitting down

TurnTableSoft is mainly used on chairs and beds to facilitate turning when sitting down.

TurnTableSoft is constructed of two round discs with a friction-reducing material on the inside, enabling the discs to rotate smoothly over each other. The material resists soiling.



TurnTableHard

For transferring a user between wheelchair and chair, toilet or bed

Many users have enough strength to stand on their legs but not the ability to move their feet. In these circumstances, TurnTableHard, which is placed on the floor under the user's feet, is an outstanding, easy-to-use aid for the more experienced caregiver and for transferring a user between chair and wheelchair, wheelchair and toilet or from wheelchair to bed.

The caregiver can guide and slow the transfer by keeping half a foot on the TurnTableHard.



Technical data ▼



SWL ▼



Diameter ▼

Unit ▼

TurnTableSoft		45 mm / 40 mm	1
TurnTable Ø38cm	135 kg	38 cm	1

DHCare Ltd.
4D/E, The Ark, 201 Talgarth Road,
Hammersmith, London, W6 8BJ, UK
www.dhcare.com

©2026 DHCare Ltd. All rights reserved. DHCare Ltd reserves the right to modify existing models at any time. Specifications may vary across regions. TurnTable_EU_29052026

