



Product Service

CERTIFICATE

No. Z1 09 11 52355 016

Holder of Certificate: **Küschall AG**
Benkenstrasse 260
4108 Witterswil
SWITZERLAND

Production Facility(ies): 52355

Certification Mark:



Product: **Active wheelchairs**

Model(s): **Ultra-Light / Compact**

| Parameters: | Ultra-Light | Compact |
|--------------------------------|-------------|-------------|
| Total length incl. footrests | 730-1100 mm | 740-1160 mm |
| Total length without footrests | ---- | 570-850 mm |
| Wheelchair width | 490-750 mm | 490-750 mm |

- Further parameters see enclosure

Tested according to: DIN EN 12182:1999
DIN EN 12183:1999

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 71358709

Date, 2009-11-25 (Torsten Zimmer)
Page 1 of 2





Product Service

Enclosure to Certificate No.: Z1 09 11 52355 016

Model designation for the active wheelchairs Ultra-Light / Compact

| | Ultra-Light | Compact |
|---------------------------------|--|---|
| Total length incl. foot rests | 730 - 1100 mm | 740 - 1160 mm |
| Total length without foot rests | ---- | 570 - 850 mm |
| Wheelchair width | 490 - 750 mm | |
| Wheelchair height | 690 - 1010 mm | 670 - 1010 mm |
| Seat width | 280, 300, 320, 340, 360, 380, 400, 420 mm (Single diagonal strut) 440, 460, 480, 500 mm (Double diagonal strut) | |
| Seat depth | 320 - 500 mm | |
| Seat height at the front edge | 430 - 530 mm | 370 - 530 mm |
| Seat angle | -2° - 16° | |
| Width of the backrest | 280 - 500 mm | |
| Height of the backrest | 300 - 510 mm | |
| Angle of the backrest | 82° - 102° | |
| Max. user weight | 100 kg (Single diagonal strut) 130 kg (Double diagonal strut) | |
| Driving wheels | | |
| - Size | 22, 24, 25 and 26 x 1" - 1 3/4" | |
| - Tyres | <input checked="" type="checkbox"/> air-filled | <input checked="" type="checkbox"/> accident proof |
| - Adjustability | <input checked="" type="checkbox"/> yes | <input type="checkbox"/> no |
| - Kind of fixation | <input checked="" type="checkbox"/> quick release | <input type="checkbox"/> rigid |
| - Adjustment of the camber | 0° / 3° | |
| Front wheels | | |
| - Size | 3, 4, 5, 6 and 7 x 3/4" - 1 1/4" | |
| - Tyres | <input checked="" type="checkbox"/> air-filled | <input checked="" type="checkbox"/> accident proof |
| - Adjustability | <input checked="" type="checkbox"/> yes | <input type="checkbox"/> no |
| - Kind of fixation | <input type="checkbox"/> quick release | <input checked="" type="checkbox"/> rigid |
| Brake system | Knee lever brakes | |
| Frame material | Aluminium | |
| Arm rests | Height adjustable and removable | |
| Foot rests | Height adjustable, foot plate(s) can be swivelled upwards | Height adjustable, can be turn aside and removable, foot plate(s) can be swivelled upwards |
| Anti tip device | existing | |

Department: MHS 5 / me
Date: 2009-11-25

Page 2 of 2