

UNIQUE

Designed by you. Built by us.



Küschall®

2024 EDITION

CONTENTS

Each personality is **UNIQUE**.

Our team of experts have been pushing the boundaries of customisation to bring you a Unique range of personalisation options, allowing you to explore endless possibilities to create the perfect chair to suit any personality.

The only limit is **YOUR** imagination !

PERFORMANCE FRAMES 3 & 4

CERAMIC BEARINGS 5

SPINERGY REAR WHEEL OPTIONS..... 5 & 6

AESTHETICS

UPHOLSTERY 7

FRAME COLOURS & FRAME 8

CASTORS & FORKS..... 9

ANODISATIONS10 & 11


EXAMPLE UNIQUE WHEELCHAIRS 12 - 15

SAMPLE EMBROIDERIES 16



Dare to express your personality through your wheelchair and adopt the style that suits you.

PERFORMANCE FRAMES



KSL ERGONOMIC SEAT

Request via Unique specials, price on application.

Please provide the full configuration of the wheelchair, along with the precise measurements of the ergo required.



K-SERIES ANGLE ADJUSTABLE FOLDING BACKREST ON THE KSL

Option SWO4444

With carbon fibre flip-up and fold-in clothesguards on the backrest

Option SWO9999 With carbon fibre removable mudguards

74°	78°	82°	86°	90°



SK FOLDING FRONT FRAME FOR K-SERIES & THE KSL

Additional weight : +600g

Subject to suitable specification. Request via Unique specials providing full configuration details to be approved.

FRAMES



KSL WELDED UNDER SEAT GRIP BAR **Option code SWO3333**

Adds additional 180g to overall weight of wheelchair



CHAMPION SL WELDED COG **& WHEEL SUPPORTERS Option code 1260**

Increased rigidity Improved driving performance
Lighter by 800gm than standard Champion

Total weight: 8.4kg Transport weight: 5.9kg
Based on 40cm x 40cm width and depth



THE KSL WITH ADJUSTABLE CENTRE OF GRAVITY **Option code 6050**

Adjusts +/- 1 position in 15mm increments
Select from 1 – 5 positions Adds additional 200g in weight

SEAT POSITION (Measured from centre of backpost to centre of axle) *

<p>DSA6100 Position 1 = 140 mm MOST ACTIVE POSITION</p> <p><input type="radio"/></p> <p>Select</p>	<p>DSA6110 Position 2 = 125 mm</p> <p><input type="radio"/></p> <p>Select</p>	<p>DSA6120 Position 3 = 110 mm</p> <p><input type="radio"/></p> <p>Select</p>	<p>DSA6130 Position 4 = 95 mm</p> <p><input checked="" type="radio"/></p> <p>Select</p>	<p>DSA6140 Position 5 = 80 mm LEAST ACTIVE POSITION</p> <p><input type="radio"/></p> <p>Select</p>
--	---	---	---	--

CERAMIC BEARINGS



CERAMIC BEARINGS FOR FRONT CASTORS Option code 5510

Ceramic bearings on castors. Only for Aluminium softroll front castors



CERAMIC BEARINGS SPINERGY WHEELS Option code 6579

Ceramic bearings on rear wheels. Only for the SpinerGY Wheels



LOOP WHEELS Maximum user weight 120kg

A range of Loop wheels are available. Request via Unique specials, price on application.

COLOURS FOR SPINERGY WHEELS

Choose from a range of colours for SpinerGY Wheels



SPINERGY REAR WHEELS



SPINERGY Lite Extreme LX

Option code 6551

All-purpose, daily use, lightweight, travel. 12 spoke, the most popular wheel in the range. 687g per wheel.



SPINERGY Lite Extreme X-Laced XLX

Option code 6552

Reinforced, all-purpose, daily use, lightweight and active. 12 spoke, with the outside spokes crossed at the hub which contributes added stiffness with a distinctive aesthetic. 821g per wheel.



SPINERGY Lite Extreme X-Laced Sport XSLX

Option code 6553

24 x 4mm thick spoked wheels. Outside spokes are cross-laced, inside spokes are radially-laced. 757g per wheel.



SPINERGY Carbon Fibre CLX

Option code 6554

Ultra lightweight, professional, high-end appearance, carbon fibre. At 530g per wheel, the Spinerger CLX is the lightest, high end wheel offered on our active range.

Choose from a range of coloured anodised options listed above. Please refer to the order form for full details.

AESTHETIC UPHOLSTERIES



UPHOLSTERY STITCHING ON BACKREST

Option codes 2241, 2242, 2243, 2244, 2245, 0002

Neon Yellow, Neon Pink, Blue, White, Red or select colour of your choice.

Only available on black Aero upholstery.

Please refer to order form for choices, or select a colour of your choice.



CUSTOM FLOCKING DESIGN Option code 0004

Design size 280mm x 400mm

Available on black Aero upholstery

Provide PDF, AI or PNG file with the design of your choice. (Subject to copyright)


CUSTOM EMBROIDERY DESIGN

Option code 0003

Design size 250mm x 250mm embroidery.

Available on black Aero upholstery

Provide PDF, AI or PNG file with the design of your choice. (Subject to copyright)



FRAME COLOURS & DESIGNS



SPECIAL FRAME PAINT

Option code 0008


Request the RAL reference for the frame colour of your choice.



CERAMIC GLAZE FINISH

Option code 4066

High quality ceramic frame coating, increases gloss and shine.



HYDRO DIPPING FRAME

Have the design of your choice transferred onto the frame for a truly bespoke wheelchair.

Request via Unique specials, price on application.
Classic varnish included. Supply file with pattern required.
(Subject to copyright)



MIRROR POLISH FRAME

Option code 4101

Classic varnish coating included. Includes frame and backrest.

CASTORS



ALUMINIUM SOFTROLL CASTORS 4" Option code 5200 to access 4" castor size

Select the colour from options 5502, 5504, 5506, 5507, 5508 or 5509.

Choose from six anodised colours; pink, orange, blue, purple, red, gold.



ALUMINIUM SOFTROLL CASTORS 5" Option code 5450 to access 5" castor size

Select the colour from options 5502, 5504, 5506, 5507, 5508 or 5509.


Choose from six anodised colours; pink, orange, blue, purple, red, gold.

FORKS



SINGLE SIDED FROGLEGS FORKS

Request via Unique specials, price on application



FROG LEGS SUSPENSION FORKS

Request via Unique specials, price on application

ANODISATIONS



ANODISED HANDRIMS

Request via Unique specials, price on application.

Choose from the following colours:

Red, purple, copper, blue, green, orange, gold, champagne, pink or mirror finish.



ANODISED CARBOLIFE HANDRIMS

Request via Unique specials, price on application Available for:

Gekko / Gekko S

Curve / Curve L

Choose from the following colours: red, purple, copper, blue, green, orange, gold, champagne, pink or mirror finish

ANODISATIONS

ANODISATION KIT 1

Option code 6960



Forks and forks head caps. Applies to all models.


Choose from the following colours: red, purple, copper, blue, green, orange, gold, champagne, pink or mirror finish.

ANODISATION KIT 2

Option code 6961

Please refer to order form for full details. Components will vary depending on model selected.

Choose from the following colours: red, purple, copper, blue, green, orange, gold, champagne, pink or mirror finish.



ANODIZATION KIT 3

Option code 6962

Please refer to order form for full details. Components will vary depending on model selected.

Choose from the following colours: red, purple, copper, blue, green, orange, gold, champagne, pink or mirror finish.



SAMPLE WHEELCHAIRS

Take a look at some samples of chairs we've already made!



THE KSL

UNIQUE FEATURES USED:

KSL welded frame

SK Champion style folding front frames

K-series type folding backrest

White backrest stitching

White seat stitching

Single-sided Froglegs castor forks

Froglegs 4" castors

Ceramic bearing castors

Ceramic bearing rear wheels

Spinergy XLX wheels (Spokes: Black, Hubs: Mirror finish, Rims: Mirror finish)

Special Premium colour : Sage Green 733

Gekko S handrims painted in frame colour, Sage Green.

Ceramic Varnish coating on frame



THE KSL UNIQUE FEATURES USED:

KSL welded frame SK Champion type folding front frames

Ergonomic seat (12 cm)

K-series style folding backrest

Blue backrest stitching

Blue seat stitching

Single-sided Froglegs castor forks

5" Aluminium softroll castors

Ceramic bearings in castors

Ceramic bearings in rear wheels

Anodised kit 3: blue

Spinergy XLX: Spokes: Black, Hubs: Blue, Rims: Blue

Special Premium Colour: Glittered Black Ceramic varnish



THE KSL

UNIQUE FEATURES USED:

KSL welded frame

SK Champion style folding front frames

Ergonomic seat (12 cm)

K-series type folding backrest

Embroidery on backrest

Fluorescent yellow backrest stitching

Special Premium Painting - White OPAL Ceramic varnish

Printed Fluorescent yellow clothes-guards

Printed fluorescent yellow footplate

Single-sided Froglegs castor forks

Aluminium softroll 5" castors

Ceramic bearing castors

Ceramic bearing rear wheels

Phosphorescent yellow on the gripping bar

Spinergy XLX wheels with personalised Spokes: White, Hubs: Purple, Rims: Purple



THE KSL

UNIQUE FEATURES USED:

KSL welded frame

SK Champion type folding front frames

Ergonomic seat (12 cm)

K-series type folding backrest

Special Premium Paint - Glittered white

Blue backrest stitching

Single-sided Froglegs castor forks

Gekko S handrims anodized: blue

Anodised kit 3: blue

Ceramic varnish

Froglegs 5" castors


Ceramic bearing castors

Ceramic bearing rear wheels

Spinergy XLX 3 wheels with personalised Spokes: white, Hubs: blue, Rims: blue

SAMPLE EMBROIDERIES

Some examples of embroidery styles available



Contact : specialsuk@invacare.com

UNIQUE
Designed by you. Built by us.

This catalogue gives examples of just some of the options available via Unique, but the possibilities are endless. Please contact your local sales manager for further details. Invacare reserves the right to make changes to the options and models presented at any time without notice. © Invacare Ltd 2024