

VICAIR

Vicair[®] Hygiene Cushion



Comfort and pressure relief



Vicair® technology

The Vicair Hygiene cushion is developed by Vicair, known for their uniquely designed pressure relief cushions that incorporate the innovative Vicair® air cell technology - SmartCell. The cushion which was designed in collaboration with Invacare is for use to support those with vulnerable tissue whilst using an Ocean Ergo* shower chair. The structure provides a dynamic surface and support. The cushion contains several hundred small tetraformed air-filled cells, SmartCells, which have very low friction and rounded corners.

The cells therefore glide easily against each other and mould to the shape of the body. The result is an optimal combination of pressure relief, support and stability.



Vicair Hygiene Cushion is made up of 4 sections with the possibility of customisation, if needed, just like the other cushions in the Vicair range.

Vicair Hygiene Cushion offers positioning for optimal redistribution of pressure and helps to counteract shear. It can be used preventively and as an effective complement to the treatment of existing pressure ulcers.

Vicair Hygiene Cushion consists of a durable low friction material and a mesh non-slip material on the underside with fastening straps in each corner.

*The cushion is designed to fit the **Invacare Ocean Ergo Range** but may be used on other brands of hygiene chairs after a risk assessment.

Features



Equipped with fastening straps in each corner that can be individually adjusted to achieve optimal placement.







The cushion is made up of 4 individual sections with mesh material on the underside. Machine washable up to 60 degrees without removing the air cells, air dry or in a drying cabinet.



The cushion is fixed with 4 fastening points.

Technical data

For more detailed information regarding this product, including the product's user manual, please visit our website www.invacare.co.uk

	Size	Total product weight	Max. user weight	Max. cleaning temperature
Vicair Hygiene	 48x45x8 VHY/4845	 0.6 Kg	 250kg	 60°



T: +44 1656 776222
uk@invacare.com
www.invacare.co.uk

Invacare UK & Ireland
Unit 4 Pencoed Technology Park
Bridgend CF35 5AQ

