## VICAIR® CENTRE RELIEF 02

## OFFLOADING PRESSURE IN THE HIGH RISK AREA















For completely offloading of the rectal area



Adjustable non-contact area



Stable positioning

The Vicair Centre Relief O2 consists of 5 individual compartments, designed to contour under loading. Each compartment is designed to influence a specific area of the pelvis, striking a unique balance between preserving skin and underlying tissue at the ischials, coccyx and trochanters, positioning the pelvis and supporting the femurs. The unique feature is the circular compartment towards the rear of the cushion. Designed especially for (partial) offloading of the coccyx/rectal area. For specific needs the filling of each of the 5 compartments can be easily adjusted, for individuals who: (1) require (partial) offloading of the rectal area, (2) have coccyx problems, (3) have localized deep pressure injuries.

## **AVAILABLE COVERS**

Vicair O2 Top Cover, Comfair Cover, Incotec Cover and Reversible Cover.

PRODUCT PERFORMANCE	
Skin protection / Anti decubitu	IS
Pressure redistribution	***
Shear force reduction	
Microclimate control	
Impact & vibration damping	***
Positioning	
Lateral stability	
Frontal stability	<b>AAAA</b>
Sliding reduction	
Allows manual adjustments	<b>AAAA</b>
Increasing endurance	
Maintenance & Safety	
Durability	***
Reliability	
Functionality when misused	
Easy to maintain	
Level of hygiene	***
Functionality & Application	
Comfort	***
Ultra lightweight	***
Allows easy transfers	
Available standard sizes	
Maximum user weight	
Haximam user weight	

TECHNICAL SPECIFICATIONS		
10 cm		
750 grams		
250 kg		
60°C max.		
	10 cm 750 grams 250 kg	

CUSHION SIZES	
Standard	
On request	

