


VicAIR celleoversigt - gram

Smartcells (mørkeblå)

Dato		Fornavn	
Pudetype		Efternavn	
Pudestr.		CPR-nr	


Vector 6 kamre, 10 cm

Str.	HMI-nr.	1+2	3	4-5	6
3540	43263	172,5	45	24	37,5
4040	43264	204	52,5	27	43,5
4045	43265	231	60	30	49,5
4050	43266	258	67,5	34,5	55,5
4343	43267	241,5	61,5	31,5	52,5
4540	43268	241,5	64,5	31,5	54
4545	43269	273	73,5	34,5	61,5
4550	43270	306	82,5	39	69
5045	43271	321	87	40,5	75
5050	43272	358,5	97,5	45	84



Vector Junior


Str.	HMI-nr.	1 + 2	3	4 + 6	5
2525	45888	120	37,5	7,5	30
3030	45889	172,5	70,5	15	57
3535	45890	285	112,5	33	90



Yes, you can.®


Twin, 6 cm

Str.	HMI-nr.	For	Bag
3540	24967	180	168
4040	19453	208,5	195
4045	19454	238,5	222
4050	24968	268,5	249
4343	24826	244,5	228
4540	24969	238,5	222
4545	19455	271,5	253,5
4550	24970	306	285
5050	24971	345	322,5




Twin, 10 cm

Str.	HMI-nr.	For	Bag
3540	24980	270	229,5
4040	19456	312	267
4045	19457	355,5	304,5
4050	24981	390	342
4343	24982	366	313,5
4540	24983	355,5	304,5
4545	19458	405	346,5
4550	24984	456	391,5
5050	24985	514,5	442,5




Academy Positioner Plus, 6 cm

Str.	HMI-nr.	For	Bag	Side
3540	24972	64,5	111	78
4040	24758	76,5	129	91,5
4045	24973	88,5	145,5	106,5
4050	24974	100,5	162	121,5
4343	24975	91,5	150	109,5
4540	24976	88,5	145,5	106,5
4545	24977	102	165	123
4550	24978	117	184,5	142,5
5045	34902	117	184,5	142,5
5050	24979	133,5	205,5	163,5



Academy Positioner Plus, 10 cm


Str.	HMI-nr.	For	Bag	Side
3540	19459	81	112,5	138
4040	24962	96	130,5	160,5
4045	24759	112,5	150	184,5
4050	24963	129	171	210
4343	24760	115,5	156	190,5
4540	19460	112,5	150	184,5
4545	19461	130,5	174	213
4550	24964	151,5	198	243
5045	34903	151,5	198	243
5050	24965	175,5	226,5	276



Yes, you can.®


Academy Adjuster, 6 cm

Str.	HMI-nr.	For	Midt	Bag
3540	30687	31,5	79,5	79,5
4040	30688	37,5	93	93
4045	30689	42	108	108
4050	30690	48	123	123
4343	30691	43,5	111	111
4540	30692	42	108	108
4545	30693	49,5	124,5	124,5
4550	30694	55,5	142,5	142,5
5045	34904	55,5	142,5	142,5
5050	30695	63	162	162




Academy Adjuster, 10 cm

Str.	HMI-nr.	For	Midt	Bag
3540	24959	46,5	126	126
4040	19462	54	147	147
4045	24957	63	168	168
4050	24958	70,5	190,5	190,5
4343	24827	64,5	174	174
4540	19463	63	168	168
4545	19464	72	193,5	193,5
4550	24960	81	219	219
5045	34905	81	219	219
5050	24961	93	249	249




Academy Multi functional Back

Str.	HMI-nr.	1+6	2+5	3+4	7	8
3645	35860	55,5	43,5	27	81	34,5
4054	35861	82,5	60	45	127,5	55,5
4554	35862	90	60	60	163,5	73,5
5054	35863	121,5	91,5	49,5	127,5	52,5



Academy Multi functional Seat


Str.	HMI-nr.	For	Midt	Bag
3752	35856	67,5	120	150
4052	35857	72	132	159
4552	35858	81	153	175,5
5052	35859	96	181,5	220,5



Yes, you can.®

Academy Back System

Str.	HMI-nr.	1+7	2+6	3+5	4	8	9
1	37298	31,5	43,5	39	15	15	15
2	37299	34,5	46,5	42	21	18	18
3	37300	37,5	49,5	45	27	21	21
4	37301	40,5	52,5	48	30	24	24
5	37302	43,5	55,5	51	33	27	27
6	37303	46,5	58,5	54	36	30	30
7	37304	49,5	61,5	57	39	33	33
8	37305	52,5	64,5	60	42	36	36



Velcroåbning på bagsiden til fyldning i nr. 9



Yes, you can.®