


Vicair celleoversigt - stk. antal

Smartcells (sølvgrå)

Dato		Fornavn	
Pudetype		Efternavn	
Pudestr.		CPR-nr	


Vicair Vector O2 9 kamre, 6 cm

Str.	HMI-nr.	1+2	3+4	5/6	7/8	9
3540	105530	37	26	12	30	29
3545	113637	42	29	13	34	33
4040	105537	43	30	14	35	34
4045	105536	48	34	15	39	38
4050	105535	53	37	17	43	42
4343	105534	49	35	16	40	39
4345	113102	49	35	16	40	39
4540	105533	48	34	15	39	38
4545	105532	54	38	17	44	42
4550	105528	60	42	19	49	47
5045	105523	60	42	19	49	47
5050	105531	67	47	21	54	52



Vicair Vector O2 9 kamre, 10 cm


Str.	HMI-nr.	1+2	3+4	5/6	7/8	9
3540	105519	49	36	21	43	35
3545	113634	55	41	23	48	39
4040	105520	56	41	23	49	40
4045	105521	63	46	26	55	45
4050	105529	70	51	29	61	50
4343	105522	65	48	27	56	46
4345	113093	68	50	28	59	48
4540	105518	63	46	26	55	45
4545	105524	71	52	30	62	50
4550	105525	79	58	33	69	56
4555	113096	87	64	36	76	61
5040	128029	70	51	29	61	50
5045	105525	79	58	33	69	56
5050	105527	88	64	37	76	62
5055	113099	97	81	40	84	68
5543	128032	83	61	35	72	54
5545	128035	87	64	36	72	61
5550	128028	97	71	40	84	68



Yes, you can.®


Vicair Vector O2 8 kamre, 6 cm

Str.	HMI-nr.	1+2	3+4	5	6+7	8
3540	63443	39	54	26	10	23
4040	63444	62	45	30	11	26
4045	63445	69	50	33	13	29
4050	63452	77	56	37	14	32
4343	63446	71	52	34	13	30
4540	63447	69	50	33	13	29
4545	63448	78	57	37	14	33
4550	63449	87	63	41	16	36
5045	63450	87	63	41	16	36
5050	63451	96	70	46	18	40




Vicair Vector O2 8 kamre, 10 cm

Str.	HMI-nr.	1+2	3+4	5	6+7	8
3540	61784	68	54	30	16	27
4040	61785	78	61	35	18	30
4045	61786	88	69	39	20	34
4050	61787	98	77	43	23	38
4343	61788	90	71	40	21	35
4540	61789	88	69	39	20	34
4545	61790	99	78	44	23	39
4550	61791	110	86	49	26	43
5045	61792	110	86	49	26	43
5050	61793	122	96	54	28	48



Vicair Active O2, 10 cm


Str.	HMI-nr.	1 + 4	2 + 3
3540	108714	52	37
4040	108715	59	43
4045	108716	66	48
4050	108717	74	54
4540	108718	66	48
4545	108719	75	54
4550	108720	83	60
5045	108721	83	60
5050	108722	92	67



Yes, you can.®



Vicair Adjuster O2, 6 cm

Str.	HMI-nr.	For	Midt	Bag
3540	106741	22	60	60
4040	106733	25	68	68
4045	106743	28	77	77
4050	106744	31	85	85
4343	106745	29	79	79
4540	106746	28	77	77
4545	106747	32	86	86
4550	106748	35	96	96
5045	106750	35	96	96
5050	106749	39	107	106



Vicair Adjuster O2, 10 cm

Str.	HMI-nr.	For	Midt	Bag
3540	106739	33	88	88
4040	106732	37	100	100
4045	106737	42	113	113
4050	106738	46	125	126
4343	106736	43	116	116
4540	106734	42	113	113
4545	106735	47	127	127
4550	106740	52	141	141
5045	106742	52	141	141
5050	106751	58	157	157




Vicair Multifunctional O2 rygpute

Str.	HMI-nr.	1+9	2+3+7+8	5	4 +6	10
3754	124838	35	40	35	33	30
4054	124839	38	43	38	35	30
4554	124840	43	49	42	40	30
5054	124841	48	54	47	44	30

Vicair Multifunctional O2, 8 cm

Str.	HMI-nr.	For	Midt	Bag
3752	35856	47	88	99
4052	35857	51	95	108
4552	35858	57	107	121
5052	35859	63	119	134



Yes, you can.®