



Product Service

CERTIFICATE

No. Z1 10 08 42097 044

Holder of Certificate: Invacare Rea AB
Org. No. 556098-9211
 Växjövägen 303
 343 71 Diö
 SWEDEN

Production Facility(ies): 42097

Certification Mark:



Product: Standard wheelchairs
Model(s): Manual wheelchair type: XLT Active, XLT Dynamic, XLT Swing, XLT Max

Parameters:

	total length incl. footrests:			
	XLT Active	XLT Dynamic	XLT Swing	XLT Max
	830–1030 mm	730–880 mm	820–1230 mm	930-1080 mm

	max. width:			
	XLT Active	XLT Dynamic	XLT Swing	XLT Max
	555–680 mm	555–680 mm	530–705 mm	815 mm

	max. height:			
	XLT Active	XLT Dynamic	XLT Swing	XLT Max
	610–1050 mm	610–1050 mm	610–1050 mm	1100 mm

Further Technical Data: see enclosure

Tested according to: EN 12183:2009

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 71373079

Date, 2010-08-03 (Torsten Zimmer)

Page 1 of 2





Enclosure to certificate no.: Z1 10 08 42097 044

Technical data for the manual wheelchair
Type: XLT Active, XLT Dynamic, XLT Swing, XLT Max

	XLT Active	XLT Dynamic	XLT Swing	XLT Max
total length incl. footrests:	830 – 1030 mm	730 – 880 mm	820 – 1230 mm	930-1080 mm
total length without footrests:	-	-	675 – 825 mm	-
max. width:	555 – 680 mm		530 – 705 mm	815 mm
max. height:	610 – 1050 mm			1100 mm
seat width:	355, 380, 405, 430, 455, 480 mm		330, 355, 380, 405, 430, 455, 480, 505 mm	505, 555, 605 mm
seat depth:	360 – 400, 410 – 450, 460 - 500 mm			460 – 500 mm, 510 – 550 mm, 560 – 600 mm
seat height:	400 – 520 mm	400 – 500 mm	390 – 520 mm	440 – 520 mm
backrest height:	200 – 490 mm			400-490 mm
seat plane angle:	0° - 14°			
backrest angle:	-10° - 10°			
max. load:	135 kg			180 kg
brake system:	knee lever brake (user brake)			
back wheels				
- size:	22 and 24 x 1 ³ / ₈ " , 24, 25 and 26 x 1"			24 x 1 ³ / ₈ " and 24 x 1"
- tyres:	air filled / accident proof			
- adjustable:	yes			
- fixation:	Quick release			
front wheels				
- size:	ø 75 x 24, ø 100 and 125 x 27, ø 140 x 37, ø 150 x 27, ø 150 x 30, ø 150 x 32, ø 180 x 40 mm			ø 140 x 37, ø 150 x 27, ø 150 x 30, ø 150 x 32 and ø 180 x 40 mm
- tyres:	air filled / accident proof			
- adjustable:	yes			
- fixation:	rigid			
frame material:	Titanium			
armrest:	height adjustable and removable			
legrest:	height adjustable	height adjustable and removable		height adjustable
anti tip device:	adjustable			

Department: MHS5 / me
 Date: 2010-08-03