



Anti corrosion statement – Invacare Mirage.

All metal components, sub assemblies and fasteners used on the above product are coated with a powder polyester paint, clear zinc & passivate, chrome plate, galvanised or anodised finish to protect them against corrosion.

Invacare will not be liable for premature corrosion of the product, if it has been used for purposes which it was not intended.

Invacare strongly recommend that the product is maintained in line with guidelines laid down in the User Guide to ensure integrity and longevity of the product.

Neil Harvey
Quality Engineering Manager.
October 2003.