



Aquatec® | **Transfer aids**



Aquatec **Trans**

- Extra large seat surface
- Safe: turntable is only unlocked by the weight of the user
- Seat surface turns and slides easily
- Soft, skin-friendly cover: easy-care, wipes clean



Aquatec **Disk**

- Extra large seat surface
- Seat surface turns very easily
- Soft, skin-friendly cover: easy-care, wipes clean
- Within seconds, the **Disk**, the shower stool **Dot** and an adapter are combined to a shower stool with a swivel seat

AQUATEC®



Yes, you can.®



Aquatec **Step**

- Stable frame
- Feet with non-slip rubber caps
- Non-slip, removable cover, machine-washable

Technical data

Order numbers	Aquatec Trans	Aquatec Disk	Aquatec Step
Turning and transfer aid, blue	4.03.002		
Turning and transfer aid, grey	4.03.008		
Rotating disk, blue		4.02.002	
Rotating disk, grey		4.02.008	
Step, blue			3.00.002
Step, white			3.00.001

Dimensions	Aquatec Trans	Aquatec Disk	Aquatec Step
Diameter	36.5 cm	36.5 cm	–
Length	53 cm	–	40 cm
Width	–	–	30 cm
Height	2.5 cm	2.5 cm	16.5 cm

More info	Aquatec Trans	Aquatec Disk	Aquatec Step
Weight	2.2 kg	1.6 kg	1.5 kg
Load capacity	130 kg	130 kg	150 kg
Colours	Blue, grey	Blue, grey	Blue, white

Invacare International Sarl
 Route de Cité-Ouest 2 - 1196 Gland
 Switzerland
 Tel: +41 (0)22 354 60 10
 Fax: +41 (0)22 354 60 11
 hqeuropa@invacare.com
 www.invacare.eu.com

© 2010 Invacare® International Sarl
 All rights reserved.

All information quoted is believed to be correct at time of print.
 Invacare® reserves the right to alter product specifications without prior consultation.

Transfer aids - EU - 09/2010



Yes, you can.®