



A tailor-made solution for heavier users

Introducing a new passive (with tilt and recline) manual wheelchair for heavier users. The *Rea Azalea Max* is designed to meet the specific needs of disabled users with extra weight.

The *Azalea Max* is a niche product and a complementary addition to the existing *Azalea* family from Invacare®. Thanks to the introduction of *Azalea Max*, dealers seeking a truly complete range of passive wheelchairs need look no further.



Yes, you can.®



Features and Options

Continued Compatibility

The *Rea Azalea Max* offers the same features as the rest of the *Rea Azalea* family such as tilt-in-space with weight-shifting and DSS (Dual Stability System) mechanism. Stability in this wheelchair originates from body length compensation in the back and leg rests, beneficial for users with limited movement and that have difficulty changing position.

High ergonomic push handle

The ergonomically shaped push handles are positioned close together for easier handling. An alternate version with higher handles offers a different position for the hands and is available as an option.



Two-way backrest

The lower and upper parts of the shaped backrest are interchangeable, giving users the option of more space in the seat or chest area.



Footplate converter

Divided into three parts, the footplate comes with a removable mid section, to facilitate transfers.



Lower-leg support

The fixed calf plate adjusts and swivels to provide optimum support for each leg.

Technical data

Azalea Max	55; 61 + 5cm	53 - 60	40 / 45	50 + 8 cm	32 - 42	33 - 50	1° - 20°
Azalea Max Assist	55; 61 + 10 cm	53 - 60	40 / 45	50 + 8 cm	32 - 42	33 - 50	1° - 20°
Azalea Max	0°-30°	SW + 25 cm	102 - 124	105 - 112	54 kg		
Azalea Max Assist	0°-30°	SW + 22 cm	102 - 124	105 - 112	54 kg		
Azalea Max	33 kg	max. 180 kg	For other <i>Rea Azalea</i> accessories, please consult the <i>Rea Azalea</i> brochure.				
Azalea Max Assist	37 kg	max. 180 kg	The <i>Rea Azalea Max</i> complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: www.invacare-rea.com under Rea wheelchairs/Technical downloads.				

Invacare International Sarl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeuropa@invacare.com
www.invacare.eu.com

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Rea Azalea Max - EU - 08/2010



Yes, you can.®

Frame colour



Pearl Grey



Grey Dartex TR23



Grey plush TR 32

Upholstery



Rea, Rea design and DSS (Dual Stability System) design are registered trademarks of Invacare® International.