

Certificate

No.: Z1 00 12 42097 002



Invacare Rea AB
Org. No. 556098-9211
Växjövägen 303
S - 343 21 Diö

with production facility(ies)
42097

is authorized to label the following products with the
TÜV Mark
concerning the enclosed attachment. See also notes overleaf.

Product: Rollstühle für Behinderte
Wheelchair for handicapped

Model: Handrim wheelchair
Type: Bellis, Bellis Simplex, Bellis Transport

Parameters:

seat width:	39 - 41 - 43 cm, 44 - 46 - 48 cm, 49 - 51 - 53 cm
seat height:	adjustable between 36 - 48,5 cm
seat depth:	40 - 52 cm
seat angle:	(-1°) - 15°
back rest height:	52 cm (adjustable to 7 cm), 0° - 30°, gas spring, mechanical
drive wheels:	20", 22", 24"
transport wheels:	16"
steering wheels:	140, 150, 180, 200 mm
total height:	110 cm
total length:	92 cm
total width:	seat width + 20,5 cm
max. load:	125 kg

The compliance of the above-mentioned product with the following essential requirements was tested on a voluntary basis:

EN 12183: 1999
TÜV PRODUCT SERVICE PPP
"mechanical wheelchairs" 03/98

Report No: 04610054473, E - (Z1 00 02 18026 031)

Released with the above certificate number by the
certification body of TÜV PRODUCT SERVICE GMBH.

Department: HANMED5 / me
Date: 12-08-2000

