



Yes, you can.®

➤ Invacare® **Docking Station**

Simple, Safe, Secure

The *Invacare*® **Docking Station** is designed to make wheelchair transportation so much easier: a simple switch of a button will safely fix and unfix the wheelchair inside a vehicle.

Users are provided peace of mind, as a mini display indicates that the connection is safe and the wheelchair locked in place.

The **Docking Station** meets the highest safety standards and conforms to ISO regulations.



NEW

Simple and convenient

The *Invacare*[®] **Docking Station** is simple to dock and easy to use. An automatic connection and guided track ensures trouble-free docking.



For added convenience, the **Docking Station** easily adapts to power wheelchairs and can be mounted in most vehicles.



ISO 'crash test' conformity

The **Docking Station** has proven stability and complies to ISO 10542-3 and ISO 7176-19 crash test requirements and related safety standards.

Features



Technical data



25 kg




700 mm - 1000 mm
(depending on PWC model)




120 mm




400 mm



Grip of T-anchor:
70 - 90 mm
above ground



260 kg
(PWC including load)



Van fixation:
floor plate, by
8 screws M10

The Docking Station varies slightly for the different PWCs.

Order numbers:	Storm	Typhoon	Typhoon II
Docking Station (car part)	I498274	I504830	I510997
Fixation Docking Station (chair part, retrofit kit)	SP1498275	SP1505049	SP1505049

All information quoted is believed to be correct at time of print. Invacare[®] reserves the right to alter product specifications without prior consultation.



INVACARE[®] Limited
South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY
Tel. +44 (0) 1656 647327 · Fax +44 (0) 1656 649016
E-mail: uk@invacare.com · www.invacare.co.uk
Sales Order E-mail: ordersuk@invacare.com

INVACARE[®] Ireland Ltd.
Unit 5 · Seatown Business Campus
Seatown Road · SWORDS · County Dublin
Tel. (353) 181 107 84 · Fax (353) 181 107 85
Email: ireland@invacare.com · www.invacare.ie