



Yes, you can.

Invacare® **EC Gantry** – Portable Hoist Systems



Invacare® **EC Gantry** Available in two versions

- Can be adjusted in 6 different heights with intervals at 7 cm.
- Easy to assemble and disassemble.
- **Gantry** with castors, height adjustable from 195 cm – 230 cm.
- **Gantry** without castors, height adjustable from 216 cm – 251 cm.

Invacare® **EC Gantry Traverse** No need for installing tracks in wall or ceiling.

- **EC - Gantry Traverse** is very easy to install and dismount
- Once mounted it is very easy to raise and lower the gantry without using tools
- **EC- Gantry** is height adjustable from 218 to 253 cm.
- Can be configured in all sizes from 2,5 m x 2,0 m to 4,0 m x 4,0 m



Technical data




200 kg



95°

Colour



white (Standard) aluminium (Optional)