



Product Service

CERTIFICATE

No. B 08 11 42940 025

Holder of Certificate: Invacare EC-Hong A/S

Ostergade 3
4270 Hong
DENMARK

Production Facility(ies):

68663

Certification Mark:



Product:

Components for hoists

Model(s):

Non rigid body support units

Parameters:

Model identification see enclosure

Remark:

The tests were performed according to ISO 10535:2006 only clause 8

Tested according to: ISO 10535:2006

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. The certification mark must not be altered in any way. See also notes overleaf.

Test report no.:

71301913-1 to -12

Date, 2008-11-27

Page 1 of 2





Enclosure to certificate no.: B 08 11 42940 025

Model identification for the non rigid body support units

Non rigid body support unit type:	Universal Low, Universal Standard, Universal High and Universal High Plus	Non rigid body support unit type:	Dress Toileting and Dress Toileting High
Size:	XS, S, M, L or XL	Size:	S, M, L or XL
SWL:	200 kg	SWL:	200 kg
Maximum washing Temperature:	95 °C	Maximum washing Temperature:	95 °C
Material:	100 % Polyester	Material:	100 % Polyester
Non rigid body support unit type:	Easy-Fit	Non rigid body support unit type:	Comfort Toileting
Size:	XS, S, M, L or XL	Size:	XS, S, M, L or XL
SWL:	200 kg	SWL:	200 kg
Maximum washing Temperature:	95 °C	Maximum washing Temperature:	95 °C
Material:	100 % Polyester	Material:	100 % Polyester
Non rigid body support unit type:	Transfer Stand Assist	Non rigid body support type:	Amputee
Size:	S, M, L or XL	Size:	S, M, L or XL
SWL:	159 kg	SWL:	200 kg
Maximum washing Temperature:	95 °C	Maximum washing Temperature:	95 °C
Material:	100 % Polyester	Material:	100 % Polyester
Non rigid body support unit type:	Standing Transfer Vest	Non rigid body support unit type:	Comfort Net
Size:	XS, S, M, L or XL	Size:	XS, S, M, L or XL
SWL:	200 kg	SWL:	200 kg
Maximum washing Temperature:	95 °C	Maximum washing Temperature:	95 °C
Material:	100 % Polyester	Material:	100 % Polyester
Non rigid body support unit type:	Groin Band for Standing Transfer Vest	Non rigid body support unit type:	Comfort Spacer
Size:	S or L	Size:	XS, S, M, L or XL
SWL:	159 kg	SWL:	200 kg
Maximum washing Temperature:	95 °C	Maximum washing Temperature:	95 °C
Material:	100 % Polyester	Material:	100 % Polyester
Non rigid body support unit type:	Stand Assist	Non rigid body support unit type:	Lifting Band
Size:	S, M, L or XL	Size:	Short or Long
SWL:	159 kg	SWL:	200 kg
Maximum washing Temperature:	95 °C	Maximum washing Temperature:	95 °C
Material:	100 % Polyester	Material:	100 % Polyester

Department: MHS5 / tz
Date: 2008-11-27