

Invacare® | *Apollo Mk IV*



Indoor practicality

With a small wheelbase, the Indoor MkIV is more manoeuvrable in and around the home. It offers great value for money with proven **Apollo** quality and is unrivalled in the market place. The **Apollo** has a fully padded seat and back as standard, and is available in 5 different seat widths and 3 seat depths to suit all client requirements.



Yes, you can.®



Features and Options

Indoor version

The **Invacare Apollo Mk IV** is the ideal indoor chair. Its narrow styling and small turning circle make it the practical choice in confined spaces, and the sensitive programmable controller gives the user the ability to manoeuvre confidently, in full control at all times.

The 30Ah batteries provide enough power to get around all day, and the handy plug-in charger in the controller allows recharging of the batteries when needed. Simply plug-in to any household mains supply to recharge.



Footrests

The adjustable footrests swing away, allowing easier access to the front of the chair for both versions.

Folding backrest

Reduces the overall size of the chair when folded for transportation or storage.

Technical data

Invacare Apollo	380 - 510 mm	410 - 460 mm	480 mm	495 mm	220 mm	330 - 430 mm
Invacare Apollo	10°	580 - 705 mm	975 mm	870 mm	620 mm	
Invacare Apollo	appr. 53 kg	max. 127 kg	30 Ah	1610 mm	max. 50 mm	
Invacare Apollo	max. 20%	160 W	appr. 12,5 km	max. 4 km/h		

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel: +353 1 8107084
Fax: +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2013 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Apollo Mk IV - UK - 08/2013 - Q99407

Frame colours



Black

More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.

