

Invacare® | **Commode chairs**



Special: Also available with a seat width of 60 cm!



H720T Special 60-48

Invacare H720T

- Stable, robust frame made from powder-coated steel tube
- High corrosion protection
- Space-saving push bar for optimal positioning over the toilet
- Swing-away armrests and swing-away, detachable footrests make conversion easy
- Stable commode seat with integrated toilet pan holder

AQUATEC®



Yes, you can.®



Invacare **Omega H750**

- Compact modern design
- Castors with brakes for safety
- Padded back, seat and arms
- Integral foot support



Invacare **H450**

- Stationary commode chair
- Space-saving: commode chair can be stacked
- Stable and corrosion-protected design with welded tubes
- Easy handling: commode pan can be removed from the side and from above



Invacare **H450LA**

- Stationary, height adjustable commode chair
- Space-saving: commode chair can be stacked
- Seat height adjustable in 7 steps from 51 - 66 cm
- Attractive design: fits well into a domestic setting
- Easy to clean thanks to the smooth surfaces



Invacare **Omega H407**





- Foldable, Space saving commode
- Ideal for infrequent or travel use
- Contoured backrest
- Complete with sanitary pan

Technical data

Order numbers	Invacare H720T	Invacare Omega H750	Invacare H450 stationary	Invacare H450LA height adjustable	Invacare Omega H407
Commode chair, mobile	1525770				
Commode chair, mobile, in special construction	H720T S				
Commode chair, mobile, lavender blue		1471043			
Commode chair, mobile, blue		1471044			
Commode chair, lavender blue			1471012	1471014	
5 pack			1471163	1471165	
Commode chair, blue			1471028	1471029	
5 pack			1471164	1471166	
Commode chair, lavender blue, 2 pack					1471010

Dimensions					
Seat width	45 cm (A)	50 cm (A)	49,5 cm (A)	49,5 cm (A)	44 cm (A)
Seat depth	45 cm (B)	46 cm (B)	45,5 cm (B)	45,5 cm (B)	42 cm (B)
Seat height	48 cm (C)	48 cm (C)	50 cm (C)	51 - 66 cm (C)	45 cm (C)
Total width	57,5 cm (D)	63,5 cm (D)	67 cm (D)	67 cm (D)	58 cm (D)
Total depth	81 cm (E)	79 cm (E)	63 cm (E)	63 cm (E)	56 cm (E)
Total height	94,5 cm (F)	92 cm (F)	89 cm (F)	90 - 105 cm (F)	91 cm (F)
Depth without footrests	54,5 cm (G)	67 cm (G)	-	-	-
Footplate depth	15,5 cm (H)	9,5 cm (H)	-	-	-
Height to armrests	75,5 cm (I)	70 cm (I)	69 cm (H)	70 - 85 cm (H)	66 cm (H)
Height, footrests to seat	34,5 - 46 cm (J)	37 cm (J)	-	-	-
Height, ground to lower seat edge:					
with holder for toilet pan	42 cm (K)	36 cm (K)	37 cm (I)	38 - 53 cm (I)	-
without holder for toilet pan	-	-	43 cm (J)	44 - 59 cm (J)	-
Height, floor to cover plate	52,5 cm (L)	50,5 cm (L)	52,5 cm (K)	53,5 - 68,5 cm (K)	47,5 cm (K)
Width between the armrests	47 cm (M)	51 cm (M)	51 cm (G)	51 cm (G)	46 cm (G)
Width, cut out in seat plate	22 cm	21,5 cm	21,5 cm	21,5 cm	21,5 cm
Depth, cut out in seat plate	26 cm	25,5 cm	25 cm	25 cm	25,5 cm
Ø Castors	5" (127 mm)	100 mm	-	-	-

More info					
Weight	14,0 kg	15,0 kg	8,3 kg	8,8 kg	8,5 kg
Load capacity	135 kg	130 kg	130 kg	130 kg	80 kg
Colour	Black anthracite	Lavender blue / Blue	Lavender blue / Blue	Lavender blue / Blue	Lavender blue

 Invacare **H720T**

 Invacare **Omega H750**

 Invacare **H450 / H450LA**

 Invacare **Omega H407**




Invacare Styxo 9630E

- Height adjustable from 39 cm to 54 cm
- Easy to clean
- Complete with seat, sanitary pan and lid


Technical data

Order numbers	<i>Invacare Styxo 9630E</i>
Toilet frame	9630E

Dimensions	
Seat width	38 cm
Seat depth	35 cm
Seat height	41-56 cm
Total width	62 cm
Total depth	49 cm
Total height	58 - 73 cm

More info	
Weight	5,5 kg
Load capacity	135 kg
Colour	white / grey

Invacare Styxo 9630E



Invacare International Sarl
 Route de Cité-Ouest 2 - 1196 Gland
 Switzerland
 Tel: +41 (0)22 354 60 10
 Fax: +41 (0)22 354 60 11
 hqeuropa@invacare.com
 www.invacare.eu.com

© 2010 Invacare® International Sarl
 All rights reserved.

All information quoted is believed to be correct at time of print.
 Invacare® reserves the right to alter product specifications without prior consultation.

Commode chairs - EU - 09/2010



Yes, you can.®

