



Easy Back

Easy Adjustability – Straight from the box!

The Easy Back provides the level of adjustability you would expect from a custom backsystem – directly from the box! This one piece back system can be fitted in seconds with just an allen key. Available in a wide range of sizes, it can be ordered with or without a headrest.

Designed to fit 7/8" & 1" tubing, this back is compatible with most of our *Invacare*® Power and Manual wheelchairs. By adding an *Invacare*® Drop Seat Board and *Flo-tech*® cushion, makes this product combination a one stop shop for all your mobility, seating and positioning needs.





Split Wings

Two independent support panels allowing individual adjustments. Ideal for clients requiring a narrower fit at the lateral area with a wider fit at the hips.




Multi Adjustable Brackets

A simple modification using just an allen key and spanner enables finely tuned adjustments to give a truly custom fit.



Reversed Dartex® Material

Two-way stretch vapour permeable, and water resistant cover. Easy to clean to aid infection control.



Easy Back Measuring Guide



Multi-size Brackets

All clamps and brackets fit both 7/8" & 1" tubing.



Quick Release Clamp Hardware

Allows quick and easy installation.

Back Size	A Dimension	B Dimension	C Dimension	D Dimension	E Dimension	F Dimension
12"/305mm	7"/178mm	8"/203mm	13"/330mm	3"/76mm	2"/51mm	7"/178mm
14"/355mm	9.5"/241mm	10.5"/266mm	15"/381mm	4"/101mm	3"/76mm	9"/228mm
16"/406mm	10.5"/266mm	11.5"/292mm	17"/432mm	4"/101mm	3"/76mm	11"/279mm
18"/457mm	11"/279mm	12.5"/317mm	19"/482mm	4"/101mm	3"/76mm	13"/330mm
20"/508mm	12"/305mm	13.5"/343mm	21"/533mm	4"/101mm	3.5"/89mm	15"/381mm
22"/559mm	13"/330mm	14.5"/368mm	23"/584mm	4"/101mm	3.5"/89mm	17"/432mm
24"/609mm	14"/355mm	15.5"/394mm	25"/635mm	4"/101mm	3.5"/89mm	19"/482mm

Part Number: 1525303

Invacare® Ltd, South Road,
Bridgend Industrial Estate, Bridgend.
CF31 3PY
Tel: 01656 647327.
Fax: 01656 649016
E-mail: ordersuk@invacare.com
Web: www.invacare.co.uk

Invacare® Ireland Limited
Unit 5 Seatown Business Campus
Seatown Road, Swords, County Dublin
Tel. (353) 18107084.
Fax: (353) 181 070 85
E-mail: ireland@invacare.com
Web: www.invacare.ie



- Crash tested to ISO16840- 4.
- Material fire tested to BS EN 1021-1 parts 1 & 2.



Yes, you can.

Freedom Designs - An Invacare® Company