

KÜSCHALL K-SERIES



Swiss exactitude

Precision is essential for küschall, and pushes us to bring more and more information to you. This is the reason why every küschall wheelchair portfolio includes the weight of each part. Thanks to these detailed data, it is easy to compare the weight of all options and make a safe choice.

Except any other indication, all weights indicated are based on seat width 400mm and seat depth 400mm. Options weights are calculated by pair when requested (e.g. for the wheels, push handles, brakes...) and with their related fixations.

All küschall options are presented in the portfolio, but offer can slightly vary according to each country. Contact your nearest küschall specialist for more information.

küschall K-Series

Content

- Frame 4
- Backrest 5
- Seating 7
- Footrest 8
- Side guards 9
- Cushions 10
- Colours 11
- Front wheels..... 12
- Rear wheels 13
- Brakes..... 17
- Safety options 18
- Accessories..... 19
- Küschall lifestyle collection..... 21

Frame

Frame material and angle

DDA0284 Aluminium frame,
DDA0285 frame angle 75° (K4) or 90° (AirLite) to the ground,
DDA0286 standard, adducted, seat depth tube +50mm
DDA0289
DDA0290
DDA0291

Total weight in the lightest configuration possible*:
approx. 7.7 kg



90° adducted frame

DDA0288 Carbon frame,
frame angle 75° (K4) to the ground

Total weight in the lightest configuration possible*:
approx. 7.5 kg



DDA0287 Titanium frame,
DDA0292 frame angle 75° (K4) or 90° (AirLite) to the ground

Total weight in the lightest configuration possible*:
approx. 7.7 kg





90° titanium frame

*All weights indicated are based on seat width 400mm and seat depth 400mm in the lightest and most compact configuration possible. Weights can change with the different configurations.

Backrest

Backrest frame

DDA0730	Aluminium	
	410 grams	
DDA0731	Titanium	
	290 grams	

Backrest upholstery




DDA0721	Light back upholstery	
	230 grams	
DDA0722	Standard black/grey backrest, tension adjustable with Velcro®	
	490 grams (including Velcro® straps)	
DDA0735	Aero backrest, tension adjustable with Velcro®	
	600 grams (including Velcro® straps)	
DDA0748	Alcantara black backrest, tension adjustable, with Velcro®	
	660 grams (including Velcro® straps)	
DDA0772	Without backrest upholstery, only for fixed backrest	

Push handles




DDA0750	Without push handles	
DDA0702	Push handles fixed, long	
	240 grams the pair (without back tubes)	
DDA0751	Push handles fixed, short	
	160 grams the pair (without back tubes)	
DDA0703	Push handles, height adjustable, rear set	
	720 grams the pair (incl. tubes and fixations)	
DDA0704	Push handles, foldable	
	560 grams the pair (incl. tubes and fixations)	
DDA0733	Push handles, height adjustable integrated	

Seating

Seat upholstery



DDA0610	Standard seat upholstery	
	220 grams	
DDA0615	Light seat upholstery	
	143 grams	
DDA0611	Tension adjustable with one central soft Velcro® band	
	560 grams	

Seat upholstery bags



DDA0616	Small bag (picture: on the right) To be fixed under the seat by Velcro® band	
	50 grams	
DDA0617	Big bag (picture: on the left) To be fixed under the seat by Velcro® band	
	70 grams	
DDA0606	Seat bag, click in and removable	
	130 grams	

Footrest



Footrest tube material

DDA1508	Aluminium	
	360 grams	
DDA1509	Titanium	
	210 grams	

Footrest tube

DDA1510	Standard footrest tube, rigid	
	370 grams	
DDA1513	Footrest cover	
	80 grams	

Footplate

DDA1570	Carbon footplate, angle adjustable	
	510 grams	
DDA1515	Footplate, flip-up to side	
	1'150 grams	

Side guards

Clothes guards / Mud guards

DDA1401 Carbon clothes-guards, flip back

Short: 420 grams the pair
Large: 540 grams the pair



DDA1408 Carbon mud guards, removable
(This option adds 2 cm to the total width of the wheelchair)

Short: 660 grams the pair
Large: 824 grams the pair



Siderests

DDA1416 Siderest height adjustable (with tools),
DDA1421 removable, short or long

Short: 1'580 grams the pair
Long: 1'700 grams the pair



DDA1411 Tubular armrest, height adjustable,
removable, swing-away, comfort
padded



Cushions

DDA0600

Seat cushion, 80 mm,
double foam with visco top

660 grams (SB38 cm*ST38 cm)



DDA0605















Alcantara cushion, black

460 grams (SB40 cm*ST36 cm)











Colours

Frame colours





Glossy black		Black matt	
Royal blue		Royal red	
Purple fire		Silver grey	
Umber glaze		Ivory white	
Smokey grey		Night shimmer	
Flashy yellow		Black grey rustic	
Mirror polished		Golden beach	

CPS-colours

Techno-anthracite		Comet red	
B-blue		Copper glaze	
Artic blue		Chrome	
Champagne		Olive	








The Colour Power Set (CPS) applies on the following parts and options: Starec & High Performance rear wheel hubs with quick release protection, 22" and 24" Starec rear wheel rims, Starec front wheels, forks, castor housing caps, front seat module caps, backrest adjustment plates.

Front wheels







DDA0800	Supersport	
	3": 180 grams the pair	
Various	Starec wheel, in various CPS colours	
	4": 364 grams the pair 5": 536 grams the pair	
DDA0860 DDA0863	Comfort wheel with black rubber	
	4": 480 grams the pair 5": 780 grams the pair	
DDA0810	Light wheel with black rubber	
	5": 624 grams the pair	
DDA0824	Skater wheel	
	4": 344 grams the pair	
DDA0801	Everyday wheel	
	4": 480 grams the pair	

Rear wheels









Rear wheel fixation and position

DDA0276	Vario-Ax, adjustable camber 1°, 3°, 7° and 10°	
	660 grams	
DDA0277	Standard axle (aluminium, black)	
	310 grams	
DDA0278	Carbon Axle	
	230 grams	
DDA0251	Rear wheel extension, 5 positions (with 2 axles on the picture)	
	800 grams the pair (without axles)	
DDA0248	Light rear wheel extension, 3 positions	
	416 grams the pair (without axles)	
DDA0630 DDA0631 DDA0632	Suspension	
	770 grams	
DDA0909	Tetra Clip	
	220 grams the pair	

Rear wheel type

DDA0915	küschall Starec wheel	
	24": 1'592 grams the pair 26": 1'700 grams the pair	
DDA0908 DDA0918	High Performance wheel (standard hub is black – if no CPS choice)	
	24": 1'356 grams the pair 25": 1'372 grams the pair	
DDA0913	Mountain Bike wheel	
	24": 1'712 grams the pair	
DDA0917	Spinergy wheel	
	24": 1'476 grams the pair 25": 1'508 grams the pair	
DDA0922	Ultra-light fibre wheel	
	24": 1'330 grams the pair	
DDA0921	Economy wheel	
	24": 3'400 grams the pair	




Rear wheels tyres

<p>DDA1001 Light wheel 7 bar</p> <p>24": 976 grams the pair 26": 1'000 grams the pair</p>	 <p>DDA1002 Profile tyre 7 bar</p> <p>24": 784 grams the pair</p> 
<p>DDA1007 Solid, light wheel</p> <p>24": 1'032 grams the pair</p>	 <p>DDA1004 Solid tyre, profile</p> <p>24": 1'140 grams the pair</p> 
<p>DDA1006 Solid tyre (KIK Type), black</p> <p>24": 904 grams the pair</p>	 <p>DDA1010 Marathon Evolution Plus</p> <p>24": 1'180 grams the pair 25": 1'235 grams the pair 26": 1'280 grams the pair</p> 
<p>DDA1011 Schwalbe speedrun</p> <p>24": 556 grams the pair 25": 564 grams the pair</p>	 <p>DDA1012, DDA1013 Tyres Schwalbe rightrun red or blue</p> <p>24": 960 grams the pair</p> 

Handrims

<p>DDA1102 Supergripp</p> <p>24": 1'160 grams the pair 26": 1'250 grams the pair</p>	 <p>DDA1104 Steel</p> <p>24": 1'460 grams the pair 26": 1'640 grams the pair</p> 
<p>DDA1106 Titanium</p> <p>24": 860 grams the pair</p>	 <p>DDA1113, DDA1114 Light Alu grey</p> <p>24": 630 grams the pair 25": 660 grams the pair 26": 760 grams the pair</p> 
<p>DDA1120 "The Surge" LT</p> <p>24": 1'020 grams the pair 25": 1'040 grams the pair</p>	 <p>DDA1117 Light Alu black</p> <p>24": 630 grams the pair 25": 660 grams the pair 26": 760 grams the pair</p> 


Spoke guards

DDA1200	Transparent	
	24": 660 grams the pair 26": 800 grams the pair	
DDA1207	«Original» Carbon	
	24": 720 grams the pair	
DDA1202	Black küsch	
	24": 680 grams the pair	

Brakes

DDA1300 Push brakes, levers bent

380 grams the pair



DDA1302 Performance brakes

420 grams the pair



DDA1303 Active brakes

300 grams the pair



DDA1307 Hub Break



DDA1310 Extension levers for standard brakes

Brakes weight doesn't change



Safety options

DDA1702
DDA1703
DDA1704

Antitipper, height adjustable, removable,
active and inactive position,
right, left or both sides

525 grams each



DDA1705

Active antitipper, active and inactive
position, removable

640 grams



DDA1710
DDA1711

Passive lights, assembled or enclosed

- Picture on left: assembled for Spinergy,
Mountain Bikes, High Performance wheels
- Picture on right: other rear wheel types
60 grams the pair



DDV1730

Seat belt

230 grams











DDA1770

Transportation kit incl. seat belt and labels



Accessories

DDA1744	Package holder	
	400 grams the pair	
DDA1745	Bottle holder	
	140 grams (including the bottle)	
DDA1746	Mobile phone holder	
	50 grams	
DDA1728	Wheel bag for 24", 25" and 26" rear wheels	
	750 grams	
DDA1700	Transit wheels 3"	
	820 grams the pair	
DDA1712 DDA1713	Küschall watch (mounted in castor housing), right or left side	
	30 grams	
DDA1715	Air pump	
	100 grams	
DDA1716	AirMan	

DDA1717

Tool kit



DDA1739

Front frame protectors (2 pieces)



80 grams the pair

Kuschall Lifestyle Collection

1531382 / 1531383 kuschall black t-shirt
1531463 / 1531464 Men or Women
1531189 / 1531379 (available in size S to XXL)
1531380 / 1531381



1513357 kuschall jacket
1513358 (available in S, M, L, XL, XXL)
1513359
1513360
1513361



1550289 kuschall sweat-shirt
1550290 (available in S, M, L, XL)
1550291
1550292



1553795 kuschall fleece bodywarmer
1553796 (available in S, M, L, XL, XXL)
1553797
1553798
1553799
1553800



1531187 kuschall baseball cap
black & silver



1486578 kuschall pen



1486577 kuschall carbon pen set



1486614 kuschall umbrella



1532047 kuschall backpack



1517781 kuschall scarf



1513362 Kuschall USB-Stick



