

# Invacare® | Pronto® M61™



## A fresh approach to mobility

The sleek and stylish *Invacare Pronto M61* is an affordable centre wheel drive power chair with 6 wheels on the ground for optimum safety and stability. It features *Invacare's Sure Step®* technology for smoother driving over lower thresholds. Designed around a compact base, the *Pronto M61* provides excellent manoeuvrability especially in confined spaces, and the electric seat provides 12 cm of elevation.



**Yes, you can.®**



Features and Options

### Optimum Comfort

The padded, shaped van seat ensures you are comfortable for day long use, while the popular configurable seating system offers a wide range of seat sizes allowing you to achieve an optimum fit. For easier transfers into or out of the **Pronto M61**, the seat will rotate 90 degrees in either direction, and the van seat armrest can flip up or alternatively the armrests can be removed altogether.



**Sure Step**® technology gives a smooth ride over low thresholds of up to 5cm.

### Upwardly mobile

With its centre wheel drive and its compact dimensions, the **Pronto M61** gives excellent manoeuvring in very small spaces and has a very tight turning circle. The simple to operate seat will electrically raise you up to 12 cm. In the elevated position it provides increased accessibility, which is ideal in an office environment or when you are out and about.



A range of seating options are available including the **Laguna** seat and back style.



A range of angle adjustable legrests and footplates available. (Configured version only)

### Technical data

<b>Pronto M61 Standard</b>	45,5 - 50 cm	41 - 46 cm	0 - 30°	89 - 94 cm
<b>Pronto M61 Configured</b>	33 - 43 cm	32 - 50 cm	0 - 12°	96 - 110 cm
<b>Pronto M61</b>	47 - 60 cm <sup>1</sup> 50 - 63 cm <sup>2</sup>	61 cm (excluding joystick)	68 cm <sup>3</sup> 97 - 102 cm <sup>4</sup>	R = 48 cm (front and rear turning radius) 5 cm
<b>Pronto M61</b>	9° / 16 %	2 x 30 Ah	2 x 180 W	16 km 6,4 km/h
		<b>Frame colours</b>		
<b>Pronto M61</b>	70 kg (91 kg with batteries)	136 kg		
			Deep blue pearl	Candy red pearl

- 1/ Seat height is measured including cushion, with 4° tilt angle
- 2/ Seat height is measured including cushion, with 9° tilt angle
- 3/ with footboard up
- 4/ with footboard down, with side mounted legrests

INVACARE® Limited  
Pencoed Technology Park - Pencoed -  
Bridgend CF35 5HZ  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: uk@invacare.com  
www.invacare.co.uk  
Sales Order E-mail:  
ordersuk@invacare.com

INVACARE® Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County  
Dublin  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: Ireland@invacare.com  
www.invacare.ie

© 2010 Invacare® International Sarl  
All rights reserved.

All information quoted is believed to be correct at time of print.  
Invacare® reserves the right to alter product specifications without prior consultation.

Pronto - UK - 11/2010  
1468372



Yes, you can.®