




PRODUCT INFORMATION UPDATE - February 2008

Power Wheelchairs - Mirage

Bulletin PU242
Date: 20/02/2008

Please find below details of the latest product changes for our Powerchair range.

Item	Description of change	New Part No. (if applicable)	Introduction date	Picture
MIRAGE	<p>In order to accommodate easier access to the Mirage tie down points, the front tie down bracketry has been moved from the inside of the frame to the outside. (See diagram).</p> <p>The serial number of the first chair with this change is :- 060812403448</p>		Feb-08	

If you require any additional information regarding any of the above changes, please contact

TECHNICAL SERVICES

TEL: 01656 647327

FAX: 01656 753299

EMAIL: uk@invacare.com