

rea |  **Azalea** /  **Azalea Assist**



Offers excellent seating comfort and many adjustment possibilities!

The Rea **Azalea** is the natural choice for the moderately active user as well as for the more passive one. The comfortable seating system allows the user to alter the seat position, ensuring maximum comfort at all times. This combined with good driving characteristics makes **Azalea** the ideal choice. **Azalea** is a self propelling wheelchair with large rear wheels, and is complimented with an alternative transit version that keeps the same frame but with small rear wheels called Rea **Azalea Assist**.





Features and Options

Tilt & recline makes life comfortable

The *Rea Azalea* and *Rea Azalea Assist* have a 21° or 25° seat tilt respectively and a 30° backrest recline. Both tilt and recline can be easily adjusted manually or electrically by the assistant and also by the user for the electrical version. With the increased seat tilt for the *Rea Azalea*

Assist an enhanced body positioning and pressure relief can be achieved for maximum comfort and rest. The backrest has a biomechanical positioning hinge, which follows the user's natural back movement when adjusting the backrest angle.



Seat size

The size of the seat is adjustable in width, depth and backrest height to meet individual needs.



Rear wheel adjustment
















The seat height and the balance can be adjusted individually.



Step less width adjustments

Individual width adjustment in the armrests and leg rests.

Technical data

<i>Rea Azalea</i>	 39-55 cm	 43 - 50 cm	 40/45 cm	 50 - 80 cm (without seat cushion)	 24 - 36 cm	 33 - 50 cm	 -1°-20° (-1°-25°)
<i>Rea Azalea</i>	 0°-30°	 SW + 22 cm (SW + 23 cm)	 90-125 cm	 95-102 cm	 34 kg (SW 44 kg)		
<i>Rea Azalea</i>	 Max 135 kg	 Transport weight: 16,5 kg	 Crash-tested ISO 7176-19	<p>The <i>Rea Azalea</i> complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: www.invacare-rea.com under Rea wheelchairs/Technical downloads.</p>			

Invacare International Sarl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeuropa@invacare.com
www.invacare.eu.com

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Rea Azalea - EU - 08/2010

Frame colour

Upholstery



Pearl Grey



Grey Dartex TR23



Grey plush TR 32



Rea, Rea design and DSS (Dual Stability System) design are registered trademarks of Invacare® International.



Yes, you can.®