



Yes, you can.®

➤ *Invacare®* **Reliant Mobile Hoists**

Multi-functional mobile hoists

Invacare® **Reliant Mobile Hoists** feature 2 models; the **Reliant 150** and **Reliant 250**, with a maximum lifting capacity of 130 kg and 175 kg. These models have been designed for home and residential use and have features that can address carer and client needs in terms of function, reliability, manoeuvrability and safety.



Warranty

Actuator	5 years
Chassis	3 years
Battery	1 year

Safe ergonomic and comfortable transfers

The *Invacare® Reliant Mobile Hoists* ensure comfort and safety during transfers. The foot operated, adjustable base facilitates ease of access around a bed or wheelchair during the transfer process. The spreader bar attachment allows the spreader bar to be used with most types of sling*. Transfers can be carried out by using a push button control to lift or lower the client to the desired height, and the large push handle allows ease of handling and control.

* Dependent on a documented risk assessment of the patient to be transferred, Plastic clip style slings can not be used.



Safety and security features

For added patient and carer safety, both mobile hoists feature a number of emergency systems, such as: automatic stop function that stops the lift in the event of the lift boom coming into contact with an obstacle. An emergency stop button, an integrated manual switch and an emergency mechanical descent system, if the remote control should fail. The multi-caster base ensures the hoist is easily manoeuvrable, even on carpeted surfaces. Both hoists can be dismantled for transport and storage.

Features



Adjustable base width, with click push brackets
Rear castors have push lock brakes.



Push-bar handle
Large handle for easier handling.



Removable battery
For ease of charging. Allows 24h use.

Technical data

Product Part Number
REL150UKPLUS
REL250UKPLUS

Product Description
Reliant 150 Mobile Hoist
Reliant 250 Mobile Hoist

Colour



Silver

Reliant 150	154 cm	38 cm	91 cm	62 cm	100 cm
Reliant 250	167 cm	40 cm	89 cm	62 cm	107 cm
Reliant 150	11 cm	40 cycles	26 kg	max. 130 kg	
Reliant 250	9 cm	40 cycles	35 kg (21 + 14)	max. 175 kg	