




Testing performed at:	DEMKO[®] Testing and Certification Lyskaer 8 DK-2730 Herlev, Denmark
-----------------------	--


Report No.:	126595/EMC-DK 99-03096	No. of pages in testreport - excl. this page	43
Applicant:	Scandinavian Mobility EC-Høng A/S		
Address:	Østergade 3, DK-4270, Høng, Denmark.		

Test item:	Nursing Bed		
Input Voltage: [Volt]	230V ~	Rated power: [Watt / VA]	70 VA
Rated Current: [Amp]		Frequency: [Hz]	50Hz
Manufacturer (Resp.):	Scandinavian Mobility EC-Høng A/S, Østergade 3, DK-4270 Høng, Denmark.		
Trademark:	Scan Bed		
Type No.	SB 750 & SB 400		

Date of receipt:	1999-09-01		
Date(s) of testing	1999-09-02	1999-09-06/07	1999-09-09

Test specification(s):	EN 60601-1-2:93 (EN 55011:91 + A1:97 + A2:96, IEC 801-3:84, EN 60801-2:93, IEC 801-4:88, EN 61000-4-5:95)
------------------------	---

Test verdict:	All tested parameters fulfil the requirements.	<small>Copy Protection</small> 
Special remarks:	The test result of type SB 750 also cover type SB 400 due to identical electrical design.	

Date:	1999-10-20	Title:	Project Engineer
Signature:		Name:	Preben Jakobsen

The test report shall not be reproduced except in full without the written approval of DEMKO[®]. The testing is only valid for the equipment tested.

DEMKO1998-1*