

# Aquatec® | Shower stools, shower chairs



Aquatec **Dot**

- Round shape, fits into every shower
- Easily disassembled: ideal for taking on travels
- Textured seat surface for a safe seat
- Rubber feet provide a secure stand




Aquatec **Disk on Dot**

- Within seconds, the **Disk**, the shower stool **Dot** and an adapter are combined to a shower stool with swivel seat
- High load capacity: 125 kg

## Technical data

Order numbers	Aquatec <b>Dot</b>	Aquatec <b>Disk on Dot</b> Seat height 55 cm	Aquatec <b>Disk on Dot</b> Seat height 50 cm
Shower stool round, seat height 50 cm	<b>10200-32</b>		
Shower stool round, seat height 55 cm	<b>10201-32</b>		
Blue Disk + Dot + Adapter		<b>9.01.002</b>	<b>9.00.002</b>
Grey Disk + Dot + Adapter		<b>9.01.008</b>	<b>9.00.008</b>
Dimensions			
Seat width	40 cm (A)	40 cm (A)	40 cm (A)
Seat depth	40 cm (A)	40 cm (A)	40 cm (A)
Seat height	50 / 55 cm (B)	57.5 cm (B)	52.5 cm (B)
Total width	46 cm (C)	46 cm (C)	46 cm (C)
Total depth	46 cm (C)	46 cm (C)	46 cm (C)
Total height	50 / 55 cm (B)	57.5 cm (B)	52.5 cm (B)
More info			
Weight	2.5 kg	4.1 kg	4.1 kg
Load capacity	125 kg	125 kg	125 kg
Colour	White	Blue, grey	Blue, grey
Accessories			
<b>Disk</b> , blue	4.02.002	-	-
<b>Disk</b> , grey	4.02.008	-	-
Adapter for <b>Disk on Dot</b>	15326AQT	-	-



**AQUATEC®**



**Yes, you can.®**



**Invacare® I-Fit 9780E**

- Height Adjustable from 38 cm to 55 cm
- Large ergonomic shaped seat for comfort
- Easy to clean
- Light weight and easy to transport
- Stylish modular design
- Hygiene recess




**Invacare I-Fit 9781E**

- Height Adjustable from 38 cm to 55 cm
- Removable backrest
- Ergonomic shaped seat with backrest for extra comfort
- Easy to clean
- Light weight and easy to transport


### Technical data

Order numbers	Invacare I-Fit 9780E	Invacare I-Fit 9781E
Shower stool	<b>9780E</b>	
Shower chair		<b>9781E</b>
<b>Dimensions</b>		
Seat width	34	34
Seat depth	37	37
Seat height	38-55	38-55
Total width	50	50
Total depth	43	43
Total height	38-55	76-93
<b>More info</b>		
Weight	3,7 kg	4,5 kg
Load capacity	180 kg	180 kg
Colour	white	white

Invacare I-Fit 9780E



Invacare I-Fit 9781E





**Aquatec Tellus**

- Stable, stackable shower stool
- Ergonomic seat surface, can be supplied with or without hygiene recess
- Height-adjustable, corrosion-protected legs
- Accessories: with the soft backrest, the shower stool can be converted into a shower chair



**Invacare® H2240/1**

- Stable shower stool
- Large ergonomic seat surface with hygiene recess makes personal hygiene easier
- Fixed armrests
- Height-adjustable, corrosion-protected legs



**Aquatec Pluto**

- Ergonomic seat surface, can be supplied with or without hygiene recess
- Height-adjustable, corrosion-protected legs
- Non-slip rubber caps (FlexiFoot) compensate height differences of up to 1 cm
- Accessories: with the soft backrest, the shower stool can be converted into a shower chair



**Aquatec Galaxy**

- Folding, height-adjustable shower chair, perfect for taking with you while travelling and compact for storage
- Ergonomic seat surface, can be supplied with or without hygiene recess
- Soft backrest for additional support


## Technical data

Order numbers	Aquatec <b>Tellus</b>	Invacare® <b>H2240/1</b>	Aquatec <b>Pluto</b>	Aquatec <b>Galaxy</b>
Shower stool without hygiene recess	<b>1470796</b>			
5 pack	<b>1470797</b>			
Shower stool with hygiene recess	<b>1470798</b>			
5 pack	<b>1470799</b>			
Accessories: Backrest, can be retrofitted.	<b>1470800</b>		<b>1470800</b>	
Shower stool with armrests and hygiene recess		<b>H2240/1</b>	<b>1470801</b>	
5 pack			<b>1470804</b>	
Shower stool with armrests, without hygiene recess			<b>1470802</b>	
5 pack			<b>1470803</b>	
Shower chair with backrest and armrests, with hygiene recess				<b>1470805</b>
Shower chair with backrest and armrests, without hygiene recess				<b>1470806</b>


Dimensions				
Seat width	40 cm (A)	41 cm (A)	40 cm (A)	40 cm (A)
Seat depth	34.5 cm (B)	40 cm (B)	34.5 cm (B)	34.5 cm (B)
Seat height	46 - 61 cm (C)	47 - 55 cm (C)	46 - 61 cm (C)	47 - 62 cm (C)
Total width	55 cm (D)	50 cm (D)	55 cm (D)	58 cm (D)
Total depth	51 cm (E)	45 cm (E)	51 cm (E)	46 cm (E)
Total height	46 - 61 cm (C)	58 - 66 cm (F)	63 - 78 cm (F)	84 - 99 cm (F)
Graduation, height adjustment	2.5 cm	2 cm	2.5 cm	2.5 cm
Height to armrests	–	–	–	63 - 78 cm (G)
Width between the armrests	–	45 cm (G)	49 cm (G)	51 cm (H)
Width, hygiene recess	–	14 cm	13 cm	13 cm
Depth, hygiene recess	–	17 cm	14 cm	14 cm

More info				
Weight	3.8 kg	5.7 kg	4.5 kg	7.5 kg
Load capacity	125 kg	120 kg	125 kg	100 kg
Colour	White	White	White	White


Aquatec **Tellus**




Invacare **H2240/1**



Aquatec **Pluto**



Aquatec **Galaxy**





Aquatec **Pico**

- Modular shower chair system
- Four configurations
- Textured seat surface for safe sitting
- Large hygiene recess, either at the front or back, due to two mounting options for the backrest
- Height-adjustable legs
- Rubber feet provide a secure stand
- Simple and easy assembly, no tools required
- Load capacity of up to 160 kg

### Options and features



**Pico**

**Pico with armrests**

**Pico with backrest**



### Technical data

Order numbers	Aquatec <b>Pico</b>
Shower stool	<b>9.20.001</b>
Shower stool with armrests	<b>9.21.001</b>
Shower chair with backrest	<b>9.23.001</b>
Shower chair complete	<b>9.22.001</b>
Backrest	<b>16469</b>
Armrests (pair)	<b>16468</b>

Dimensions	
Seat width	43 cm (A)
Seat depth	42 cm (B)
Seat height	42,5 - 57,5 cm (C)
Total width	57,5 cm (D)
Total depth	52 cm (E)
Total height	82,5 - 97,5 cm (F)
Height to armrests	60,5 - 75,5 cm (G)
Width between the armrests	43 cm (H)
Width, hygiene recess	18 cm
Depth, hygiene recess	14,5 cm
Graduation, height adjustment	2,5 cm

More info	
Weight, shower stool	3,2 kg
Weight, shower stool with armrests	4,1 kg
Weight, shower chair with backrest	4,0 kg
Weight, shower chair complete	4,8 kg
Load capacity	160 kg!
Colour	White

Aquatec **Pico**





Aquatec **Pico 3 in 1**

- Stationary, height adjustable multifunctional shower chair
- Versatile: shower chair, commode or commode support frame
- Easy handling: toilet pan can be removed from all four sides
- Absolutely corrosion-free
- Armrests, backrest, toilet pan and toilet pan holder are included in the scope of supply
- Seat height adjustable in 6 steps from 42.5 – 57.5 cm
- Very high ILoad capacity of 160 kg

### Options and features



**Pico 3 in 1** with one armrest



**Pico 3 in 1** as commode support frame



**Pico 3 in 1** with universal soft seat



### Technical data

Order numbers	Aquatec <b>Pico 3 in 1</b>
Shower chair complete	<b>1525887</b>
Universal soft seat	<b>1526032</b>

Dimensions	
Seat width	43 cm (A)
Seat depth	42 cm (B)
Seat height	42,5 - 57,5 cm (C)
Total width	57,5 cm (D)
Total depth	52 cm (E)
Total height	82,5 - 97,5 cm (F)
Height to armrests	60,5 - 75,5 cm (G)
Width between the armrests	43 cm (H)
Graduation, height adjustment	2,5 cm
Height, ground to lower seat edge:	
with holder for toilet pan	34 - 49 cm
without holder for toilet pan	37,5 - 52,5 cm
Width, cut out in seat plate	21,5 cm
Depth, cut out in seat plate	25 cm
Graduation, height adjustment	2,5 cm

More info	
Weight, shower chair complete	4,8 kg
Load capacity	160 kg!
Colour	White

Aquatec **Pico 3 in 1**



Invacare International Sarl  
Route de Cité-Ouest 2 - 1196 Gland  
Switzerland  
Tel: +41 (0)22 354 60 10  
Fax: +41 (0)22 354 60 11  
hqeurope@invacare.com  
www.invacare.eu.com

© 2010 Invacare® International Sarl  
All rights reserved.

All information quoted is believed to be correct at time of print.  
Invacare® reserves the right to alter product specifications without prior consultation.

Shower stools - EU - 09/2010



**Yes, you can.®**