











XLT DYNAMIC

														
SHf		SHr						L=75	L=110	L=150	mm			
		1	2	3	4	5	6	*						
22"														
39					39	38			2			75 (3")		
40				40	39	38			1			75 (3")		
41				40	39	38			2			100 (4")		
42			41-42	40	39	38			1			100 (4")		
43		43	41-42	40	39	38				2		75 (3")		
43		43	41-42	40	39	38				4		125 (5")		
43		43	41-42	40	39	38				3		100 (4")		
43		43	41-42	40	39	38	* D3 (SD40)		1			125 (5")		
44		43	41-42	40	39	38				1		75 (3")		
44		43	41-42	40	39	38	* D2 (SD40)			3		125 (5")		
44	44	43	41-42	40	39	38				2		100 (4")		
45	44	43	41-42	40	39	38				1		100 (4")		
45	44	43	41-42	40	39	38	* D1 (SD40)			2		125 (5")		
46	44	43	41-42	40			* D1 (SD40)			1		125 (5")		
47	44	43	41-42	40							2	75 (3")		
47	44	43	41-42	40							3	100 (4")		
47	44	43	41-42								4	125 (5")		
48	44	43	41-42								1	75 (3")		
48	44	43	41-42								2	100 (4")		
48	44	43	41-42								3	125 (5")		
49	44	43	41-42								1	100 (4")		
49	44	43	41-42								2	125 (5")		
50	44	43									1	125 (5")		
24" / 25" Low profile														
40						40			1			75 (3")		
41						40			2			100 (4")		
42					41-42	40			1			100 (4")		
43				43	41-42	40	* D5 (SD40)			2		75 (3")		
43				43	41-42	40	* D1-6 (SD40)			4		125 (5")		
43				43	41-42	40	* D5 (SD40)			3		100 (4")		
43				43	41-42	40	* D1-6 (SD40)		1			125 (5")		
44			44	43	41-42	40	* D4-5 (SD40)			1		75 (3")		
44			44	43	41-42	40	* D1-6 (SD40)			3		125 (5")		
44			44	43	41-42	40	* D4 (SD40)			2		100 (4")		
45		45	44	43	41-42	40	* D3-4 (SD40)			1		100 (4")		
45		45	44	43	41-42	40	* D3-6 (SD40)			2		125 (5")		
46	46	45	44	43	41-42	40	* D2-6 (SD40)			1		125 (5")		
47	46-47	45	44	43	41-42	40	* D2-4 (SD40)				2	75 (3")		
47	46-47	45	44	43	41-42	40	* D2-4 (SD40)				3	100 (4")		
47	46-47	45	44	43	41-42	40	* D2-5 (SD40)				4	125 (5")		
48	46-47	45	44	43	41-42	40	* D1-4 (SD40)				1	75 (3")		
48	46-47	45	44	43	41-42	40	* D1-4 (SD40)				2	100 (4")		
48	46-47	45	44	43	41-42	40	* D1-5 (SD40)				3	125 (5")		
49	46-47	45	44	43	41-42		* D1-4 (SD40)				1	100 (4")		
49	46-47	45	44	43	41-42		* D1-4 (SD40)				2	125 (5")		
50	46-47	45	44	43			* D1-6 (SD40)				1	125 (5")		

* = Not applicable
SD = seat depth

XLT DYNAMIC

											
SHf	SHr						L=75	L=110	L=150	mm	
	1	2	3	4	5	6	*				
26" Low profile											
42						41	* D1-6 (SD40)	1		100 (4")	
43					42-43	41	* D1-6 (SD40)		2	75 (3")	
43					42-43	41	* D1-6 (SD40)		4	125 (5")	
43					42-43	41	* D1-6 (SD40)		3	100 (4")	
43					42-43	41	* D1-6 (SD40)	1		125 (5")	
44				44	42-43	41	* D1-6 (SD40)		1	75 (3")	
44					42-43	41	* D1-6 (SD40)		3	125 (5")	
44				44	42-43	41	* D1-6 (SD40)		2	100 (4")	
45			45	44	42-43	41	* D1-6 (SD40)		1	100 (4")	
45			45	44	42-43	41	* D1-6 (SD40)		2	125 (5")	
46		46	45	44	42-43	41	* D1-6 (SD40)		1	125 (5")	
47	47	46	45	44	42-43	41	* D1-6 (SD40)			2 75 (3")	
47	47	46	45	44	42-43	41	* D1-6 (SD40)			3 100 (4")	
47	47	46	45	44	42-43	41	* D1-6 (SD40)			4 125 (5")	
48	47-48	46	45	44	42-43	41	* D1-6 (SD40)			1 75 (3")	
48	47-48	46	45	44	42-43	41	* D1-6 (SD40)			2 100 (4")	
48	47-48	46	45	44	42-43	41	* D1-6 (SD40)			3 125 (5")	
49	47-48	46	45	44	42-43	41	* D1-6 (SD40)			1 100 (4")	
49	47-48	46	45	44	42-43	41	* D1-6 (SD40)			2 125 (5")	
50	47-48	46	45	44	42-43		* D1-6 (SD40)			1 125 (5")	

* = Not applicable
SD = seat depth