

Invacare® | **Softform® Odstock Wedge**



Pressure reducing wedge

The *Invacare Softform Odstock Wedge* was developed in partnership with clinicians to assist in the reduction of interface pressures on the sacrum and heels without significantly increasing pressure on the thighs or calves. It is designed to improve and maintain patient posture in the recumbent position, and to reduce the tendency for the client to slip down the bed, thereby reducing shearing forces on the sacrum. The **Odstock Wedge** is ideal for improving posture for those with back discomfort.

Softform®



Yes, you can.®



Features and Options

Reduces pressure on the heels

Clinical research has effectively demonstrated that the degree of knee flexion necessary to achieve substantial levels of pressure relief is relatively small. The *Invacare Softform Odstock Wedge* was designed by clinicians to reduce pressure on the heels through re-positioning the hip and elevating the knee.

Reduces lumbar lordosis

By tilting the pelvis, the **Odstock Wedge** also reduces lumbar lordosis – desirable following back surgery and ideal for improving the sleeping posture of patients with back problems.



High frequency welded seams

Ensures the core is protected from fluid ingress, increasing product integrity.



High quality construction

The two-way stretch cover, is waterproof and vapour-permeable. The foam core is constructed from, high density foam.



Designed to be used in conjunction with the *Invacare Heelpad*

Used in conjunction, these two products combine to significantly reduce pressure on the vulnerable heel area.

Technical data



Width: 570 mm



Length: 865 mm



Height: 115 mm



Weight: 1.5 kg



Cleaning: 80°

Softform Odstock Wedge

Refer to instructions for use and relevant infection control guidelines.

Manufactured to European Standards:
Interior Foam: BS 5852: Crib 5,
Cover: BS7175: Crib 5

Colours



Fuchsia (Standard)



Blue (Option)

Warranty: 1 year

