

Spare parts **AQUATEC ORCA**

Pos.	Bezeichnung	Article#
1	Suction cups with thread, 2 tags (pair)	15071
2	Suction cups with toggle, 1 tag (pair)	15072
3	Bolt, ø 10 mm	16679
4	Bolt for gliding shoe	15105
5	Suction cup lid	15121
6	Thrust washer for scissors	15344
7	Flat spring	15415
8	Gliding shoe	16446
9	Bushing latch	16448
10	Latch	16449
11	Holding device (excluding screws)	16455
12	Bolt right	1470579
13	Screws phone jack , 40x22	15348
14	Lens head screw DIN 7985, M5 x 30	16781
15	Nut for actuator fastening	16779
16	Pipe clamp for actuator fastening	16778
17	Washer A2 DIN 9021	16780
18	Axle	16793
19	Locking left	16791
20	Locking right	16792
21	Pressure spring for locking	16794
22	Pulley	16450
23	Retaining ring A2	16456
24	Actuator	16782
25	Cover seat, white	16606
26	Cover backrest, white	16604
27	Battery charger	16429
28	Hand control complete	1470623
29	Sideflap with filmhinges	15892
	Cover set seat and backrest, white	1470576
	Suction cup set, ø 8 cm (pair with toggle, 1 tag/ pair with thread, 2 tags)	15353
12	Bolt, left	1470580
	User manual	16451

Exploded view **AQUATEC ORCA**

