


Customer :
 Customer No. :
 Phone :

Contact :
 Commission :
 Serial No. :




Item	Qty	Part-No.	Part-Designation	
-	1		SWING OUT FRONT LIGHT SYSTEM 24 V, CPL.	
1	1	1420413	HEADLAMP LH 24V CPL.	
2	1	1420412	HEADLAMP RH 24V CPL.	
1	1	SP1420413	HEADLAMP LH 12V CPL. (10km/h)	->05/02
2	1	SP1420412	HEADLAMP RH 12V CPL. (10km/h)	->05/02
3	2	P31108	LIGHTBULB 24V/10W	
4	2	1419271	LIGHTBULB 24V/21W, YELLOW	
3	2	1423458	LIGHTBULB 12V/10W	->05/02
4	2	1423457	LIGHTBULB 12V/21W, YELLOW	->05/02
15	1	1433737	HEADLIGHT LENS	
5	4	F20110	NYLOC NUT M10 DIN 985	
6	4		SELF-TAPPING CROSS HEAD SCREW 4,8 X 16	
8	4		BOLT M5 X 12 DIN 7500	
9	4		BOLT M10 X 75 DIN 7380	
11	1	1418547-9005	BEARING PLATE LH	
12	1	1418544-9005	BEARING PLATE RH	
13	1	1428082	SWIVEL MECHANISM RH	
14	1	1428081	SWIVEL MECHANISM LH	
16	1	SP1422822	FIXING SET FRONT HEADLIGHT FOR SB AB 43CM	
10	2		SPACER	
7	4		BOLT 4,8 X 25 DIN 7981	

Storm 3



Customer :
 Customer No. :
 Phone :

Contact :
 Commission :
 Serial No. :




Item	Qty	Part-No.	Part-Designation	
-	1	SP1419335	LIGHT SUPPORT WITHOUT LIGHT & LUGGAGE CARRIER CPL.	
-	1	SP1419336	LIGHT SUPPORT WITH LIGHT WITHOUT LUGGAGE CARRIER CPL.	
-	1	SP1419337	LIGHT SUPPORT WITHOUT LIGHT FOR LUGGAGE CARRIER CPL.	
-	1	SP1419338	LIGHT SUPPORT WITH LIGHT FOR LUGGAGE CARRIER CPL.	
1	1	1419274	REFLECTOR, RED	
2	2		BOLT M5x16 DIN 7380 SV	
3	1	1419333	REAR LIGHT SUPPORT LH	
3	1	1419334	REAR LIGHT SUPPORT RH	
3	1	1421815	LIGHT SUPPORT LH WITH LUGGAGE CARRIER CUTOUT	
3	1	1421814	LIGHT SUPPORT RH WITH LUGGAGE CARRIER CUTOUT	
4	2		GETO WELL NUT M5	
5	1		WASHER 5,3 X 15 X 1,6	
6	1	1419331	LAMP CASING, CLEAR LH	
6	1	1419332	LAMP CASING, CLEAR RH	
6	1	1421738	LAMP CASING, SILVER LH	->06/03
-	1	SP1421738	LAMP CASING, NO LIGHTS	07/03->
6	1	1421739	LAMP CASING, SILVER RH	->06/03
-	1	SP1421739	LAMP CASING, NO LIGHTS	07/03->
7	1	1419935	REAR LIGHT 24V CPL.	
8	1	1419271	LIGHTBULB 24V/21W, YELLOW	
9	1	P31191	LAMP 24V/5W	
-	1	SP1419935	REAR LIGHT 12V CPL. (10km/h)	->05/02
7	1	1419935	REAR LIGHT 24V CPL.	
8	1	1423457	LIGHTBULB 12V/21W, YELLOW	->05/02
9	1	1423459	LIGHTBULB 12V/5W (SOFITTE)	->05/02

Storm 3



Customer :
 Customer No. :
 Phone :

Contact :
 Commission :
 Serial No. :



Item	Qty	Part-No.	Part-Designation
-	-		ONLY FOR STORM³ TRUE TRACK
-	1	SP1433547	LIGHT SUPPORT WITHOUT LIGHT & LUGGAGE CARRIER
-	1	SP1433556	LIGHT SUPPORT WITH LIGHT WITHOUT LUGGAGE CARRIER
-	1	SP1433557	LIGHT SUPPORT WITHOUT LIGHT FOR LUGGAGE CARRIER
-	1	SP1433558	LIGHT SUPPORT WITH LIGHT FOR LUGGAGE CARRIER
1	2	1419935	REAR LIGHT 24V CPL.
2	2	1419274	REFLECTOR, RED
3	2	1432553	WARNING STICKER FOR RIM GEARLESS MOTORS
4	4		HEX NUT M5
5	2		HEX HEAD SCREW M6 X 12
6	4		BOLT M5x16 DIN 7380 SV
7	1	1429342	LIGHT SUPPORT LH
7	1	1430143	LIGHT SUPPORT LH WITH LUGGAGE CARRIER CUTOUT
8	1	1429343	LIGHT SUPPORT RH
8	1	1430144	LIGHT SUPPORT RH WITH LUGGAGE CARRIER CUTOUT
9	4	1422052	GETO WELL NUT M5
10	2	1419337	RIVET
11	2		WASHER 6,4 X 12,5 X 1,6
12	6		WASHER 5,3 X 15 X 1,6
13	1	1419331	LAMP CASING, CLEAR LH
14	1	1419332	LAMP CASING, CLEAR RH
15	1	1419271	LIGHTBULB 24V/21W, YELLOW
16	1	P31191	LAMP 24V/5W

Storm 3



Release
 07 / 2004
 Chapter 6
 Sheet 3