

Invacare® | **Stratos Pro**



Simply faster

The new *Invacare Stratos Pro* is a unique new piston aerosol compressor presented in a premium package. Combined with **Medjet Pro** nebuliser and accessories, it provides fast and effective long-term aerosol treatment. Compact, reliable and simple, **Stratos Pro** is economical and easy to use.



Yes, you can.®



Options and Features

Faster and more effective

The *Invacare Stratos Pro* patented one-block high-pressure compressor pump design prevents flow losses thus enhancing the system's efficiency. Coupled to the *Medjet Pro*, it nebulises 2 ml of NaCl 0,9%: in approximately three minutes.

Flexible and easy to clean

Both *Stratos Pro* and *Medjet Pro* can nebulise a large array of medications down to <5 µm droplets. The inhalation valve is activated with minimal respiratory flow, making *Medjet Pro* the ideal choice for children. Whether for home or hospital use *Medjet Pro* can be easily sterilised in boiling water, with cold disinfecting liquids or with an autoclave.



Easy to carry

Lightweight premium design with integrated handle and carrying bag.

Efficient

Single block compressor provides efficient continuous nebulisation and short inhalation times.

Optimum Performance

- Residual volume is reduced down to 0,5 ml.
- Inhalation and exhalation valve reduce drug loss
- Up to 15 ml easy-to-read graduated scale

Technical data



17,2 cm



10,1 cm



17,3 cm



1,7 kg



230 V, 50 Hz

Stratos Pro

Medjet Pro

MMAD: 4,25 µm
 Mass Fraction < 5 µm: 55,1 %
 Nebulisation Rate: 0,55 ml/min average

Invacare International Sarl
 Route de Cité-Ouest 2 - 1196 Gland
 Switzerland
 Tel: +41 (0)22 354 60 10
 Fax: +41 (0)22 354 60 11
 hqeuropa@invacare.com
 www.invacare.eu.com

© 2010 Invacare® International Sarl
 All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Stratos Pro - EU - 08/2010



Yes, you can.®