

REA DAHLIA

Altezze possibili in combinazione con diametro ruote

REA Dahlia 30°

Seat Hight	Rear Wheel Size	Seat Tilt as delivered	Backrest recline ad delivered
330 mm	16"	30	30
	20"	10	30
350 mm	16"	30	30
	20"	20	20
	22"	10	30
380 mm	20"	20	10
400 mm	20"	30	20
	22"	20	10
430 mm	22"	30	20
	24"	20	10
450 mm	16"	30	30
	24"	30	20
480 mm	16"	30	30
500 mm	20"	30	30

REA Dahlia 45°

Seat Hight	Rear Wheel Size	Seat Tilt as delivered	Backrest recline ad delivered
350 mm	16"	45	30
	20"	30	30
380 mm	16"	45	30
	20"	45	10
	22"	30	30
400 mm	20"	45	30
430 mm	20"	45	30
	22"	45	30
450 mm	22"	45	30
	24"	45	30
480 mm	16"	45	30
	24"	45	30
500 mm	16"	45	30