


# KÜSCHALL COMPACT

Technical Sheet

Legenda misure



Compatibilità freni

Scissor brakes	Seat depth	Rear wheel position					
		1	2	3	4	5	6
ST32-38		Red	Red	Red	Red	Green	Green
ST40		Red	Red	Red	Red	Green	Green
ST42		Red	Red	Red	Green	Green	Green
ST44		Red	Red	Green	Green	Green	Green
ST46		Red	Green	Green	Green	Green	Green
ST48		Green	Green	Green	Green	Green	Green
ST50		Green	Green	Green	Green	Green	Green

Posizionamento montaggio ruotine anteriori

Front wheel / Front wheel attachment	High mounted	Standard mounted	Low mounted
4"	41 - 48	44 - 51	47 - 53
5"	44 - 49	47 - 52	50 - 53
6"	46 - 51	49 - 53	52 - 53
7"	48 - 52	51 - 53	54 - 55*