

INVACARE® TOP END®

Force® CC



Yes, you can.®

INVACARE® TOP END® Force® CC



The *Invacare Top End Force CC* is a handcycle you can ride on the trails and is lightweight yet robust with extreme climbing gears, sturdy 559-26" bike wheels equipped with mountain bike tires and a stop-on-a-dime front disk brake!

The design is a cross between our popular *Invacare Top End Force RX* handcycle and the *Invacare Top End Force G* handcycle with enhanced ground clearance and aggressive back angle adjustments that put you in charge of those hills, rocks and rough terrain.

Now you can force the pace....cross country or off road!

FEATURES AND OPTIONS



Front disk brake standard.



Optional dual rear disk brakes.



Brake lever for optional rear disk brakes.



Welded footrest clamp.



TECHNICAL DATA



15" or 17"
(approx. 381 or 432 mm)



Wheelbase:
60"
(approx. 1524 mm)



Front and rear:
26"



Max. 250 lbs
(approx. 113 kg)

Please refer to order form for all technical specifications.