



Product Service

CERTIFICATE

No. Z1 10 07 32457 011

Holder of Certificate: INVACARE France Operations SAS

Route de St. Roch
37230 Fondettes
FRANCE

Production Facility(ies): 32457

Certification Mark:



Product: Beds for disabled persons

Model(s): Alegio

Parameters:	dimensions (length x width)	2100 x 1010 mm
	mattress support platform	2000 x 900mm
	height adjustment	365 to 810 mm
	voltage supply	230 V ±10%, 50 Hz
	max. current input	0,8 A
	voltage output	24 V
	operating time	10% max . 6 minutes/hour
	applied part	type B
	degree of protection	IP 66
	protection class	II
	sound level	45 dB (A)
	safe working load	170 kg
	max. patient weight	135 kg
	angle A	70°
	angle B	21°
	angle C	90°
	angle D	17°

Tested according to: EN 1970/A1:2005

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 71362564

Date, 2010-07-29

(Torsten Zimmer)



Page 1 of 1