



Product Service

CERTIFICATE

No. Z1 09 11 52355 016

Holder of Certificate: **Küschall AG**
Benkenstrasse 260
4108 Witterswil
SWITZERLAND

Production Facility(ies): 52355

Certification Mark:



Product: **Active wheelchairs**

Model(s): **Ultra-Light / Compact**

Parameters:	Ultra-Light	Compact
Total length incl. footrests	730-1100 mm	740-1160 mm
Total length without footrests	----	570-850 mm
Wheelchair width	490-750 mm	490-750 mm

- Further parameters see enclosure

Tested according to: DIN EN 12182:1999
DIN EN 12183:1999

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 71358709

Date, 2009-11-25 (Torsten Zimmer)
Page 1 of 2





Product Service

Enclosure to Certificate No.: Z1 09 11 52355 016

Model designation for the active wheelchairs Ultra-Light / Compact

	Ultra-Light	Compact
Total length incl. foot rests	730 - 1100 mm	740 - 1160 mm
Total length without foot rests	----	570 - 850 mm
Wheelchair width	490 - 750 mm	
Wheelchair height	690 - 1010 mm	670 - 1010 mm
Seat width	280, 300, 320, 340, 360, 380, 400, 420 mm (Single diagonal strut) 440, 460, 480, 500 mm (Double diagonal strut)	
Seat depth	320 - 500 mm	
Seat height at the front edge	430 - 530 mm	370 - 530 mm
Seat angle	-2° - 16°	
Width of the backrest	280 - 500 mm	
Height of the backrest	300 - 510 mm	
Angle of the backrest	82° - 102°	
Max. user weight	100 kg (Single diagonal strut) 130 kg (Double diagonal strut)	
Driving wheels		
- Size	22, 24, 25 and 26 x 1" - 1 3/4"	
- Tyres	<input checked="" type="checkbox"/> air-filled	<input checked="" type="checkbox"/> accident proof
- Adjustability	<input checked="" type="checkbox"/> yes	<input type="checkbox"/> no
- Kind of fixation	<input checked="" type="checkbox"/> quick release	<input type="checkbox"/> rigid
- Adjustment of the camber	0° / 3°	
Front wheels		
- Size	3, 4, 5, 6 and 7 x 3/4" - 1 1/4"	
- Tyres	<input checked="" type="checkbox"/> air-filled	<input checked="" type="checkbox"/> accident proof
- Adjustability	<input checked="" type="checkbox"/> yes	<input type="checkbox"/> no
- Kind of fixation	<input type="checkbox"/> quick release	<input checked="" type="checkbox"/> rigid
Brake system	Knee lever brakes	
Frame material	Aluminium	
Arm rests	Height adjustable and removable	
Foot rests	Height adjustable, foot plate(s) can be swivelled upwards	Height adjustable, can be turn aside and removable, foot plate(s) can be swivelled upwards
Anti tip device	existing	

Department: MHS 5 / me
Date: 2009-11-25

Page 2 of 2