

# Certificate

No.: Z1 01 09 42097 022



Invacare Rea AB  
Org. No. 556098-9211  
Växjövägen 303  
343 21 Diö  
Sweden

with production facility(ies)  
42097

is authorized to label the following products with the  
**TÜV Mark P.**  
Please see notes overleaf.

**Product:** Rollstühle für Behinderte  
Wheelchair for handicapped

**Model:** Manual wheelchair  
Type: WHEELER TRANSPORT

**Parameters:**

seat width:	260, 290, 320, 350 mm
seat depth:	260 to 380 mm
seat height:	390 to 500 mm
backrest height:	260 to 380 mm
max. width:	570 mm
max. depth:	800 mm
frame material:	aluminium, colour coated
tyres, front:	140 x 40 mm (massive)
tyres, rear:	16" (air-filled)
brake system:	drum brake
max. load:	60 kg

The compliance of the above-mentioned product with the following essential requirements was tested on a voluntary basis:

EN 12182: 1999  
TÜV PRODUCT SERVICE PPP 31010  
"manual wheelchairs" 01/00

Report No: 70019817

Released with the above certificate number by the  
certification body of TÜV PRODUCT SERVICE GMBH.

Department: HANMED5 / ut  
Date: 24-Sep-2001

