

Certificate



No.: Z1 01 04 42097 020

Invacare Rea AB
Org. No. 556098-9211
Växjövägen 303
343 21 Diö
Sweden

with production facility(ies)
42097

is authorized to label the following products with the
TÜV Mark P.

Please see notes overleaf.

Product: Rollstühle für Behinderte
Wheelchair for handicapped

Model: Manual wheelchair
Type: WHEELER

Parameters:

seat width:	260, 290, 320, 350 mm
seat depth:	260 to 380 mm
seat height:	390 to 500 mm
backrest height:	260 to 380 mm
max. width:	570 mm
max. depth:	800 mm
frame material:	aluminium, colour coated
tyres, front:	150 x 27 mm (massive)
tyres, rear:	24" (air-filled)
brake system:	knee lever brake or drum brake
max. load:	60 kg

The compliance of the above-mentioned product with the following essential requirements was tested on a voluntary basis:

EN 12182: 1999
TÜV PRODUCT SERVICE PPP 31010
"manual wheelchairs" 01/00

Report No: 70012550

Released with the above certificate number by the
certification body of TÜV PRODUCT SERVICE GMBH.

Department: HANMED5 / me
Date: 03-Apr-2001



Handwritten signature

ZERTIFIKAT ♦ CERTIFICATE ♦ CERTIFICADO ♦ CERTIFICAT ♦ 認 証 証 書 ♦ CERTIFICATE ♦ ZERTIFIKAT ♦