

Dyna-Form® Jupiter

Static mattress range



A cost-effective solution

Single-sided castellated foam



Affordable quality mattress

Comfort & Support

- Castellated foam with tear duct cut: smaller cells at the heels, larger cells sacrally for targeted pressure redistribution.
- Castellated foam with independent cells offer greater immersion and envelopment.

Multi-Directional Stretch

- Multi-stretch, vapour permeable fabric impregnated with an anti-fungal agent reduces shear and friction forces.
- Fully sewn or high frequency welded cover prevents fluid ingress and supports infection control and decontamination.
- Launderable up to 80 °C.

Quality Assured

- CME foam and fabrics for durability and long-term value.
- Fully replaceable parts extend product life and reduce total cost of ownership.

Technical data

Height	Width	Length	Max. User Weight	Total weight	Max. Cleaning Temperature
120, 140, 150 mm	760-1200 mm	1790-2050 mm	165 kg	11 kg	80 °C

Manufactured to: BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177 (ignition Source 5)

Essential nursing care is pivotal in pressure injury management. This mattress contributes to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries.

Contact us:
Invacare International GmbH
Neuhofweg 51
4147 Aesch BL, Switzerland
globalcomms@dhcare.com
www.dhcare.com

Legal Manufacturer:
DHG Ltd.
Withey Court
Western Industrial Estate
Caerphilly, United Kingdom
CF83 1BF

Invacare International GmbH is a member of the DHCare group of companies. Invacare International GmbH reserves the right to modify existing product models at any time. Product specifications may vary across regions. ©2026 Invacare International GmbH. All rights reserved. Dyna-Form_Jupiter_EU_29102026